



● Comprehensive non-targeted chemical fingerprinting of coffee silverskin extracts with MRMS

Coffee silverskin is a major coffee bean roasting by-product, which is currently underutilized and mainly discarded as industrial waste. However, silverskin is a rich source of polyphenols and other bioactive ingredients, and thus a potential feed-stock for pharmaceutical, cosmetics and food sectors.

Abstract

Coffee silverskin, a by-product of the coffee roasting process, contains a plethora of bioactive compounds, such as caffeine, lipids, chlorogenic acids, and

melanoidins, which possess a considerable potential in several industrial applications. In this study, a comprehensive non-targeted chemical fingerprinting of solvent extracts of pelletized coffee silverskin residue was

performed by using ultrahigh-resolution MRMS technology, giving access to the hundreds of chemical constituents, including organic acids, polyphenols, sugars, and nitrogen-containing heterocycles.

Keywords:
coffee, silverskin,
polyphenol, MRMS,
metabolomics




Figure 1: Selected APPI spectra in positive ion mode of coffee silverskin extracts using different extraction solvents.

Introduction

Coffee and other caffeinated drinks are the most consumed beverages in the world. The highest consumption per capita occurs in Finland, where about 10 kg of roasted coffee beans are consumed per person a year. Coffee is enjoyed mainly due to its taste, habit and its stimulating effect, caused mainly by caffeine, while it contains no essential nutrients. Moderate coffee consumption has also been associated with the reduced risk of neurodegenerative diseases, like Parkinson's disease, type II diabetes, and some types of cancer [1]. Coffee silverskin (CS) is the thin outermost layer of the coffee bean and the only by-product obtained from the coffee roasting process. Over 10 million tons of coffee is roasted globally every year [2], leading to the estimated CS production of about 200,000 tons.

Despite the huge availability of this feedstock, its current utilization is very limited, and it is mostly used as a solid fuel or soil fertilizer. However, CS is a rich source of polyphenols, lipids, and other bioactive compounds, and thus its revalorization has gained more attention recently [3,4]. The most interesting compounds include chlorogenic acids (CGAs) and melanoidins, which could find use in pharmaceutical, cosmetics, food, and techno-chemical sectors [5-7]. Compounds from CS can be recovered by various methods such as hot water or solvent extraction. Due to the chemical complexity of CS extracts, rapid and sensitive analysis methods for their characterization are desired. In this work, we employed ultrahigh-resolution Magnetic Resonance Mass Spectrometry (MRMS) technology for comprehensive, non-targeted chemical fingerprinting of the solvent extracts of CS.

Materials and Methods

Solvent extraction

The coffee silverskin pellets, kindly provided by Meira roastery (Helsinki, Finland), were subjected to continuous Soxhlet extraction. The solvents used were chloroform, dichloromethane, hexane, toluene, acetone, acetonitrile, methanol, ethanol, and water (HPLC grade). Both, non-polar and polar solvents were used to assess their ability to extract different types of compounds from the silverskin pellets. Prior to the mass spectrometric analysis, the obtained extracts were further diluted with methanol (for negative-ion ESI) or a methanol/toluene mix (10:1, v/v) (for positive-ion APPI), to the approximate concentration of 50 - 100 $\mu\text{g}/\text{mL}$.

MS analysis

All mass spectrometric analyses were performed on a solarix 12T XR MRMS instrument (Bruker Daltonik GmbH, Bremen, Germany), equipped with a dynamically harmonized ICR cell (ParaCell). All samples were directly infused to an Apollo-II ESI and APPI ion source by using a syringe pump (flow rate of 6.7 $\mu\text{L}/\text{min}$ for APPI and 2 $\mu\text{L}/\text{min}$ for ESI). The ions were detected at the m/z range of 92 - 2000 with a mass resolving power of $\sim 530,000$ at m/z 300 (transient length 1.05 s). Two hundred single scans were co-added for each mass spectrum and processed in magnitude mode. The instrument was externally calibrated by sodium trifluoroacetate (NaTFA) clusters. Furthermore, internal mass recalibration was accomplished with a custom-made calibration list, containing known analytes from different compound classes. For the elemental formula search (SmartFormula), the following parameters were used: mass error ≤ 0.8 ppm; relative intensity $\geq 0.1\%$;

signal-to-noise (S/N) ratio ≥ 5.0 ; H/C ratio ≤ 3 ; DBE ≤ 80 ; elemental formula $^1\text{H}_{1-200}$ $^{12}\text{C}_{1-100}$ $^{16}\text{O}_{1-25}$ $^{14}\text{N}_{1-5}$ $^{32}\text{S}_{1-2}$.

Data processing and structural annotations

The mass spectra were processed using DataAnalysis 5.0 (Bruker Daltonik GmbH, Bremen, Germany). Further analysis and structural annotations were accomplished by using MetaboScape 5.0 (Bruker Daltonik GmbH, Bremen, Germany) with CompoundCrawler database search engine.

MS analysis

Ultrahigh-resolution MRMS represents an unparalleled analysis tool for non-targeted chemical fingerprinting of natural extracts and other complex mixtures, giving access to hundreds or even thousands of analytes in a single measurement. When combined with different ionization techniques such as electrospray ionization (ESI) or atmospheric pressure photoionization (APPI), both

polar and non-polar analytes can be detected.

Based on the data, all CS extracts were highly complex with up to ~ 4600 and ~ 2200 spectral features detected with ESI negative ion mode and APPI positive ion mode, respectively. Figure 1 shows selected APPI spectra for three different solvents. The most abundant compounds detected with APPI in positive ion mode included different acids, di- and triterpenoids, sterols, phenolic acids, and nitrogen heterocycles (Figure 2). The most abundant compound detected with APPI was caffeine ($\text{C}_8\text{H}_{10}\text{N}_4\text{O}_2$), which was present in both non-polar and polar





Figure 3: Example of nitrogen-containing alkaloids found in coffee silverskin extracts: caffeine (left), trigonelline (right).



m/z meas.	M meas.	Ions	Name	Molecular For...	Annotations	AQ	Boxplot	Chloroform	DCM_Coffe...	Toluene_Ext...	HEXANE_Co...	Ethanol_Cof...	Methanol_C...	Water_Coffe...	aCETONE_C...
1	330.21892	330.21947	16-O-methylcafestol	C ₂₁ H ₃₀ O ₃	M										
2	399.36232	398.35505	24-methylenecholesterol	C ₂₈ H ₄₈ O	M										
3	153.05462	152.04734	3',4'-dihydroxyacetophenone	C ₉ H ₈ O ₃	M										
4	208.07307	208.07361	3,4-dimethoxycinnamic acid	C ₁₁ H ₁₂ O ₄	M										
5	139.07528	138.06800	3-ethylcatechol	C ₉ H ₁₀ O ₂	M										
6	159.06784	159.06838	6-methoxyquinoline	C ₁₀ H ₉ NO	M										
7	411.36212	410.35484	7-dehydrostigmaterol	C ₂₉ H ₄₈ O	M										
8	430.38057	430.38112	alpha-tocopherol	C ₂₉ H ₅₀ O ₂	M										
9	321.20610	320.19882	atractyligenin	C ₁₉ H ₃₂ O ₄	M										
10	316.20336	316.20391	cafestol	C ₂₉ H ₄₈ O ₃	M										
11	554.43265	554.43320	cafestol palmitate	C ₃₉ H ₇₀ O ₄	M										
12	181.04959	180.04232	caffeic acid	C ₉ H ₈ O ₄	M										
13	194.07986	194.08041	caffeine	C ₈ H ₁₀ N ₄ O ₂	M										
14	111.04403	110.03675	catechol	C ₆ H ₆ O ₂	M										
15	291.23196	290.22469	coffeeol	C ₁₉ H ₃₀ O ₂	M										
16	426.38551	426.38606	cycloartenol	C ₃₀ H ₅₀ O	M										
17	298.19269	298.19324	dehydrocafestol	C ₂₀ H ₃₂ O ₂	M										
18	296.17705	296.17760	dehydrokahweol	C ₂₀ H ₃₂ O ₂	M										
19	183.06516	182.05788	dihydrocaffeic acid	C ₉ H ₁₀ O ₄	M										
20	289.25260	288.24532	ent-kaur-16-en-19-ol	C ₃₀ H ₅₂ O	M										
21	194.05737	194.05791	ferulic acid	C ₁₉ H ₁₈ O ₄	M										
22	314.18782	314.18836	kahweol	C ₂₀ H ₃₂ O ₃	M										
23	552.41691	552.41746	kahweol palmitate	C ₃₉ H ₇₀ O ₄	M										
24	275.27339	274.26611	kaurane	C ₂₉ H ₄₈	M										
25	210.07477	210.07532	liberine	C ₉ H ₁₀ N ₄ O ₃	M										
26	124.03923	123.03196	nicotinic acid	C ₆ H ₇ NO ₂	M										
27	414.38559	414.38614	sitosterol	C ₂₉ H ₅₀ O	M										
28	396.37506	396.37561	stigmastan-3,5-diene	C ₂₉ H ₄₈	M										
29	412.36999	412.37054	stigmaterol	C ₂₉ H ₄₈ O	M										
30	181.07200	180.06473	theophylline	C ₇ H ₈ N ₄ O ₂	M										
31	138.05483	137.04755	trigonelline	C ₇ H ₉ NO ₂	M										

Figure 2: A table showing the most abundant compounds detected in the coffee silverskin extracts with APPI in positive ion mode and their relative abundance variation across different samples (different extraction solvents).




Figure 4: A table showing the most abundant compounds detected in the coffee silverskin extracts with ESI in negative ion mode and their relative abundance variation across different samples (different extraction solvents).




Figure 5: PCA analysis (scores plot) of the coffee silverskin extracts based on ESI data in negative ion mode. Non-polar and polar solvents are clearly separated in the scores plot while water stands out of the two groups.

solvents. The other nitrogen-containing heterocycles were 6-methoxyquinoline, liberine, theophylline and trigonelline, all naturally occurring alkaloids in coffee beans (Figure 3). The other abundant compounds

included cafestol, kahweol, and their dehydro-forms, which are diterpenoids that have been associated with a variety of pharmacological effects of coffee. This confirms that CS is a rich source of valuable

bioactive compounds. Since APPI does not efficiently ionize some of the more polar, aliphatic or alicyclic compounds, complementary data were acquired from the extracts by using ESI in negative ion mode. The

most abundant compounds detected with ESI included different acids, carbohydrates, and their derivatives (Figure 4). Among fatty acids, linoleic (C18:2), palmitic (C16:0), oleic (C18:1), and arachidic (C20:0) acids were detected. A plethora of chlorogenic acids (i.e. quinic and caffeic acids and their esters) were also observed. Carbohydrates (e.g. galactose and arabinose) were more enriched in the polar solvents.

To assess the overall impact of different solvents on the chemical composition of extracts, two-dimensional principal component analysis (PCA) of the ESI data was performed in MetaboScape. Figure 5 depicts the PCA scores plot, showing that non-polar as well as polar solvents are grouped and clearly

separated from each other, while water stands out of the two solvent types, mainly due to higher content of nitrogen-containing analytes. Therefore, by choosing an appropriate solvent, specific types of compounds can be recovered from coffee silverskin for possible further applications.

Acknowledgements

Meira Ltd. (Helsinki, Finland) is thanked for providing the coffee silverskin sample for this study.

Conclusion

- Ultrahigh-resolution MRMS represents a powerful tool for non-targeted chemical fingerprinting of complex mixtures, exemplified here for coffee silverskin extracts. Very simple sample preparation protocols, unparalleled data acquisition speed, and confident assignment of chemical formulae for hundreds or even thousands of analytes in a single mass spectrum are the key analytical benefits of this technology.



Learn More

You are looking for further Information?
Check out the link or scan the QR code for more details.

www.bruker.com/solarix



References

- [1] Grosso G, Godos J, Galvano F, Giovannucci EL (2017). *Coffee, Caffeine, and Health Outcomes: An Umbrella Review*. *Annu. Rev. Nutr.*, **37**, 131–56.
- [2] International Coffee Organization (www.ico.org)
- [3] Narita Y, Inouye K. (2014). *Degradation Kinetics of Chlorogenic Acid at Various pH Values and Effects of Ascorbic Acid and Epigallocatechin Gallate on Its Stability under Alkaline Conditions*. *Food Res. Int.*, **61**, 16-22.
- [4] Toschi TG, Cardenia V, Bogana G, Mandrioli M, Rodriguez-Estrada MT (2014). *Coffee Silverskin: Characterization, Possible Uses, and Safety Aspects*. *J. Agric. Food Chem.*, **62**, 10836–10844.
- [5] Bessada SMF, Alves RC, Oliveira MBPP (2018). *Coffee Silverskin: A Review on Potential Cosmetic Applications*. *Cosmetics*, **5**, 5.
- [6] Tores de la Cruz S, Iriondo-DeHond A, Herrera T, Lopez-Tofiño Y, Galvez-Robleño C, Prodanov M, Valezquez-Escobar F, Abalo R, Dolores Del Castillo M (2019). *An Assessment of the Bioactivity of Coffee Silverskin Melanoidins*. *Foods*, **8**, 68.
- [7] Barbero-López A, Monzó-Beltrán J, Virjamo V, Akkanen J, Haapala A (2020). *Revalorization of Coffee Silverskin as a Potential Feedstock for Antifungal Chemicals in Wood Preservation*. *Int. Biodeterior Biodegradation*, in press.

For Research Use Only. Not for Use in Clinical Diagnostic Procedures.

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0

Bruker Scientific LLC

Billerica, MA · USA
Phone +1 (978) 663-3660

ms.sales.bdal@bruker.com – www.bruker.com