

A GENERIC KIT-BASED APPROACH FOR LC-MS/MS QUANTIFICATION OF URINARY ALBUMIN FOR CLINICAL RESEARCH

Mary E. Lame, Caitlin M. Dunning, and Erin E. Chambers
Waters Corporation 34 Maple Street Milford, MA 01757

INTRODUCTION

Albumin (~MWT 66.5 kDa) is the most abundant protein in blood.¹ Its presence in urine is often one of the first signs of kidney damage. Thus, it has become an important biomarker for renal disease.²⁻⁴ As a result, accurate measurement in urine is of high interest in drug discovery and clinical research.

Under normal kidney function, urinary albumin levels are quite low (< 30 mg/day), but following renal injury albumin levels in urine can exceed 300 mg/day.²⁻⁴ Existing affinity-based methods for urinary albumin quantification include: immunoturbidimetric, ELISA and radioimmunoassay.⁵⁻⁸ While these immunoassay (IA) methods are sensitive and simple to execute, poor reagent reproducibility, lack of standardization, cross-reactivity, limited linear dynamic range, and other shortcomings have led to increased interest in LC-MS based methods. With its many benefits (e.g., multiplexing, selectivity, dynamic range and fast method development), LC-MS is widely accepted for protein quantification. However, this workflow can be complex and time consuming, often taking 24 hours to achieve analytically sensitive and accurate quantification.

This work describes the accurate quantification of urinary albumin over 3.5 orders of magnitude (0.1—500 µg/mL), which is completed in <4 hours using commercially available digestion and peptide purification kits and generic protocols from only 15 µL of sample. This first standardized analytical method is shorter (9X) and more sensitive (50X) than published methods.^{3,8}

METHODS

Sample Preparation

Human serum albumin (HSA) was spiked into normal human urine. Urine samples (15 µL) were prepared for LC-MS analysis using the ProteinWorks eXpress Digest Kits and Protocols*. After digestion, peptides were cleaned-up using the ProteinWorks µElution SPE Clean-up Kit and Protocol.

LC-MS Conditions

LC-MS/MS quantification of signature peptides was performed using a Waters Xevo TQ-XS triple quadrupole MS (ESI+). Chromatographic separation was achieved using an ACQUITY UPLC system with an ACQUITY UPLC Peptide BEH C18, 300Å, 1.7 µm, 2.1 mm x 150 mm column and 0.1% formic acid in water and acetonitrile mobile phases. A total of 11 HSA tryptic peptides (highlighted in blue, Figure 1), and 2 MRM transitions per peptide, were monitored for quantification. MS conditions are summarized in Table 1.




Figure 1. Amino acid sequence of human serum albumin⁹; Tryptic peptides used for quantification are highlighted in blue.

Peptide	Precursor Charge State	MRM Transition	Collision Energy (eV)	Product Ion Identification
YLVEIAR	[M+2H] ²⁺	464.25 > 277.16 464.25 > 651.35	22	[1H+] ¹ /b2 [1H+] ¹ /y5
FQNALLVR	[M+2H] ²⁺	480.79 > 276.13 480.79 > 685.44	22	[1H+] ¹ /b2 [1H+] ¹ /y6
LVNEVTEFAK	[M+2H] ²⁺	575.31 > 694.38 575.31 > 937.46	25	[1H+] ¹ /y6 [1H+] ¹ /y8
VFDEFKPLVEEPQNLIK	[M+3H] ³⁺	682.37 > 712.44 682.37 > 970.52	25	[1H+] ¹ /y6 [1H+] ¹ /y8
AEFAEVSK	[M+2H] ²⁺	440.42 > 201.09 440.42 > 680.36	22	[1H+] ¹ /b2 [1H+] ¹ /y6
LCTVATLR	[M+2H] ²⁺	467.25 > 274.12 467.25 > 660.40	22	[1H+] ¹ /b2 [1H+] ¹ /y6
DLGEENFK	[M+2H] ²⁺	476.23 > 229.12 476.23 > 723.33	22	[1H+] ¹ /b2 [1H+] ¹ /y6
TYETTTLEK	[M+2H] ²⁺	492.75 > 265.12 492.75 > 720.38	23	[1H+] ¹ /b2 [1H+] ¹ /y6
QATALVELVK	[M+2H] ²⁺	500.81 > 488.31 500.81 > 587.38	22	[1H+] ¹ /y4 [1H+] ¹ /y5
LVAASQAALGL	[M+2H] ²⁺	507.30 > 189.12 507.30 > 712.40	22	[1H+] ¹ /y2 [1H+] ¹ /b8
RPCFSALEVDETYVPK	[M+3H] ³⁺	637.64 > 244.17 637.64 > 961.46	22	[1H+] ¹ /y2 [1H+] ¹ /b8

Table 1. Final MS conditions for HSA tryptic peptides, including precursor and fragment ions; Primary tryptic peptides used for quantification are highlighted in blue.

RESULTS

Peptide	Curve (µg/mL)	Weighting	Linear Fit (R ²)	Mean % Accuracy	% Accuracy Range
YLVEIAR	0.1-500	1/X	0.998	100.0	92.3-114.2
FQNALLVR			0.999		91.4-112.9
LVNEVTEFAK			0.997		91.8-111.8
VFDEFKPLVEEPQNLIK			0.998		93.1-113.0

Table 2. Linear dynamic range and standard curve statistics for the 4 primary albumin tryptic peptides: YLVEIAR, FQNALLVR, LVNEVTEFAK, and VFDEFKPLVEEPQNLIK used for quantification. Urine samples were digested and extracted protein quantification digestion and peptide purification kit.




Figure 2. UPLC chromatographic separation of albumin tryptic peptides, digested in human urine (50 µg/mL).

Peptide	Albumin QC overspike concentration (µg/mL)	Urine Lot #1		Urine Lot #2		Urine Lot #3	
		Mean (N=3) % Accuracy	% CV	Mean (N=3) % Accuracy	% CV	Mean (N=3) % Accuracy	% CV
YLVEIAR	0.000	100.0	0.5	100.0	4.4	100.0	4.0
	0.150	96.6	1.7	101.9	2.9	103.0	1.7
	0.400	92.2	0.4	103.1	3.5	101.7	1.2
	2.000	87.2	1.4	104.9	2.0	102.7	4.4
	20.000	85.3	0.7	97.0	8.1	99.8	4.0
	80.000	94.6	3.9	101.0	5.3	99.6	4.1
FQNALLVR	0.000	100.0	5.0	100.0	1.9	100.0	4.0
	0.150	99.2	0.6	101.7	2.3	101.0	1.0
	0.400	93.6	1.2	104.3	5.7	103.8	0.7
	2.000	88.2	2.0	101.8	4.4	100.8	1.5
	20.000	85.8	0.3	98.8	6.4	97.8	5.7
	80.000	96.3	1.6	101.7	2.4	99.3	1.7
LVNEVTEFAK	0.000	100.0	1.0	102.8	3.9	99.3	2.3
	0.150	100.0	4.7	100.0	4.6	100.4	3.2
	0.400	95.0	2.0	102.5	2.8	103.8	1.9
	2.000	91.5	1.5	105.3	1.4	103.1	3.1
	20.000	86.7	0.8	106.4	0.4	107.5	0.5
	80.000	85.6	0.6	98.0	6.1	98.4	3.0
VFDEFKPLVEEPQNLIK	0.000	100.0	2.0	100.0	2.4	100.0	4.3
	0.150	94.9	3.7	103.3	3.0	102.0	1.3
	0.400	90.1	1.6	106.3	3.3	102.9	1.3
	2.000	86.5	1.0	107.1	2.9	103.9	3.5
	20.000	86.7	0.5	98.9	7.5	101.2	2.5
	80.000	95.1	4.9	104.2	2.7	98.7	2.6

Table 3. Summary of the QC sample statistics (precision and accuracy) in Urine Lots #1-3 for the 4 primary tryptic peptides used to quantify albumin.

Table 4. Calculated endogenous urinary albumin concentrations in 3 lots of human urine using the YLVEIAR, FQNALLVR, LVNEVTEFAK, and VFDEFKPLVEEPQNLIK tryptic peptides of HSA.

Urine Lot #	Peptide	Mean (N=3) Calculated Endogenous Concentration (µg/mL)	Mean (N=3) Calculated Endogenous Concentration (nM)	% CV
1	YLVEIAR	2.520	37.9	0.5
	FQNALLVR	2.069	31.2	5.0
	LVNEVTEFAK	2.428	36.5	4.7
	VFDEFKPLVEEPQNLIK	1.486	22.4	2.0
	Mean Albumin Concentration (4-peptides)	2.126	32.0	3.0
2	YLVEIAR	11.544	173.8	4.4
	FQNALLVR	11.129	167.5	1.9
	LVNEVTEFAK	11.159	168.0	4.6
	VFDEFKPLVEEPQNLIK	10.241	154.2	2.4
	Mean Albumin Concentration (4-peptides)	11.018	165.9	3.3
3	YLVEIAR	6.668	100.4	4.0
	FQNALLVR	6.193	93.2	4.0
	LVNEVTEFAK	6.383	96.1	3.2
	VFDEFKPLVEEPQNLIK	5.715	86.0	4.3
	Mean Albumin Concentration (4-peptides)	6.240	93.9	3.9




Figure 3. Representative QC chromatograms for the 4 primary HSA peptides: YLVEIAR (A), FQNALLVR (B), LVNEVTEFAK (C), and VFDEFKPLVEEPQNLIK (D) used to quantify albumin in human urine.

DISCUSSION

Using a generic kit-based approach (with simple step-wise protocols and standardized, pre-measured reagents) for digestion and subsequent peptide purification, accurate albumin quantification in urine was achieved.

- Total sample preparation time, including SPE, was <4 hours.
- UPLC chromatographic separation of albumin tryptic peptides, digested in human urine, is illustrated in Figure 2. Peak widths for all peptides were <4.5 seconds wide.
- Through direct digestion (no affinity purification) of 15 µL of urine and subsequent peptide purification, quantification limits of 0.1—500 µg/mL was readily achieved for the four primary HSA signature tryptic peptides (Table 2). In addition, standard curves were linear over 3.5 orders of magnitude with mean accuracies ranging from 91.4—114.2%.
- Employing a mixed-mode SPE clean-up step, effectively removing buffer salts and excess digestion reagents, while also concentrating the sample, significantly improved accuracy and precision of the QCs across all lots of urine tested.
- The accuracy and precision for HSA quantified in human urine (3 lots) was excellent with accuracies ranging from 85.3—107.5% and CVs <8.1%. HSA QC statistics are highlighted in Table 3, and illustrated in Figure 3 for the 4 primary HSA peptides.
- Endogenous urinary albumin concentrations were accurately quantified in three individual lots of urine and are summarized in Table 4.

CONCLUSION

Using a generic kit-based approach for digestion and subsequent peptide purification, accurate and reproducible quantification of endogenous urinary albumin was achieved. The analytical sensitivity (0.1 µg/mL), broad linear dynamic range, and selectivity of this LC-MS method reliably quantifies both low endogenous and elevated urine levels that would be expected in normal and disease populations while approaching the analytical sensitivity of immunoassays.

References

- Wikipedia contributors, Human serum albumin. *Wikipedia, The Free Encyclopedia*. https://en.wikipedia.org/w/index.php?title=Human_serum_albumin&oldid=750699486 (accessed 27Feb2017).
- Tesch GH, Review: Serum and urine biomarkers of kidney disease: A pathophysiological perspective. *Nephrology*. 2010; 15: 609-616.
- Beasley-Green A, Burris NM, Bunk DM, Phinney KW, Multiplexed LC-MS/MS Assay for Urine Albumin. *J. Proteome Res*. 2014; 13: 3930-3939.
- deJong PE, Gansvoort RT, Bakker SJ, Macroalbuminuria and microalbuminuria: Do both predict renal and cardiovascular events with similar strength? *Nephrol*. 2007; 20:375-380
- Polkinghorne KR, Detection and measurement of urinary protein. *Curr. Opin Nephrol. Hypertens*. 2006;15: 625-630
- Choi S, Choi EY, Kim HS, Oh SW, On-site quantification of human urine albumin by a fluorescence immunoassay. *Clin Chem*. 2004; 50:1052-1055
- Seegmiller, JC, Sviridov D, Larson TS, Borland TM, Hortin GL, Liske JC, Comparison of urine albumin quantification by immunoturbidimetry, competitive immunoassay, and protein-cleavage liquid chromatography-tandem mass spectrometry. *Clin. Chem*. 2009; 55:1991-1994
- Cell Biolabs, Inc. (San Diego, CA) Human Albumin ELISA Kit (STA-383) Product manual
- UniProt, (16May2014). UniProtKB - P02768 (ALBU_Human). Primary accession number:P02768 accessed 27Feb2017 from <http://www.uniprot.org/uniprot/P02768>