## **CORRELATING DIFFERENTIAL ENDOGENOUS METABOLITE PROFILES WITH THE PHARMACOKINETICS OF GEFITINIB** AND ITS ASSOCIATED DRUG METABOLITES USING AN ION MOBILITY BASED APPROACH

Lee A. Gethings<sup>1</sup>, Martin Palmer<sup>1</sup>, Adam King<sup>1</sup>, Lauren Mullin<sup>1</sup>, Ian D. Wilson<sup>2</sup>, Robert S. Plumb<sup>3</sup> <sup>1</sup>Waters Corp., Wilmslow, UK; <sup>2</sup>Imperial College, London, UK; <sup>3</sup>Waters Corp., Milford, MA, USA

## INTRODUCTION

Gefitinib (iressa), a selective inhibitor of the epidermal growth factor receptor (EGFR; HER1) tyrosine kinase, used for treating non-small cell lung cancer, is extensively metabolized in animals and humans. However, the methods used to determine its metabolic fate employed HPLC separations of 15-45 minutes duration, which are incompatible with modern high-throughput drug discovery. However, technologies such as UPLC, ion mobility (IM) spectrometry and high-resolution MS (HRMS), when combined, can achieve high-throughput without compromising quality. Thus, UPLC provides excellent separation efficiency allowing rapid analysis, IM offers an additional separation, resolving co-eluting interferences, and structurally relevant CCS values, whilst HRMS provides excellent data for metabolite identification. Here we applied UPLC-IMS-MS to characterize and identify Gefitinib-related metabolites and endogenous metabolite changes in the mouse urine



Urine analysis using for the orally dosed group (PO) using the HILIC-IMS-MS methodology allowed for the identification of Gefitinib-related metabolites and endogenous species. An example drug metabolite originating from Gefitinib is at m/z 449, providing the characteristic pharmacokinetic profile. The box-whisker plot (A) shows the elevated intensity over time, reaching a C<sub>max</sub> at 3-8 hrs before returning to comparative initial levels at 8-24 hrs. Unsupervised PCA plots (pareto scaling applied) representing the PO group with Gefitinib and its related metabolites subtracted from the data (B). The trend observed for changes at the endogenous metabolite level reflects that of the expected pharmacokinetic profile as shown with (A) (i.e. deviation from the predose level before returning close to original, predose levels).





- 2. Giles K, Ujma J, Wildgoose J, et al. A Cyclic Ion Mobility-Mass Spectrometry System. Anal Chem. 2019;91(13):8564-8573
- 3. Chong, J et al. Using MetaboAnalyst 4.0 for Comprehensive and Integrative Metabolomics Data Analysis. 2019, Current Protocols in Bioinformatics 68, e86

/aters THE SCIENCE OF WHAT'S POSSIBLE.™