

Lipid Nanoparticle Impurity Monitoring Using Single Quadrupole Mass Detection for Regulated Environments

Authors: *Duanduan Han, Kellen DeLaney, Robert E. Birdsall, Ying Qing Yu*
Affiliation: Waters Corporation, Milford, MA

Introduction

Lipid nanoparticles (LNPs) represent a novel solution to deliver gene-based therapeutics that is challenging to execute. LNPs are a multi-component biologic comprised of lipid species with varying properties impacting chromatography and detector response. New approaches are needed for characterization and monitoring to reduce the associated risks and ensure active pharmaceutical ingredients are delivered in a safe and efficacious manner. We present an LC-based optical/MS dual-detector platform with improved sensitivity/diagnostic power for LNP workflows. A single quadrupole mass detector provides complementary mass data for raw material impurity screening, MS spectral library-based lipid confirmation, and stability monitoring. Supported by compliant-ready informatics, this workflow is easy to deploy in both non-regulated and regulated environments alike for efficient method development and migration.

Lipid Nanoparticles

Composition

Ionizable lipid (SM102)
• Binds & protects RNA
• Facilitates endosomal escape
• Neutral charge at pH 7.4

50%

38.5%

Cholesterol
• Important for endocytosis

Phospholipid (DSPC)
• Structural stability

10%

1.5%

PEGylated lipid (DMG-PEG2000)
• Stabilizes particle
• Protect from opsonization

Figure 1. Four lipid components of LNP. Cholesterol, DSPC and DMG-PEG2000 provide structural stability and are typically outsourced. The ionizable lipids are often proprietary and produced in-house to protect intellectual property.

LC: ACQUITY™ UPLC™ Premier System

Sample(s):

Cholesterol (CHO)
Distearoylphosphatidylcholine (DSPC)
PEGylated lipid DMG-PEG 2000*
SM-102 (ionizable lipid)
Dlin-MC3-DMA (ionizable lipid)

Injection volume: 3 µL

Sample temp.: ambient

Column: ACQUITY Premier CSH™ Phenyl-Hexyl Column

1.7 µm, 2.1 mm X 50 mm (p/n:186009474)

Column temp.: 50 °C

MS: ACQUITY QDa™ Mass Detector

Scan range: 150-840 m/z @ 5Hz

Capillary voltage: 1.5 kV, cone voltage: 15 V, probe temp.: 600 °C

Methods

MP A: H₂O, 0.4% formic acid v/v (MS-grade)
MP B: MeCN, 0.6% formic acid v/v (MS-grade)

Gradient:

Time	Flow (mL/min)	%A	%B	Curve
Initial	0.4	40	60	Initial
6.0	0.4	10	90	6
8.0	0.4	10	90	6
8.5	0.4	40	60	6
12.0	0.4	40	60	6

Detection




Figure 2. Detector type and response. Detector response is shown across 3 detector types for a panel of lipids representative of components used in the production of lipid nanoparticles. Peaks 1-5 are SM102, DOTMA, cholesterol, DMG-PEG2000, and DSPC. SM102, cholesterol, DMG-PEG2000, and DSPC were prepared in a molar ratio of 50:38.5:1.5:10 representative of a formulated LNP sample. DOTMA was spiked-in at half the concentration relative to SM102.

Results

Enhanced analysis capabilities with complementary mass data

Increased sensitivity




Figure 3. ELSD meets the analytical needs of labs, detecting lipid mass loads as low as 15 ng. Selected ion recording (SIR) drives down LOD to 1.5 pg and extends LOD by 4 orders which is suitable to assess trace impurities. Green shaded concentrations represent the dynamic range. Sample: ionizable lipid DOTMA.

Drug formulation and stability





Figure 4. Labile bonds such as esters are subject to degradation. An example of which is shown here for the hydrolysis of DSPC under high pH. MS provides mass data for tentative assignment of degradant peaks as the methylated fatty acid methyl stearate (299.2 m/z) and the polar head group (258.0 m/z).

Product / process related impurities



Compound	m/z	Proposed cause
[MC3 -2H +H] ⁺	640.6	Desaturation
[MC3 +H] ⁺	642.6	Main peak
[MC3 +2H +H] ⁺	644.6	Saturation
[MC3 -2H +O +H] ⁺	656.7	Desaturation/Oxidation
[MC3 +O +H] ⁺	658.8	Oxidation
[MC3 +2H +O +H] ⁺	660.8	Saturation/Oxidation
[MC3 +2O +H] ⁺	674.6	Saturation/Double oxidation

Product consistency






Figure 6. MS data reveal volatile impurities (up to 5% peak area) which otherwise would have gone undetected in the ELSD, requiring 9 times the mass load to be detected. Sample: DSPC, typically outsourced by LNP manufacturers.

Leverage MS to reduce risks in process development and manufacturing workflows

a) ELSD/MS workflow



Empower™ 3 Software




Cholesterol 369.3 m/z

SM102 710.7 m/z


DMG-PEG2000 741.8 m/z

DSPC 790.7 m/z

MS confirmation



b) Raw material screening



c) Compositional analysis




Figure 7. Complementary mass data reduces error in process development and manufacturing workflows through spectral confirmation. a) Spectral libraries based on retention time and mass can be utilized to confirm peak identity in raw material screening and compositional analysis. b) To demonstrate this, ELSD and mass data were acquired in the same acquisition run for the DSPC lipid spiked with methyl stearate to act as a low-level impurity present in a raw material shipment. In this example, we can view both the ELSD and QDa information from the results window with pertinent peak information. The table shows best match (and structural information if available) found in the library for the spectra associated with each peak seen in the ELSD chromatogram. c) This workflow is also applicable in a compositional analysis setting that may occur in formulation or release assays as part of QA/QC testing.

Conclusion

- Multi-detector array for efficient development and manufacturing
- MS-support for increased sensitivity and complementary mass data
- Spectral library support for identification and structure assignment
- Compliant-ready software solution to facilitate easy method migration

