

CONVERSION AND INTEGRATION OF OMICS DATA FROM A PROTOTYPE, BENCHTOP MULTI-REFLECTING TIME-OF-FLIGHT (MRT) PLATFORM WITH THIRD-PARTY INFORMATIC WORKFLOWS

Lee A. Gethings, Ian Morns, Pete Reay, Simon Jones, Nyasha Munjoma, Jayne Kirk, Richard Lock
Waters Corp., Wilmslow, Cheshire, United Kingdom

INTRODUCTION

OMIC studies typically generate large and complex data sets, derived from various sample types such as biofluids. As mass spectrometry (MS) technology rapidly advances and with data acquisition methods such as data independent analysis (DIA) progressing, the ability to delve deeper into the metabolome and lipidome is significantly improved, whilst at the same time generating high-dimensional datasets. In this study, we describe a data processing pipeline which converts MS data (DIA and DDA based) that were collected using a novel, multi-reflecting ToF mass spectrometry into a generic file format, that can be readily processed by a variety of third-party informatic tools. Data were primarily processed via Skyline¹ and MS-Dial² for lipidomic datasets, whilst metabolomic datasets were processed using XCMS³.




Figure 1. waters_connect™ and routes of connectivity to third party in-formatic packages.

- mzML output options available within the LC-MS method created within AME. Users have the option of tailoring the mzML format for compatibility to the various third party informatic packages available:
- Automatic generation of mzML during the course of the acquisition to ensure
 - Centroiding function for those applications requiring centroided data. Peak detection is performed with lockmass correction applied.
 - Compression to minimise the size of mzML's (binary data), particularly suitable for larger sample sets and samples of high complexity.
 - Writing index out to file, whereby each mass spectrum and chromatogram are indexed (referenced). An option which is dependent on the third party application used for processing.
 - m/z precision available as 32 or 64-bit versions (dependent on the instrument used for data acquisition).

Once created, the mzML is available in the defined subdirectory folder specified by the user.




Figure 2. Data convert integration into waters_connect. Flexibility is demonstrated with the option of wither generating mzML automatically during data acquisition.




Figure 3. LC-MS analysis of biological samples (urine and plasma) for polar metabolomics and lipidomics. Analytes were chromatographically separated using HILIC (metabolites) or reversed-phase (RP) (lipids) and data collected using the benchtop Xevo MRT mass spectrometer. mzML files were automatically generated via AME during acquisition and processed with a variety of third party informatic tools.




Figure 4. Lipidomic data processing via a variety of third party informatic tools, including Skyline and MS-Dial. Data represented by the Skyline example is based on the EquiSPLASH spiked components, whereas an endogenous lipid (PC(36:1)) is shown when processing with MS-Dial.

RESULTS

Lipidomic data resulting from the plasma extracts were subsequently processed using Skyline and MS-Dial, whilst polar metabolite data relating to the urine samples were processed via XCMS. Exemplar data for each software package is illustrated for the lipidomics (below) and metabolomics (right).

Lipidomics

NIST plasma extracts were spiked with EquiSPLASH™ (Avanti® Polar Lipids, Alabama, USA). Data were imported as mzML format into Skyline and screened against a library consisting of the EquiSPLASH components, highlighting consistency, linear dynamic range and high mass accuracy. The MS-Dial workflow provides a discovery approach with identification of endogenous components, indicating high mass accuracy at the precursor and fragment ion level.

Metabolomics

NIST urine extracts from each of the three levels were pooled to form a QC, which was periodically acquired throughout the LC-MS analysis. Data as mzML format were then uploaded to XCMS for peak picking, statistical interrogation and compound identification.




Figure 5. Representative NIST urine processed through XCMS. Data corresponds to normal human urine collected from nonsmokers with no exposure to tobacco smoke (L1), nonsmokers exposed to tobacco smoke (secondhand) (L2) and smokers who consume at least one pack of cigarettes (L3).

CONCLUSION

- Data Convert provides a flexible and easy approach to generating mzML files from DDA and DIA acquisitions.
- Integration of Data Convert into waters_connect allows a seamless workflow of automatically generating mzML during data acquisition, in addition to having flexibility for users to export mzML from projects at any point in time.
- mzML files can be tailored for the intended third party informatics package (i.e., centroided, compressed etc).
- Data from the benchtop Xevo MRT has demonstrated the workflow from data acquisition through to data processing with third party tools, such as Skyline, MS-Dial and XCMS using the waters_connect/Data Convert pipeline.
- Metabolomic and lipidomic datasets highlight the high quality data that can be generated from the Xevo MRT mass spectrometer:
 - Superior mass accuracy (sub 1 ppm)
 - High mass resolution which can be maintained with fast, high throughput LC gradients
 - High levels of sensitivity/dynamic range

References

- Tsugawa, H. *et al.*, MS-DIAL: data-independent MS/MS deconvolution for comprehensive metabolome analysis. *Nat Methods* **12**, 523-526 (2015).
- Adams, K.J., *et al.*, Skyline for Small Molecules: A Unifying Software Package for Quantitative Metabolomics. *J Proteome Res.* **19**(4), 1447-1458 (2020).
- Tautenhahn R. *et al.*, XCMS Online: A Web-Based Platform to Process Untargeted Metabolomics Data. *Anal. Chem.* **84**(11), 5035-5039 (2012).
- Barupal *et al.*, Generating the Blood Exposome Database Using a Comprehensive Text Mining and Database Fusion Approach. *Environ Health Perspect.* **127**(9), Doi: 10.1289/EHP4713.