

Waters



## ACQUITY UPLC™ Columns

**C<sub>18</sub> • Shield RP<sub>18</sub> • C<sub>8</sub> • Phenyl  
BEH Column Chemistries**

**Acquity**  
Ultra Performance LC™

# Ultra Performance LC™

Waters changed the landscape and future of chromatography with the ACQUITY Ultra Performance LC™ (UPLC™) system. Chromatographers need no longer choose between the speed of short columns and the resolution of long columns. Separations scientists can now enjoy the best of both worlds: speed and resolution with the bonus of increased sensitivity. UPLC™ delivers more information faster—without compromising data integrity.


Waters ACQUITY UPLC™ systems are holistically designed to dramatically improve resolution, sample throughput and sensitivity. Key innovations include:

- Small, pressure tolerant particles
- High pressure fluidic modules
- Minimized system volume
- Negligible carryover
- Reduced cycle times
- Fast response detectors
- Integrated system software and diagnostics


## UPLC™—Speed, Sensitivity and Resolution—Without Compromise

Conditions  
 HPLC Column: XTerra® MS C<sub>18</sub> 2.1 x 150 mm, 5 µm  
 UPLC™ Columns: ACQUITY UPLC™ BEH C<sub>18</sub> 2.1 x 50 mm, 1.7 µm  
 ACQUITY UPLC™ BEH C<sub>18</sub> 2.1 x 100 mm, 1.7 µm  
 Mobile Phase A: 0.1% HCOOH  
 Mobile Phase B: ACN  
 Isocratic Mobile Phase Composition: 75% A; 25% B  
 Flow Rates: HPLC: 0.20 mL/min  
 UPLC™: 0.50 mL/min  
 Injection Volume: 5 µL  
 Sample Diluent: 0.1% HCOOH  
 Sample Concentration: 17 µg/mL  
 Temperature: 50 °C  
 Detection: UV @ 270 nm  
 Instrument: ACQUITY UPLC™ with ACQUITY UPLC™ 2996 PDA

Compounds  
 1. Quercetin  
 2. Kaempferol  
 3. Isohammetin



The differences in the speed, sensitivity and resolution of the UPLC™ vs HPLC separations shown below is graphically represented here.



Whether your objective is maximum speed (with the same resolution) or maximum resolution (in less time), UPLC™ allows you to achieve both – in addition to increased sensitivity – without compromising chromatographic fidelity.


# BEH Technology™

In 1999, Waters launched the XTerra® family of HPLC columns featuring first generation hybrid particle technology (HPT). Waters patented\* HPT enabled XTerra® columns to become the most successful column product in the history of Waters. In HPT, the best properties of inorganic (silica) and organic (polymeric) packings are combined to produce a material that has superior mechanical strength, efficiency, high pH stability and peak shape for bases.



The first generation methyl hybrid particles of XTerra® columns did not possess the mechanical strength or efficiency necessary to realize fully the potential speed, sensitivity and resolution capabilities of UPLC™ technology. Therefore, a new pressure-tolerant particle needed to be created. A new, second generation hybrid material was developed which utilizes a bridged ethyl-siloxane/silica hybrid (BEH) structure. Compared to the first generation methyl hybrid particle of XTerra® columns, the BEH particle of ACQUITY UPLC™ BEH columns exhibits improved efficiency, strength and pH range. BEH Technology™ is a key enabler of the speed, sensitivity and resolution of UPLC™ separations.

## ACQUITY UPLC™ Column Performance—Enabled by BEH Technology™



\*Patent No. 6,686,035 B2

# Why Additional ACQUITY UPLC™ BEH Chemistries?

The fundamental resolution equation for isocratic separations states that resolution ( $R_s$ ) is proportional to the square root of column efficiency ( $N$ ). According to chromatographic theory, column efficiency ( $N$ ) is inversely proportional to particle size ( $d_p$ ). Thus, smaller particles provide higher resolution. The highly efficient 1.7  $\mu\text{m}$  BEH particles allow chromatographers to maximize the efficiency ( $N$ ) of their separation when used in the ACQUITY UPLC™ system.

$$R_s = \frac{\sqrt{N}}{4} \left( \frac{\alpha - 1}{\alpha} \right) \left( \frac{k}{k + 1} \right)$$

System Efficiency      Selectivity      Retentivity





$$N \propto \frac{1}{d_p}$$

*The improved efficiency ( $N$ ) of a UPLC™ separation makes small changes in selectivity ( $\alpha$ ) a very powerful method development tool. These changes are now further enabled by a new expanded range of ACQUITY UPLC™ BEH chemistries.*

However, the fundamental resolution equation also states that by changing the selectivity ( $\alpha$ ) or retentivity ( $k$ ), resolution can also be improved. Lower efficiency HPLC separations required chromatographers to try at random one of the many different types of column chemistries in order to obtain an adequate separation. Since the efficiencies of UPLC™ separations are 2-3 times higher than those of HPLC separations, only a few ligands need to be explored to provide the desired separation.

The chemistries of the ACQUITY UPLC™ BEH column family were carefully chosen to produce the ideal combination of ultra-performance efficiencies, wide pH range and complementary selectivities. When combined with the increased efficiencies of UPLC™ technology, the ACQUITY UPLC™ BEH  $C_{18}$ ,  $C_8$ , Shield RP<sub>18</sub> and Phenyl column chemistries enable the rapid development of faster and more robust separations.

## ACQUITY UPLC™ BEH Chemistry and Particle Information

ACQUITY UPLC™ BEH Chemistry	$C_{18}$	$C_8$	Shield RP <sub>18</sub>	Phenyl
				
<b>Bonded Phase</b>				
<b>Ligand Type</b>	Trifunctional $C_{18}$	Trifunctional $C_8$	Monofunctional Embedded Polar Group	Trifunctional $C_6$ Phenyl
<b>Ligand Density*</b>	3.1 $\mu\text{mol}/\text{m}^2$	3.2 $\mu\text{mol}/\text{m}^2$	3.3 $\mu\text{mol}/\text{m}^2$	3.0 $\mu\text{mol}/\text{m}^2$
<b>Carbon Load*</b>	18%	13%	17%	15%
<b>Endcap Style</b>	Proprietary	Proprietary	TMS	Proprietary
<b>pH Range</b>	1-12	1-12	2-11	1-12
<b>BEH Particle</b>				
<b>Pore Diameter*</b>	130Å	130Å	130Å	130Å
<b>Pore Volume*</b>	0.7 mL/g	0.7 mL/g	0.7 mL/g	0.7 mL/g
<b>Surface Area*</b>	185 $\text{m}^2/\text{g}$	185 $\text{m}^2/\text{g}$	185 $\text{m}^2/\text{g}$	185 $\text{m}^2/\text{g}$
<b>90/10 Ratio*</b>	1.5	1.5	1.5	1.5

\* Expected or Approximate Values




# ACQUITY UPLC™ BEH C<sub>18</sub> and C<sub>8</sub> Columns

The vast majority of reversed-phase LC separations take place on columns that contain C<sub>18</sub> or C<sub>8</sub> bonded stationary phases due to their stability, retentivity and reproducibility. In addition, these hydrophobic ligands most often provide the desired separation. ACQUITY UPLC™ BEH C<sub>18</sub> and C<sub>8</sub> columns were designed to be the universal columns of choice for most UPLC™ separations by providing the widest pH range.

ACQUITY UPLC™ BEH C<sub>18</sub> and C<sub>8</sub> columns incorporate trifunctional ligand bonding chemistries which produce superior low pH stability and ultra-low column bleed. This low pH stability is combined with the high pH stability of the 1.7 μm BEH particle to deliver the widest usable pH operating range. In addition, these new chemistries also utilize new, proprietary endcapping processes which produce outstanding peak shape for bases. These bonding chemistries and particle synthesis innovations produce the sharpest peaks, highest efficiencies and maximum MS sensitivities.

## Explosives Separations on ACQUITY UPLC™ BEH Chemistries





## Ordering Information


ACQUITY UPLC™ BEH Chemistry		Description	Packaging	Part Number
<b>C<sub>18</sub></b>		1.0 x 50 mm, 1.7 μm	Column	186002344
		1.0 x 100 mm, 1.7 μm	Column	186002346
		2.1 x 50 mm, 1.7 μm	Column	186002350
		2.1 x 100 mm, 1.7 μm	Column	186002352
		1.0 x 50 mm, 1.7 μm	3-Pack	176000861
		1.0 x 100 mm, 1.7 μm	3-Pack	176000862
		2.1 x 50 mm, 1.7 μm	3-Pack	176000863
		2.1 x 100 mm, 1.7 μm	3-Pack	176000864
<b>C<sub>8</sub></b>		1.0 x 50 mm, 1.7 μm	Column	186002875
		1.0 x 100 mm, 1.7 μm	Column	186002876
		2.1 x 50 mm, 1.7 μm	Column	186002877
		2.1 x 100 mm, 1.7 μm	Column	186002878
		1.0 x 50 mm, 1.7 μm	3-Pack	176000882
		1.0 x 100 mm, 1.7 μm	3-Pack	176000883
		2.1 x 50 mm, 1.7 μm	3-Pack	176000884
		2.1 x 100 mm, 1.7 μm	3-Pack	176000885
<b>Shield RP<sub>18</sub></b>		1.0 x 50 mm, 1.7 μm	Column	186002851
		1.0 x 100 mm, 1.7 μm	Column	186002852
		2.1 x 50 mm, 1.7 μm	Column	186002853
		2.1 x 100 mm, 1.7 μm	Column	186002854
		1.0 x 50 mm, 1.7 μm	3-Pack	176000874
		1.0 x 100 mm, 1.7 μm	3-Pack	176000875
		2.1 x 50 mm, 1.7 μm	3-Pack	176000876
		2.1 x 100 mm, 1.7 μm	3-Pack	176000877
<b>Phenyl</b>		1.0 x 50 mm, 1.7 μm	Column	186002882
		1.0 x 100 mm, 1.7 μm	Column	186002883
		2.1 x 50 mm, 1.7 μm	Column	186002884
		2.1 x 100 mm, 1.7 μm	Column	186002885
		1.0 x 50 mm, 1.7 μm	3-Pack	176000905
		1.0 x 100 mm, 1.7 μm	3-Pack	176000906
		2.1 x 50 mm, 1.7 μm	3-Pack	176000907
		2.1 x 100 mm, 1.7 μm	3-Pack	176000908

# ACQUITY UPLC™ BEH Shield RP<sub>18</sub> Columns

Embedded polar group (Shield) RP columns contain stationary phases that combine the hydrophobicity of an alkyl ligand with the hydrophilicity of an embedded polar group. Features of embedded polar group columns include alternate selectivity to that of alkyl reversed-phase columns, excellent peak shape for bases and aqueous mobile phase compatibility. The unique selectivity of Shield RP phases, especially for polyphenolic compounds, has been attributed to the embedded polar groups acting as hydrogen bond acceptors. ACQUITY UPLC™ BEH Shield RP<sub>18</sub> columns are designed to provide selectivities that complement the ACQUITY UPLC™ BEH C<sub>18</sub> and C<sub>8</sub> phases.

ACQUITY UPLC™ BEH Shield RP<sub>18</sub> columns combine Waters patented Shield Technology™ with BEH Technology™ by incorporating an embedded carbamate group into the bonded phase ligand. Waters shield ligands consistently demonstrate increased retention of phenolic compounds versus straight chain alkyl columns. The alternate selectivity and excellent peak shape from the embedded polar group ligand, when combined with the wide pH range and ultra-efficiency of the 1.7 µm BEH particle, provide a necessary and powerful tool for UPLC™ method development.

## Fast Analgesics Separations




# ACQUITY UPLC™ BEH Phenyl Columns

Phenyl ligand-containing reversed-phase columns can provide complementary selectivities when compared to other straight-chain alkyl stationary phases, especially for analytes that contain aromatic rings. Traditional weaknesses of phenyl ligands include pH stability, reproducibility and peak shape. ACQUITY UPLC™ BEH Phenyl columns were intelligently designed to overcome these weaknesses and provide complementary selectivities, pH stability and excellent peak shape for all compounds.

ACQUITY UPLC™ BEH Phenyl columns utilize a trifunctional C<sub>6</sub> alkyl tether between the phenyl ring and the silyl functionality. This ligand, combined with the same proprietary endcapping processes as the ACQUITY UPLC™ BEH C<sub>18</sub> and C<sub>8</sub> columns, provides ultra-low column bleed, long column lifetimes and excellent peak shape. This unique combination of ligand and endcap on the 1.7 μm BEH particle creates a new dimension in selectivity and efficiency for challenging UPLC™ separations.

## Separations of Caffeic Acid Derivatives in Echinacea Purpurea



The active components in various Echinacea preparations can be divided into three major groups: caffeic acid derivatives, polysaccharides, and lipophilic components. Studies have shown that there is a variation in the quality of popular Echinacea-containing products. Thus, scientific monitoring of Echinacea's active ingredients is desired. Shown here are sub-four-minute separations of the primary caffeic acid derivatives—polyphenolic compounds that have been shown to possess antioxidant properties.



# Long UPLC™ Column Lifetimes

The innovation of ACQUITY UPLC™ BEH columns does not stop with the development of rugged and efficient stationary phases. Another major focus was the production of stable UPLC™ columns that provide the longest possible lifetimes under UPLC™ conditions. New column hardware was designed to minimize band broadening and ensure leak-free connections. An innovative frit design was implemented in order to keep the 1.7 µm particles in the column and out of

the detector or MS source. New proprietary column-packing stations and processes were developed to ensure a stable packed column bed and long, reproducible column lifetimes.


The result: Waters ACQUITY UPLC™ BEH columns provide column lifetimes under UPLC™ conditions which meet and/or exceed HPLC column lifetimes run under HPLC conditions.

## Long Column Lifetimes

Conditions  
Column: ACQUITY UPLC™ BEH Shield RP<sub>18</sub> 2.1 x 50 mm, 1.7 µm  
Mobile Phase A: 0.1% HCOOH in Water  
Mobile Phase B: 0.1% HCOOH in ACN  
Flow Rates: Low Pressure: 0.5 mL/min (~6,000 psi)  
High Pressure: 1.25 mL/min (~13,000 psi)  
Isocratic Mobile Phase Conditions: Low Pressure: 67% A; 33% B  
High Pressure: 62% A; 38% B  
Injection Volume: 5 µL  
Sample Diluent: MeOH  
Temperature: 65 °C  
Detection: UV @ 225 nm  
Instrument: ACQUITY UPLC™ with ACQUITY UPLC™ TUV

Compounds  
1. Ketorolac (208 µg/mL)  
2. Naproxen (13 µg/mL)

Isocratic Column Aging Procedure:  
1. Install and equilibrate column.  
2. Set flow rate to 0.5 mL/min (6,000 psi), inject and record data.  
3. Increase flow rate to 1.25 mL/min to achieve backpressure of 13,000 psi.  
4. Inject 200 to 400 times.  
5. Reduce flow rate to 0.5 mL/min (6,000 psi), inject and record data.  
6. Cycle and repeat steps 2 – 5.



**Over 4,000  
Injections at  
13,000 psi  
and 65 °C**

No demonstrated losses in efficiency, peak shapes or retention times after 4,000 injections at 13,000 psi and 65 °C. The column was removed after 4,000 injections—still fully functional.


Data courtesy of Dr. Ken Wehmeyer, The Proctor & Gamble Company

# Excellent Reproducibility

Waters continues to set the industry standard for column-to-column and batch-to-batch reproducibility. Beginning with the Symmetry® brand of HPLC columns in 1994 and continuing with the X Terra®, Atlantis® and SunFire™ brands, Waters columns provide consistent results, excellent peak shape and long column lifetimes. ACQUITY UPLC™ BEH

columns are manufactured in the same cGMP, ISO 9002-certified facilities that produce these industry-leading HPLC column brands. Method development scientists can be assured that the UPLC™ separation produced this year can be reproduced year after year.

## Excellent Batch-To Batch Reproducibility



# ACQUITY UPLC™ BEH Column Innovations

## Engineering

- New column hardware
- Low band broadening
- Innovative frit design
- Leak-free connections

## Bulk Synthesis

- Efficient 1.7  $\mu\text{m}$  particles (BEH technology™)
- Most technologically advanced porous particles ever created
- Combination of highest efficiencies, widest pH range and superior mechanical strength
- Stable column beds at UPLC™ pressures


## Software

- Paperless tracking of column history with new eCord™ technology

## Column Packing

- Columns must withstand UPLC™ operating pressures
- New proprietary packing methods
- New test instruments

ACQUITY UPLC™ columns feature eCord™ technology



Columns optimized for UPLC™ separations require innovative hardware and manufacturing processes that are not apparent when simply looking at the outside of the column. Like the holistically designed ACQUITY UPLC™ system, attention to every detail was critical for success. ACQUITY UPLC™ columns were designed to be an integral part of the low bandspread UPLC™ system. The ACQUITY UPLC™ column and system could not be developed by simply “designing down” or re-engineering an existing HPLC column or system. Typical HPLC system extra-column volumes and pressure limits would severely compromise the performance of UPLC™ columns. Like the BEH particle, the ACQUITY UPLC™ column is the most technologically advanced LC column ever created.

ing and sizing commercial quantities of 1.7  $\mu\text{m}$  BEH particles had to be developed. New packing stations and methods had to be invented and implemented since UPLC™ columns are packed and tested quite differently than HPLC columns. Additionally, UPLC™ instrumentation is necessary to test these columns, something that no other manufacturer possesses. Lastly, since UPLC™ technology was created with the future in mind, ACQUITY UPLC™ columns incorporate eCord™ technology — a step towards the paperless laboratory. Besides storing each column’s unique Certificate of Analysis, the eCord™ tracks column usage such as date of installation, number of injections, number of sample sets, maximum temperature and pressure and the date that the column was last used. All of this information travels with the column and is easily printable. The eCord™ is permanently attached to the column and the data cannot be erased.

A great deal of cutting edge manufacturing technology creates ACQUITY UPLC™ columns. Procedures for reproducibly produc-

## Sales Offices

Austria and European Export (Central South Eastern Europe, CIS and Middle East) (43) 1 8771807	Denmark (45) 46 598080	Japan (81) 3 3471 7191	Switzerland (41) 62 889 2030
Australia (61) 2 9933 1777	Finland (358) 9 506 4140	Korea (82) 2 820 2700	Taiwan (886) 2 2543 1898
Belgium (32) 2 7261000	France (33) 1 30 48 72 00	Mexico (5255) 5524 76 36	United Kingdom (44) 208 238 6100
Brazil (55) 11 5543 7788	Germany (49) 6196 40 06 00	The Netherlands 31 (0)76-50 87 200	All other countries: Waters Corporation U.S.A. 508 478 2000/800 252 4752
Canada 800 252 4752	Hong Kong (852) 2964 1800	Norway (47) 63 84 60 50	
China (8610) 84518918	Hungary (36) 1 350 5086	Poland (48) 22 833 4400	
CIS/Russia (7) 095 336 7000	India and India Subcontinent (91) 80 28371900	Puerto Rico (787) 747 8445	
Czech Republic (42) 02 6171 1384	Ireland (353) 1 448 1500	Singapore (65) 6278 7997	
	Italy (39) 02 27421 1	Spain (34) 93 600 93 00	Buy our products online at <a href="http://www.waters.com">www.waters.com</a> — click on shop
		Sweden (46) 8 555 11500	



The quality management system of Waters' manufacturing facilities in Taunton, Massachusetts and Wexford, Ireland complies with the International Standard ISO 9001:2000 Quality Management and Quality Assurance Standards. Waters' quality management system is periodically audited by the registering body to ensure compliance.

# Waters

©2005 Waters Corporation.  
Waters, ACQUITY UPLC, ACQUITY, UPLC, Ultra Performance LC, ACQUITY Ultra Performance LC, eCord, BEH Technology, XTerra, Symmetry, Atlantis, and SunFire are trademarks of Waters Corporation. All other trademarks are acknowledged.