

¹For complete experimental details, refer to full technical brief [720005252EN](#) at [waters.com](#)

INTRODUCTION

When transferring a method to a different HPLC/UHPLC instrument, the new instrumentation must provide the same separation and meet all of the system suitability requirements of the original method/instrument. For gradient separations, the impact of dwell volume can be dramatic and should be considered when transferring a method across chromatographic instruments.

LC CONDITIONS

Sample:	Abacavir and related substances
Sample diluent:	Water
Systems:	1. Agilent 1100 Series LC System 2. ACQUITY UPLC® H-Class System with ACQUITY UPLC PDA Detector and CH-A
Transfer:	Gradient delay was entered using SmartStart Technology
Wavelength:	245 nm
Column temp.:	36 °C
Column:	CORTECS® C ₁₈ , 2.7 µm, 4.6 x 75 mm
Mobile phase A:	0.1% Trifluoroacetic acid in water
Mobile phase B:	85% Methanol in water
Flow rate:	1.85 mL/min
Injection volume:	5 µL

Gradient:	Time (min)	%A	%B	%C
	Initial	95	5	Initial
	6.38	70	30	6.38
	10.37	10	90	10.37

Click on the [underlined blue text](#) for details on the products used in this application.

Waters Corporation

34 Maple Street, Milford, MA 01757 U.S.A.
T: 1 508 478 2000 F: 1 508 872 1990
[www.waters.com](#)


Waters, The Science of What's Possible, ACQUITY UPLC, and CORTECS are registered trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

CONCLUSIONS

An assay for the analysis of abacavir and related substances was successfully transferred from an Agilent 1100 Series LC System to an ACQUITY UPLC H-Class System. Using gradient SmartStart Technology in the instrument method, the differences in system volume were factored into the method. No changes to the gradient table were required.

Name	t _r Agilent 1100	t _r ACQUITY UPLC H-Class
1. Descyclopropyl abacavir	4.23	4.35
2. Abacavir	6.34	6.4
3. 1R, 4R Trans abacavir	6.64	6.7
4. O-(4-Chloro-2,5-diaminopyrimidinyl)-abacavir	6.93	6.99
5. O-t-Butyl-abacavir	7.70	7.73

Agilent 1100 Series LC System



ACQUITY UPLC H-Class System

