

Comprehensive Drug Analysis in Hair Samples: Extraction & UHPLC-MS/MS Quantification

Waters™

Authors: Jonathan Danaceau¹, Lisa J. Calton²
Affiliations: ¹Waters Corporation, Milford, MA, USA , ²Waters Corporation, Stamford Avenue, Wilmslow, Cheshire, U.K.

INTRODUCTION

The use of hair as a matrix for forensic toxicology continues to increase in popularity. As a specimen, hair offers several benefits. Sample collection is simple and easily supervised, and once collected, hair can be easily transported and stored at room temperature prior to analysis. Hair also provides an extended window of detection for drug exposure, enabling detection months and even years after use. This makes analysis of illicit drugs and pharmaceuticals in hair useful for situations where other samples may not be appropriate or available for analysis such as determining drug exposure in post-mortem toxicology, drug facilitated sexual assault (DFSA), or for other forensic testing in which long-term monitoring is desired.

The objective of this work was to develop, optimize and validate a method for the extraction and quantification of a comprehensive panel of drugs in hair to satisfy the confirmation cut-off values recommended by the society of Hair Testing (SoHT)¹. This was accomplished by optimizing pulverization, incubation conditions, the analytical workflow and solid phase extraction (SPE). The sample preparation protocol resulted in consistent recoveries and well-controlled matrix effects. The resulting method was linear, accurate and precise for all target compounds and easily met the SoHT cut-off criteria for all target analytes.

METHODS

Chemicals: Certified reference materials and internal standards were from Millipore Sigma and Cayman Chemical.

Reagents: Reference standards were used to prepare working multi-analyte calibrator and QC solutions in methanol. External quality control samples were acquired from Comedical (Italy) and consisted of authentic hair with drugs incorporated into the keratin matrix with assigned values.

Sample preparation: Samples were decontaminated by sequential washing with aqueous buffer and solvents. Bulk hair samples were pulverized using a Precellys Tissue Homogenizer and 2 mL CKMix Lysing Kits (Bertin Technologies, Montigny-le Bretonneux, FR) for 6 x 6400 rpm for 40 seconds each. Figure 1. shows a schematic of the pretreatment workflow. Samples were then extracted using Waters Oasis™ MCX 30 mg Plates. Figure 2 summarizes the SPE procedure

UHPLC-MS/MS analysis: A Waters ACQUITY™ UPLC™ I-Class (FTN) System was interfaced with a Xevo™ TQ Absolute Tandem Mass Spectrometer detector. Chromatography was performed using a Waters UPLC BEH™ C18 Column (1.7 μm, 2.1 x 100 mm) with a column temp. of 40°C. Mobile phase A (MPA) was 0.1% formic acid in water and mobile phase B (MPB) was 0.1% formic acid in LC-MS grade acetonitrile; flow rate was 0.6 mL/min. The LC gradient started at 2% B rising linearly to 67% MPB by 3.33 mins, then 90% B by 3.5 minutes before returning to 2% B from 3.6 to 4.0 mins. Injection volume was 2 μL. Mass spectrometer conditions were: source temperature 150°C, capillary voltage 1.0 kV desolvation gas (at 1000 L/h, 500°C) and cone gas (at 10 L/h). Data was processed with MassLynx™ Software and QUAN Review Application in the waters_connect™ Platform.




Figure 1. Graphical schematic of the pretreatment workflow



Figure 2. Solid Phase Extraction procedure

| Analyte | RT | R ² | Analyte | RT | R ² |
|----------------------|------|----------------|----------------------------|------|----------------|
| Morphine | 0.89 | 0.998 | Meperidine | 1.82 | 0.999 |
| Oxymorphone | 0.94 | 0.999 | Norbuprenorphine | 1.89 | 0.996 |
| Hydromorphone | 1.02 | 0.999 | Chloriazepoxide | 1.92 | 0.998 |
| Dihydrocodeine | 1.20 | 0.997 | Trazodone | 1.96 | 1.000 |
| Naloxone | 1.2 | 0.998 | Cocaethylene | 1.99 | 0.999 |
| Codeine | 1.23 | 0.999 | Phencyclidine | 2.06 | 0.999 |
| Noroxycodone | 1.29 | 0.997 | N-Pyrrolidino Etionitazine | 2.17 | 0.994 |
| Amphetamine | 1.31 | 0.995 | Fentanyl | 2.12 | 0.996 |
| Naltrexone | 1.26 | 0.994 | α-Hydroxymidazolam | 2.12 | 0.996 |
| Oxycodone | 1.27 | 0.999 | Midazolam | 2.15 | 0.998 |
| 6-acetyl morphine | 1.28 | 0.996 | Etionitazine | 2.21 | 0.996 |
| MDA | 1.29 | 0.995 | Flurazepam | 2.20 | 0.999 |
| Metadesnitazine | 1.36 | 0.998 | Buprenorphine | 2.24 | 0.996 |
| Hydrocodone | 1.33 | 0.997 | EDDP | 2.34 | 0.999 |
| O-desmethyl Tramadol | 1.32 | 0.999 | Methadone | 2.55 | 0.996 |
| Methamphetamine | 1.35 | 0.998 | α-Hydroxylprazolam | 2.55 | 0.998 |
| MDMA | 1.36 | 0.997 | α-Hydroxytriazolam | 2.55 | 0.999 |
| Phentermine | 1.41 | 0.996 | Nitrazepam | 2.56 | 0.999 |
| MDEA | 1.47 | 0.999 | Oxazepam | 2.63 | 0.999 |
| Ritalinic acid | 1.48 | 0.998 | Lorazepam | 2.70 | 0.998 |
| Norfentanyl | 1.53 | 0.999 | Clonazepam | 2.69 | 0.999 |
| Benzoyllecgonine | 1.52 | 1.000 | Alprazolam | 2.72 | 0.999 |
| 7-aminoclonazepam | 1.71 | 0.999 | 2-hydroxy-ethylflurazepam | 2.72 | 0.999 |
| Tramadol | 1.67 | 0.999 | Nordiazepam | 2.73 | 0.999 |
| N-desmethyiltramadol | 1.68 | 0.999 | Triazolam | 2.77 | 0.998 |
| Methylphenidate | 1.69 | 0.999 | Desalkylflurazepam | 2.82 | 0.999 |
| 7-aminoflunitrazepam | 1.71 | 0.999 | Flunitrazepam | 2.83 | 0.998 |
| Cocaine | 1.79 | 1.000 | Temazepam | 2.91 | 0.998 |
| Normeperidine | 1.80 | 0.999 | Diazepam | 3.10 | 0.999 |

Table 1. Target analytes, retention times and R² values from calibration curves

RESULTS

Recovery and Matrix Effects




Figure 3. Recovery and matrix effects for all compounds in the multi-drug panel.

Recoveries and matrix effects (ME) for all analytes are shown in Figure 3. Recoveries ranged from 6-79% with 49/58 > 40%. All %RSDs were under 17%. Amine stimulants and EDDP had lower recoveries than other analytes, but they were consistent and enabled accurate quantification even at the lowest concentrations. Matrix effects ranged up to –59% with only 6 compounds exceeding 40% ion suppression. They were also consistent with all S.D. values <10%. Phentermine and metadesnitazine were subject to some endogenous interferences which interfered with their quantification. All other compounds were well controlled.

CONCLUSIONS

- A method for the extraction and quantitative analysis of multiple classes of drugs from hair was developed and optimized to balance extraction efficiency and the stability needs of all compounds.
- The resulting method readily passed quantitative validation criteria for all compounds (excepting phentermine and metadesnitazine)
- The sensitivity requirements of SoHT were met for all compounds
- All analytes were positively identified in external quality control samples, with good quantitative correlation
- An efficient and timely workflow enables extraction within a 3-4 hour timeframe, enabling same-day batch analysis

References: 1—Cooper et al. 2012. Forensic Sci. Int 281, 20-24.

Quantitative Analysis

Calibration curves ranged from 0.01-1.0 ng/mg for most drugs, with a few exceptions. Norfentanyl and 6-acetyl morphine ranged from 0.002-0.2 ng/mg and fentanyl ranged from 0.001-0.1 ng/mg. Table 1 lists R² values for all analyte calibration curves. Table 2 lists the intra-batch statistics for all the analytes in the panel. With the exception of phentermine and metadesnitazine, all compounds met validation criteria for accuracy and precision, both for intra-batch results (shown) and inter-batch results. Table 2 also shows the cut-offs recommended by SoHT. All compounds satisfied the designated cut-off concentrations listed.

Table 2. Mean accuracy and precision for within batch QC results, together with lower limits of quantification (LLOQ) and SoHT cut-offs

| | Within Batch Statistics | | | | | | SoHT Cut-Off |
|----------------------|-------------------------|------|-------|------|-------|-----|--------------|
| | Low | Med | High | LLOQ | | | |
| | Mean | %CV | Mean | %CV | Mean | %CV | |
| Morphine | 101.3 | 4.2 | 105.1 | 5.8 | 93.4 | 5.6 | 0.01 |
| Oxymorphone | 98.6 | 2.7 | 100.4 | 1.1 | 96.4 | 3.3 | 0.01 |
| Hydromorphone | 98.2 | 2.6 | 101.5 | 1.2 | 97.3 | 4.2 | 0.01 |
| Dihydrocodeine | 98.1 | 2.5 | 101.5 | 1.8 | 97.2 | 4.8 | 0.01 |
| Naloxone | 91.6 | 5.1 | 100.8 | 4.5 | 94.2 | 4.4 | 0.01 |
| Codeine | 92.0 | 8.5 | 98.4 | 0.7 | 94.5 | 4.3 | 0.01 |
| Noroxycodone | 94.7 | 7.7 | 99.5 | 5.8 | 94.0 | 4.3 | 0.01 |
| Amphetamine | 117.6 | 10.9 | 101.5 | 2.5 | 95.8 | 2.1 | 0.01 |
| Naltrexone | 102.5 | 4.7 | 98.8 | 3.3 | 91.6 | 4.5 | 0.01 |
| Oxycodone | 96.4 | 2.5 | 103.6 | 1.9 | 97.4 | 3.6 | 0.01 |
| 6-acetyl morphine | 93.3 | 5.4 | 103.1 | 5.6 | 94.8 | 2.4 | 0.002 |
| MDA | 99.2 | 10.0 | 102.7 | 3.3 | 96.4 | 5.3 | 0.03 |
| Metadesnitazine | 95.9 | 4.1 | 97.6 | 2.3 | 86.9 | 3.4 | — |
| Hydrocodone | 98.2 | 4.1 | 103.5 | 1.6 | 94.5 | 3.6 | 0.01 |
| O-desmethyl Tramadol | 96.8 | 4.1 | 103.1 | 0.8 | 96.2 | 3.5 | 0.01 |
| Methamphetamine | 96.5 | 3.0 | 103.7 | 3.5 | 94.6 | 5.7 | 0.01 |
| MDMA | 96.2 | 6.1 | 102.0 | 5.1 | 97.8 | 6.9 | 0.01 |
| Phentermine | ND | ND | 44.0 | 85.5 | 92.3 | 7.2 | — |
| MDEA | 96.7 | 8.0 | 101.8 | 3.8 | 97.8 | 4.6 | 0.01 |
| Ritalinic acid | 82.7 | 9.0 | 94.7 | 7.1 | 94.3 | 5.1 | 0.01 |
| Norfentanyl | 91.5 | 4.3 | 101.4 | 2.8 | 105.8 | 5.1 | 0.002 |
| Benzoyllecgonine | 100.5 | 3.2 | 100.9 | 1.5 | 94.5 | 2.8 | 0.01 |
| 7-aminoclonazepam | 100.8 | 3.0 | 97.5 | 2.8 | 92.0 | 2.7 | 0.01 |
| Tramadol | 97.8 | 1.5 | 103.2 | 2.4 | 98.3 | 2.6 | 0.01 |
| N-desmethyiltramadol | 96.5 | 5.5 | 105.2 | 3.4 | 96.0 | 3.0 | 0.01 |
| Methylphenidate | 99.3 | 1.7 | 101.0 | 1.6 | 95.5 | 2.7 | 0.01 |
| 7-aminoflunitrazepam | 96.4 | 4.2 | 100.9 | 1.6 | 95.3 | 3.1 | 0.01 |
| Cocaine | 99.7 | 3.2 | 101.3 | 2.3 | 96.0 | 4.2 | 0.01 |
| Normeperidine | 97.2 | 4.1 | 102.4 | 1.6 | 96.7 | 3.2 | 0.01 |


| | Within Batch Statistics | | | | | | SoHT Cut-Off |
|----------------------------|-------------------------|------|-------|------|-------|-----|--------------|
| | Low | Med | High | LLOQ | | | |
| | Mean | %CV | Mean | %CV | Mean | %CV | |
| Meperidine | 97.4 | 2.7 | 101.7 | 1.1 | 95.2 | 3.9 | 0.01 |
| Norbuprenorphine | 102.8 | 7.7 | 97.3 | 4.3 | 97.6 | 3.1 | 0.01 |
| Chloriazepoxide | 99.9 | 4.1 | 113.9 | 0.6 | 102.3 | 2.4 | 0.01 |
| Trazodone | 100.2 | 2.5 | 105.2 | 2.2 | 99.3 | 3.6 | 0.01 |
| Cocaethylene | 99.0 | 2.3 | 101.4 | 1.9 | 96.5 | 3.3 | 0.01 |
| Phencyclidine | 98.0 | 2.1 | 101.5 | 1.0 | 95.7 | 3.7 | 0.01 |
| N-Pyrrolidino Etionitazine | 98.6 | 3.7 | 100.1 | 2.9 | 90.5 | 3.2 | 0.01 |
| Fentanyl | 99.5 | 2.8 | 101.5 | 2.3 | 94.2 | 2.6 | 0.001 |
| α-Hydroxymidazolam | 108.7 | 7.9 | 105.8 | 3.2 | 93.6 | 3.3 | 0.03 |
| Midazolam | 94.6 | 4.3 | 102.0 | 1.7 | 96.6 | 1.9 | 0.01 |
| Etionitazine | 96.4 | 5.8 | 100.3 | 4.1 | 91.8 | 4.2 | 0.01 |
| Flurazepam | 98.6 | 3.6 | 103.5 | 2.5 | 99.8 | 2.1 | 0.01 |
| Buprenorphine | 96.6 | 6.5 | 97.4 | 4.8 | 93.4 | 4.3 | 0.01 |
| EDDP | 100.7 | 4.0 | 102.2 | 1.9 | 96.5 | 3.5 | 0.01 |
| Methadone | 99.3 | 3.9 | 103.7 | 2.1 | 97.3 | 2.9 | 0.01 |
| α-Hydroxylprazolam | 92.9 | 9.1 | 101.7 | 2.6 | 97.1 | 2.1 | 0.01 |
| α-Hydroxytriazolam | 94.8 | 9.9 | 101.8 | 2.2 | 99.3 | 3.1 | 0.01 |
| Nitrazepam | 96.5 | 4.9 | 101.3 | 2.1 | 100.2 | 3.8 | 0.01 |
| Oxazepam | 103.7 | 11.2 | 107.6 | 6.0 | 94.1 | 5.0 | 0.02 |
| Lorazepam | 97.8 | 10.1 | 103.6 | 2.1 | 95.9 | 4.9 | 0.02 |
| Clonazepam | 104.6 | 5.4 | 106.0 | 4.4 | 93.3 | 3.3 | 0.03 |
| Alprazolam | 97.8 | 5.9 | 100.3 | 1.2 | 95.7 | 3.5 | 0.01 |
| 2-hydroxyethylflurazepam | 101.8 | 8.5 | 102.6 | 7.3 | 96.7 | 4.0 | 0.01 |
| Nordiazepam | 98.7 | 2.7 | 101.6 | 2.8 | 96.5 | 4.0 | 0.01 |
| Triazolam | 102.5 | 3.4 | 103.0 | 3.9 | 94.1 | 2.2 | 0.01 |
| Desalkylflurazepam | 94.7 | 2.5 | 103.0 | 3.0 | 96.6 | 2.8 | 0.01 |
| Flunitrazepam | 93.9 | 1.7 | 101.6 | 2.1 | 95.2 | 3.1 | 0.01 |
| Temazepam | 96.0 | 4.0 | 101.6 | 2.3 | 97.2 | 0.9 | 0.01 |
| Diazepam | 96.0 | 2.1 | 98.5 | 2.3 | 96.2 | 2.6 | 0.01 |

External Assessment

Table 3. Observed mean values (N=5), nominal concentrations and reference ranges of external control samples. All included compounds were positively identified. Fentanyl and 6-acetyl morphine had values beyond the range of this assay. For the remaining compounds, 18/22 (82%) were within the control limits assigned by the manufacturer and all had %CVs <10%.

| | EQC Result | | | | | |
|-------------------|--------------|------|---------------|-------------|-------------|------------|
| | Mean (ng/mg) | %CV | Nominal Conc. | Lower Limit | Upper Limit | Acceptable |
| Morphine | 0.473 | 3.31 | 0.48 | 0.310 | 0.650 | y |
| Dihydrocodeine | 0.583 | 4.35 | 0.47 | 0.310 | 0.630 | y |
| Codeine | 0.253 | 3.41 | 0.27 | 0.180 | 0.360 | y |
| Amphetamine | 0.464 | 4.07 | 0.51 | 0.330 | 0.690 | y |
| Oxycodone | 0.532 | 5.51 | 0.47 | 0.300 | 0.630 | y |
| 6-acetyl morphine | 0.841 | 2.68 | 0.57 | 0.370 | 0.770 | n* |
| MDA | 0.394 | 6.58 | 0.42 | 0.270 | 0.570 | y |
| Methamphetamine | 0.516 | 2.57 | 0.60 | 0.390 | 0.810 | y |
| MDMA | 0.661 | 9.37 | 0.74 | 0.480 | 1.000 | y |
| MDEA | 0.653 | 7.35 | 0.63 | 0.410 | 0.850 | y |
| Benzoyllecgonine | 0.500 | 4.07 | 0.36 | 0.230 | 0.490 | n |
| Tramadol | 0.434 | 2.18 | 0.52 | 0.340 | 0.700 | y |
| Cocaine | 0.862 | 3.61 | 0.90 | 0.590 | 1.220 | y |
| Norbuprenorphine | 0.091 | 9.17 | 0.06 | 0.039 | 0.081 | n |
| Cocaethylene | 0.876 | 3.80 | 0.58 | 0.380 | 0.780 | n |
| Fentanyl | 0.224 | 3.17 | 0.14 | 0.090 | 0.190 | n* |
| Buprenorphine | 0.065 | 4.78 | 0.05 | 0.033 | 0.068 | y |
| EDDP | 0.546 | 4.17 | 0.41 | 0.270 | 0.550 | y |
| Methadone | 0.866 | 8.91 | 0.58 | 0.380 | 0.780 | n |
| Lorazepam | 0.198 | 3.18 | 0.27 | 0.180 | 0.360 | y |
| Alprazolam | 0.255 | 4.23 | 0.23 | 0.150 | 0.310 | y |
| Nordiazepam | 0.230 | 2.83 | 0.33 | 0.210 | 0.450 | y |
| Temazepam | 0.177 | 3.27 | 0.25 | 0.160 | 0.340 | y |
| Diazepam | 0.224 | 1.89 | 0.27 | 0.180 | 0.360 | y |

Figure 4. Correlation between the nominal values and our observed values. Each point represents a different analyte, the overall slope of 1.05 indicates excellent overall agreement with the assigned values.



TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

For Forensic Toxicology use only

Waters, ACQUITY, UPLC, Xevo, BEH, Oasis, waters_connect and MassLynx are all trademarks of Waters Technologies Corporation

©2025 Waters Corporation
720009134EN