

# Identification of Novel Per-and-Polyfluoroalkyl Substance (PFAS) Isoform in Textile Using a Multi-Reflecting Time-of-Flight Mass Spectrometer Technology

Hania Hollins<sup>1</sup>, Sara Beverley<sup>1</sup>, Martin Palmer<sup>1</sup>, Jayne Kirk<sup>1</sup>, Richard Lock<sup>1</sup>, Sean Wu<sup>2</sup>

<sup>1</sup>Waters Corporation Wilmslow, Cheshire, UK, <sup>2</sup>Waters Corporation, Milford, MA, USA

# Waters™

720009415EN

## 1 Introduction

As awareness of the health risks associated with per- and polyfluoroalkyl substances (PFAS) grows, there is an increasing demand to determine total PFAS content to help reduce human exposure. Traditional targeted methods provide only a narrow insight of the PFAS content, making non-targeted screening with high-resolution mass spectrometry an effective approach to bridge quantitative analyses and discovery workflows.

In this study, high-resolution mass spectrometry (HRMS) combined with a non-targeted screening (NTS) workflow was applied to characterize PFAS profiles in textiles. Textile reference material (LCG) was analyzed using data dependent acquisition mode (DDA) and processed with waters\_connect™ Software Platform. Different PFAS isoforms and halogenated contaminants were identified with high confidence levels based on accurate mass measurement of the precursor and product ions.

## 2 Experimental




Figure 1. Illustration of PFAS extraction from textile reference material (TRM). Briefly, PFAS were extracted from milled textile material derived from outdoor clothing using 3x0.5 g of TRM according to EN 17681-1:2025<sup>1</sup>.

## 3 LC-MS conditions

### Chromatographic conditions:

- **UPLC™ system:** Waters ACQUITY Premier modified with PFAS solution installation kit (p/n: 176004548).
- **Column:** ACQUITY Premier BEH™ C18, 1.7 μm, 2.1 x 50 mm, 90 Å (p/n:186009452)
- **Column temp.:** 40°C
- **Sample temp.:** 10°C
- **Injection volume:** 5 μL

- **Mobile phase A:** 95:5 water:methanol + 2 mM ammonium acetate
- **Mobile phase B:** 100% methanol + 2 mM ammonium acetate

### Chromatographic gradient

Time (min)	Flow rate (mL/min)	%A	%B	Curve
0	0.5	100	0	Initial
0.5	0.5	80	20	6
3.0	0.5	55	45	6
6.5	0.5	20	80	6
7.0	0.5	5	95	6
8.5	0.5	5	95	6
9.0	0.5	100	0	1
11.0	0.5	100	0	1

### MS acquisition parameters

	Parameters	Value
Source	Capillary voltage	0.5 kV
	Cone voltage	10 V
	Source temp.	100°C
	Desolvation temp.	250°C
	Cone gas	100 L/h
Transmission Settings	StepWave RF	100 V
	Body gradient	5 V
Acquisition mode	Data dependent (DDA)	Top 5
Collision energy	Low mass ramp	20 - 40 V
	High mass ramp	35 - 75 V

## 4 Rational

A merged database composed from CompTox<sup>2</sup> (14735 compounds) and EPA<sup>3</sup> (4969 compounds) was created. Common compounds were removed. To all unique 19610 compounds, the number of "F" atoms and the mass defect were calculated.




Figure 2 left is a plot of the Mass defect as function of the number of "F" atoms. Right is the frequency of mass defect in the created library.

## 5 Results

Using Mass defect directed DDA reduced the number of MS/MS spectra by ~25 %




Figure 3. Number of MS/MS spectra triggered without (pink) and with (blue) mass defect in three consecutive injections of a standard mix at 5 ng/mL. The number of DDA triggered was reduced by ~25% in each of the individual injections.

## 6 Results




Figure 4. The frequency of triggered MS/MS using a DDA method without (pink) and with Mass Defect Derived DDA.

The frequency of DDA triggered on PFAS-like compounds was improved. Figure 4 demonstrated that without MD derived DDA (pink), compounds with Mass defects between 0.15-0.85 were selected from the survey scan and fragmented

➔ MD derived DDA improved the number of PFAS-like compounds selected for MS/MS

## 7 Identification of PFAS in Textile Reference Material




Figure 5. A: Extracted ion chromatogram of N-ethyl perfluorooctane sulfonamidoethanol (N-EtFOSE) in textile extract. B (Top) Survey scan and B (bottom) MS/MS spectra of N-EtFOSE. N-EtFOSE is identified at m/z 630.02551 with mass measurement accuracy 1.1 ppm.

## Conclusions

- Mass defect driven data dependent acquisition reduced the number of unambiguous MS/MS spectra.
- Improved the number of identified PFAS and PFAS-like compound.
- DDA driven by mass defect can be used for other types of persistent contaminants such as polychlorinated and polybrominated biphenyls.

## References

1. EN 17681-1:2025 Textiles and Textile products – Per and polyfluoroalkyl substances (PFAS) – Part 1: Analysis of an alkaline extract using liquid chromatography and tandem mass spectrometry
2. U.S. Environmental Protection Agency. CompTox Chemicals Dashboard. <<https://comptox.epa.gov/dashboard/chemical-lists/pfasmaster>> accessed November 2023.
3. Suspect List of Possible Per- and Polyfluoroalkyl Substances (PFAS). <<https://data.nist.gov/od/id/nds2-2387>> Accessed March 2023

ACQUITY, Xevo, UNIFI, UPLC, BEH, waters\_connect are trademarks of Waters Technologies Corporation.