

# Improving Detection and Selectivity of PFAS Molecules Using Cyclic Ion Mobility Wideband Enhancement

Michael McCullagh<sup>1</sup>, Yuriy Pyatkivskyy<sup>2</sup>, Brian McCullough<sup>1</sup>, Darren Hewitt<sup>1</sup>, Zoltan Laszlo<sup>1</sup> and Chris Hughes<sup>1</sup>  
<sup>1</sup>Waters Corporation, Stamford Avenue, Altrincham Road, Wilmslow, SK9 4AX. UK. <sup>2</sup>Waters Corporation, Milford, Massachusetts, USA.

# Waters™

## Introduction

Polyfluorinated alkyl substances (PFAS) are emerging environmental contaminants that have been linked to increased risks of cancer and other diseases in humans. Monitoring PFAS concentrations in human biofluids is essential for understanding exposure levels and pathways. However, accurately identifying PFAS at extremely low concentrations (ng/L) remains challenging. To address this, a Cyclic™ IMS P20 Instrument equipped with Wideband Enhancement (WBE) can be programmed to target characteristic PFAS drift time trendlines and synchronize them with the time-of-flight (ToF) pusher (see Figures 1 and 2).<sup>1</sup> This acquisition mode improves the ToF duty cycle and lowers detection limits. Analysis using Cyclic ion mobility with WBE provides enhanced selectivity, sensitivity, and collision cross section (CCS) measurements, which serve as complementary descriptors for PFAS identification.<sup>2,3</sup>

## Experimental

LC-Cyclic-IMS-MS analyses of a dilution series of Native PFAS: PFAC30PAR (Wellington laboratories) were performed using a quadrupole-Cyclic-IMS-time-of-flight mass spectrometer (Cyclic IMS resolution R ~65–145). Reversed-phase LC separation employed mobile phases of: A (95:5 water: methanol, with 2 mM ammonium acetate) and B (methanol with 2 mM ammonium acetate) using a 22 min gradient at 0.3 mL/min. An ACQUITY™ UPLC™ BEH™ C18 column (100 mm × 2.1 mm, 1.8 μm) at 35 °C with a 10 μL injection volume was used. The modified ACQUITY i-Class System included a PFAS-free conversion kit and Atlantis™ Premier BEH C18 AX Isolator Column, 2.1 x 50mm, 5 μm.<sup>4</sup>

## Results

The environmental impact of PFAS exposure is influenced by isomeric distribution, as both linear (L-PFAS) and branched (br-PFAS) forms have been associated with specific health risks, including increased infant burden during pregnancy. Enhanced sensitivity in the WBE analytical method for precursor ions, enables the detection of PFAS at lower concentrations compared to conventional HDMSE acquisition. The WBE acquisition strategy allows targeting of mass/drift time trendlines, with potential to reduce false detections caused by isobaric biological interferences. PFAS molecules exhibit distinct ion mobility drift time trendlines due to their unique CCS-to-mass ratios, a result of their fluorinated structures. This property causes PFAS to appear “denser” than hydrocarbon-based molecules, shifting their drift time trendline downward relative to molecules of similar *m/z*.<sup>5</sup> For a dilution series of PFAC30PAR we compared LC-Cyclic-IMS-MS (HDMSE) and WBE LC-Cyclic-IMS-MS, LC-Cyclic-IMS-MSMS, LC-Cyclic-IMS-DIA-MS responses for PFAS. The WBE mode was calibrated to optimize drift time alignment (See Figure 3) for PFAS constituents of PFAC30PAR, the improvement in sensitivity resulting from enhanced selectivity and ToF duty cycle is shown in Figure 4 and Figure 5. An RMS S/N 48 for 0.4 ng/L for the PFOS homologous series have been achieved as well as improved detection limits for branched PFOS byproducts, P5MHpS and P6MHpS, compared to using conventional HDMSE (see Figure 6).




Figure 1. Cyclic IMS P20 System.




Figure 2. Schematic representation of WBE Cyclic ion mobility selectivity combined with ToF pusher synchronization.




Figure 3. *m/z* versus drift time plots demonstrating WBE: (a) WBE off, showing selection of singly charged PFAS characteristic drift times for calibration; (b) WBE on, illustrating improved ion mobility mass spectrometry selectivity for PFAS detection.




Figure 4. Comparison of WBE on/off conditions: base peak intensity (BPI) for PFAC30PAR (4 ng/L), and PFOS (*m/z* 499) assessed using extracted mass chromatograms (EMC) and RMS signal-to-noise (S/N).

## Conclusions

- Integration of Cyclic ion mobility and WBE improves analytical selectivity and sensitivity for PFAS.
- Application of WBE with Cyclic ion mobility HDMSMS enabled PFOS detection at 0.4 ng/L.
- LC-Cyclic-IMS-MS with WBE provides a flexible platform for both targeted and non-targeted analysis of known and emerging PFAS in complex matrices.

### References

1. Hughes C. Enhanced Glycopeptide Identification Using the SYNAPT XS Q-ToF with Ion Mobility Enabled Wideband Enhancement. 2020. Waters Application Note 720006858.
2. Enhanced Identification Confidence and Specificity for PFAS Analysis Using Cyclic Ion Mobility Mass Spectrometry Collision Cross Sections. 2024. Waters Application Note 720008536.
3. McCullagh, M., Kass, I., Lioupi, A., Theodoridis, G., Plumb, R., Dowd, S., & Adams, S. (2024, June). The utility of cyclic ion mobility to improve selectivity and analysis efficiency of environmental PFAS contamination and exposure [Poster presentation]. ASMS Conference on Mass Spectrometry and Allied Topics, Anaheim, CA, United States.
4. PFAS Analysis Kit for ACQUITY UPLC Systems User Guide 720006689 < <https://www.waters.com/webassets/cms/support/docs/720006689en.pdf>.
5. Dodds, J. N., Hopkins, Z. R., Knappe, D. R. U., & Baker, E. S. (2020). Rapid characterization of per- and polyfluoroalkyl substances (PFAS) by ion mobility spectrometry-mass spectrometry (IMS-MS). *Analytical Chemistry*, 92(6), 4427–4435. <https://doi.org/10.1021/acs.analchem.9b05036>




Figure 5. Comparison of br-PFOS and L-PFOS WBE responses on and off at 4 ng/L, demonstrating a ~10 × increase in signal intensity and improved branched isomer chromatographic peak fidelity.




Figure 6. Detection limits and RMS S/N obtained for PFOS using HDMSMS with WBE.

“For Research Use Only. Not for use in diagnostic procedures.”

Waters, SELECT SERIES, Cyclic, Atlantis, ACQUITY, BEH, UPLC are trademarks of Waters Corporation or its affiliates.”

**Conflict of Interest Disclosure:** Yuriy Pyatkivskyy, Brian McCullough, Darren Hewitt, Zoltan Laszlo, Chris Hughes and Michael McCullagh are employees of Waters Corporation and are presenting the poster on behalf of the company.