

ANALYSIS OF FREE INOSITOL STEREOISOMERS IN DIETARY SUPPLEMENTS AND FOODS BY HYDROPHILIC LIQUID CHROMATOGRAPHY

Jinchuan Yang, Stephanie Harden, and Paul Rainville (contact: Jinchuan_yang@waters.com)
Waters Corporation, 34 Maple St., Milford MA 01757

INTRODUCTION

Inositols are cyclic carbohydrates with nine stereoisomers (Figure 1), of which myo- and D-chiro-inositol are the most abundant in humans, while others (e.g., scyllo-, neo-, and epi-) are also present at lower levels. Inositols are associated with important health benefits, including the management of insulin resistance, metabolic syndrome, PCOS, and gestational diabetes. Adequate intake is important for babies whose digestive and metabolic systems are still developing, and for individuals with inositol deficiency. With increasing incorporation of inositol into supplement and food products, there is growing demand for analytical methods capable of profiling inositol stereoisomers across various food products.

HPLC and GC are commonly used for inositol analysis. While GC provides high resolution, it requires derivatization and long analysis times. HPLC enables faster analysis without derivatization, but standard methods are typically limited to individual isomers. There is no standard method for simultaneous determination of inositol stereoisomers in supplement and food products.

OBJECTIVES

This study aimed to develop HILIC-MS/MS methods for the simultaneous separation and quantification of free inositol stereoisomers in dietary supplements and food products.






Figure 1. Structures of inositol stereoisomers

EXPERIMENTAL

Sample Prep.:



LC-MS/(MS) Methods:

For dietary supplements		For infant formula, milk, and grain				
Software: Empower™ 3 CDS	ACQUITY QDa II Mass Detector setting:	Software: MassLynx™	XevoTQ-S cronos MS setting:			
Ionization mode: ESI-	Capillary voltage: 0.5 kV	Ionization mode: ESI-	Capillary voltage: 2.50 kV			
Probe Temp.: 600 °C	Cone voltage: 7.0 V	Source temp.: 150 °C	Desolvation temp.: 350 °C			
SIR [M]: 179 Da, 185 Da (ISTD)		Cone gas flow: 40 L/Hr	Desolvation gas flow: 650 L/Hr			
LC conditions:	ACQUITY UPLC™ BEH™ Amide Column, or ACQUITY Premier BEH Amide VanGuard™ FIT Column (1.7 µm, 2.1 × 150 mm), Col. temp.: 25 °C	MRM transitions and detection parameters:				
MP A: ACN:water (90:10 v/v) with 0.01% NH ₃		Analyte	Dwell (secs)	MRM Transitions	Cone Voltage (V)	Collision Energy (eV)
MP B: ACN:water (50:50 v/v) with 0.01% NH ₃ , 20 mM NH ₄ HCO ₃		inositol	0.1	179.0 > 86.9 (Quan)	20	12
Gradient elution table:		inositol-d ₆	0.1	179.0 > 98.9 (Qual)	20	14
Time	Flow Rate (mL/min)	A%	B%	Curve		
Initial	0.35	91%	9%	Initial		
5.7	0.35	91%	9%	6		
13.1	0.35	79%	21%	6		
17.4	0.35	79%	21%	6		
18.0	0.35	91%	9%	6		
23.0	0.35	91%	9%	6		

RESULTS AND DISCUSSION

Chromatographic separation optimization

- Column temperature
- Gradient elution conditions
- Injection volume
- QDa II optimization (capillary voltage and probe temperature)




Figure 2. SIR chromatograms of inositol standards (allo-, muco-, epi-, chiro-, neo-, myo-, and scyllo-inositol) and an ISTD (myo-inositol-d₆). Resolution for allo-/muco-, chiro-/neo-, and myo-/scyllo- are 1.6, 2.3, and 2.0, respectively.

Selectivity

- no interference from common monosaccharides (glucose, fructose, galactose, and mannose)

Ruggedness

- No statistically significant difference across 3 columns (different batches and format)

Table 1. Comparison of results across three columns (different batches and formats).

	myo-Inositol results (µg/mL)			D-chiro-Inositol results (µg/mL)		
	Col. #1	Col. #2	Col. #3	Col. #1	Col. #2	Col. #3
Measurement 1	0.586	0.614	0.623	0.568	0.627	0.624
Measurement 2	0.616	0.608	0.634	0.629	0.596	0.600
Measurement 3	0.640	0.600	0.633	0.593	0.606	0.645
Average	0.614	0.607	0.630	0.597	0.610	0.623
SD	0.027	0.007	0.006	0.031	0.016	0.023
One-way ANOVA P-value	0.30			0.45		
Comments	The difference among the averages are not statistically significant at α=0.05			The difference among the averages are not statistically significant at α=0.05		

Linearity and sensitivity

Table 2A. Calibration results and limit of quantitation (LOQ) for dietary supplements.

Analyte	Equation	R ²	Conc. Range (ppm)	LOQ (ppm)
Scyllo-Inositol	Y = -3.73 x 10 ⁻³ X ² + X - 5.32x10 ⁻³	0.9993	0.2 - 15	0.2
Myo-Inositol	Y = -9.50 x 10 ⁻⁴ X ² + 0.616 X - 8.02 x 10 ⁻²	0.9995	0.3 - 15	0.3
Neo-Inositol	Y = -1.30 x 10 ⁻³ X ² + 1.16 X - 3.46 x 10 ⁻¹	0.9993	0.9 - 15	0.9
Chiro-Inositol	Y = -6.27 x 10 ⁻³ X ² + 1.47 X - 5.00 x 10 ⁻³	0.9981	0.1 - 15	0.1
Epi-Inositol	Y = -9.41 x 10 ⁻⁴ X ² + 0.891 X - 1.36	0.9978	0.9 - 15	0.9
Muco-Inositol	Y = -6.02 x 10 ⁻³ X ² + 1.94 X - 7.25 x 10 ⁻¹	0.9990	0.9 - 15	0.9
Allo-Inositol	Y = -1.81 x 10 ⁻² X ² + 2.19 X - 1.11 x 10 ⁻¹	0.9993	0.09 - 15	0.09

Table 2B. Calibration results and LOQ for foods.

Analyte	Equation	R ²	Conc. Range (µg/mL)	LOQ* (µg/mL)	LOQ in matrix** (µg/mL)
Myo-Inositol	Y = 4.62 x 10 ⁻³ X ² + 1.66 X - 2.59 x 10 ⁻²	0.9998	0.15 - 15	0.09	0.3
D-Chiro-Inositol	Y = -4.54 x 10 ⁻² X ² + 4.82 X + 1.13 x 10 ⁻²	0.9991	0.05 - 15	0.05	0.2

Note: * LOQ estimated in solvent; ** LOQ estimated in whole milk.

Accuracy and precision

Table 3A. Summary of accuracy and precision results in spiking experiments for dietary supplements.

Dietary Supplement	myo-Inositol					
	Mean Conc. (µg/mL)	RSD (n=3)	Spiked Conc. (µg/mL)	Spiked level*	RSD (n=3)	Recovery
DS-1	2.71	0.5%	4.56	168%	0.3%	105.0%
DS-2	3.72	0.4%	7.60	204%	0.1%	101.1%
DS-3	2.40	0.6%	7.60	317%	0.1%	103.4%
DS-4	2.35	0.3%	1.52	65%	0.2%	106.0%

Dietary Supplement	D-chiro-Inositol					
	Mean Conc. (µg/mL)	RSD (n=3)	Spiked Conc. (µg/mL)	Spiked level*	RSD (n=3)	Recovery
DS-1	0.00	-	4.479	-	0.7%	99.9%
DS-2	0.10	1.1%	7.465	-	0.7%	97.9%
DS-3	0.06*	7.4%	7.465	-	0.8%	97.9%
DS-4	0.03*	6.9%	1.493	-	0.7%	99.3%

Note: * spiked level relative to the native level; * close to LOQ.

Table 3B. Summary of accuracy and precision results in spiking experiments for foods.

	myo-Inositol			D-chiro-Inositol		
	Spike Conc. (µg/mL)	RSD (n=3)	Recovery	Spike Conc. (µg/mL)	RSD (n=3)	Recovery
Solvent	0.61	4.4%	101%	0.61	5.2%	98%
Soybean milk	2.43	2.0%	96%	2.39	6.7%	102%
	0.61	3.6%	107%	0.61	2.5%	103%
Almond milk	1.22	2.2%	107%	1.19	1.0%	104%
Oat milk	1.22	3.2%	110%	1.19	7.5%	100%
Whole milk	1.22	2.2%	108%	1.19	2.5%	96%
	2.43	2.2%	102%	2.39	4.7%	107%
Cornmeal	0.61	3.2%	103%	0.60	3.5%	99%
Infant Formula	0.61	0.6%	89%	0.60	3.0%	97%
	0.61	1.7%	117%	0.61	5.0%	109%

Sample analysis results

Table 4A. Summary of dietary supplement analysis results.

Dietary Supplement	myo-Inositol			D-chiro-Inositol		
	Label value (mg/serving)	Mean (mg/serving)	Rel. to label	Label value (mg/serving)	Mean (mg/serving)	Rel. to label
DS-1	2000	1965	98.2%	50	0.0	0%
DS-2	2000	1893	94.6%	50	48.6	97%
DS-3	2000	2158	107.9%	50	58.2	116%
DS-4	2000	2049	102.4%	50	24.1	48%

Table 4B. Summary of food analysis results.

	myo-Inositol		D-chiro-Inositol	
	Mean Conc. (mg/100 g) ^a	RSD (n=3)	Mean Conc. (mg/100 g)	RSD (n=3)
Soybean milk	4.16	3.0%	3.10	2.3%
	4.08	3.8%	3.22	3.7%
	4.22	1.4%	3.29	1.8%
	Average (interday, n=9)		Average (interday, n=9)	
	4.16	2.9%	3.21	3.5%
Almond milk	2.46	0.6%	N.D.	
Oat milk	0.64	7.2%	N.D.	
Whole milk	3.13	1.7%	N.D.	
	2.88	3.2%		
	Average (interday, n=6)			
	3.01	5.1%		
Cornmeal	34.72	2.2%	N.D.	
Infant formula	83.45	1.1%	N.D.	
	77.18	0.2%		
	Average (interday, n=6)			
	80.32	4.3%		

Note: ^a mg/100mL for cornmeal and infant formula.

CONCLUSION

- We have developed fit-for-purpose solutions for analysis of inositol stereoisomers in dietary supplements and foods.
- Both analytical solutions demonstrated excellent linearity, sensitivity, specificity, accuracy, precision and ruggedness.
- These methods were successfully applied to commercial products, demonstrating their suitability for routine quantitation of free inositol stereoisomers across a variety of dietary and food matrices.

Reference

- 1) Yang, J., Harden, S., and Rainville, P. HILIC-MS/MS Analysis of Free Inositol Stereoisomers in Foods, Application Note 720009200, Waters Corporation, 2026.
- 2) Yang, J., Harden, S., and Rainville, P. Analysis of Free Inositol Stereoisomers in Dietary Supplements by Hydrophilic Liquid Chromatography using the Arc™ Premier System and ACQUITY QDa II Mass Detector, Application Note 720009186, Waters Corporation, 2025.

Please scan the QR Codes on the right to download the application notes.

