6.1052.4X0 Metrosep C Supp 1 - XXX/4.0

6.1052.410 Metrosep C Supp 1 - 100/4.0

6.1052.420 Metrosep C Supp 1 - 150/4.0

6.1052.430 Metrosep C Supp 1 - 250/4.0

DE

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EN

Column material

Polyvinyl alcohol with carboxyl groups, particle
size 5 µm

Dimensions

6.1052.410: 100 x 4.0 mm

6.1052.420: 150 x 4.0 mm

6.1052.430: 250 x 4.0 mm

pH range

1 to 12 (eluent)

1 to 12 (sample)

Temperature range

20 to 40 °C

Recommended standard temperature: 40 °C (for
all lengths)

Maximum pressure

15 MPa (150 bar)

Maximum flow rate

1.5 mL/min (for all lengths)

Recommended standard flow rate: 1.0 mL/min
(for all lengths)

Application

Determination of monovalent and bivalent cati-
ons with suppression and conductivity detection.

Eluent

- Standard eluent:
5 mmol/L nitric acid
50 µg/L rubidium (from RbNO₃ salt)
- Permitted eluent additives and sample addi-
tives:
0 to 30% acetone, 0 to 50% acetonitrile, no
methanol

Preparation

- Rinse the column with eluent for 2 h.
- To avoid high backpressure, we recommend
rinsing the column after installation at a low
flow rate (0.40 mL/min) for approx. 20 min.

Guard column

Metrosep C Supp 1 Guard/4.0 (6.1052.500)

Metrosep C Supp 1 S-Guard/4.0 (6.1052.510)

Storage

Store the column in ultrapure water at 4 to 8 °C.

Do not store the column below 0 °C.

NOTE

To avoid performance losses during storage,
we recommend keeping the column in a
refrigerator.

Regeneration

NOTE

Ensure that the maximum pressure is never exceeded during regeneration.

If the pressure becomes too high, reduce the flow rate.

1. Disconnect the column outlet from the attached function units such as suppressor or detector. Instead, collect the flow of liquid in a beaker.
2. Rinse the column with ultrapure water before and after regeneration.

Depending on the type of contamination, proceed in accordance with the following instructions:

Organic contamination

Regenerate the column in the direction opposite to the flow at a flow rate of 1 mL/min with the following solutions in succession:

1. 1 h with ultrapure water
2. 1 h with acetonitrile-water mixture (30:70)
3. 1 h with ultrapure water

Inorganic contamination

1. Add 30% acetonitrile to the standard eluent.
2. Regenerate the column in the direction opposite to the flow at a flow rate of 1 mL/min for 1 h.

General notes

- Sample solutions must be microfiltered (0.45 µm filter).
- The column is not resistant to alcohol additives in the eluent.
- The eluent has only a limited shelf life due to bacterial and fungal growth. We therefore recommend preparing fresh eluent regularly.
- To protect the separation column, we recommend using the pulsation absorber

(6.2620.150) to reduce the injector pressure surges.

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