Adjusting a sample rack in tiamo / Probenrack justieren in tiamo

EN

Adjusting the sample rack

It could happen that the rack position openings of the 874 Oven Sample Processor do not precisely match the opening of the oven. The sample vial could become jammed as a result.

This lack of exactitude can be rectified by adjusting the rack. The rack can be adjusted by correcting the rack offset in *tiamo*[™].

Moving to a sample position

Proceed as follows:

- Open the manual control in *tiamo*[™].
- Under **Sample changers ▶ 874_1**, click on **Tower 1**.
- On the tab **Move**, enter the **Target position 1** and click on **[Start]**.

If the opening of rack position 1 is now located considerably to the left or right of the middle of the oven, then the value of the rack offset must be corrected in the rack configuration.

Note


Carry out the correction of the rotation angle of the rack (**Rack offset $-5.0^\circ \dots +5.0^\circ$**) in *tiamo*.

1 mm deviation of the rack to the left requires a correction of the rotation angle by -0.3° .

1 mm deviation of the rack to the right requires a correction of the rotation angle by $+0.3^\circ$.

Correcting the rack offset

- In *tiamo*[™], select the **Configuration**.
- Double-click on the device name **874_1**.
- Click the tab **Rack**.
- Click **[Rack data]** to open the settings of the rack parameters.

Rack data

Rack name:

Rack code:

Number of positions:

Rack parameters Lift positions Special beakers

Beaker radius samples: mm

Beaker sensor:

Rack offset: °

- Correct the value for **Rack offset** according to the observed deviation.
- Close the rack data configuration with **[OK]**.
- Click on **[Initialize rack]**.

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