

# Oxygenated hydrocarbons

## Separation of DME, MTBE and ETBE

### Application Note

Energy & Fuels

#### Authors

Agilent Technologies, Inc.

#### Introduction

Gas chromatography using an Agilent CP-SilicaPLOT column separates DME, MTBE and ETBE in 30 minutes.




**Agilent Technologies**

## Conditions

Technique : GC-wide-bore  
Column : Agilent CP-SilicaPLOT, 0.32 mm x 30 m fused silica PLOT (df = 4  $\mu$ m) (Part no. CP8567)  
Temperature : 90 °C (10 min)  $\rightarrow$  250 °C, 10 °C/min  
Carrier Gas : He, 150 kPa (1.5 bar, 21 psi)  
Injector : Split,  
T = 250 °C  
Detector : ECD  
T = 275 °C  
Solvent Sample : pentane  
Courtesy : J. Kuipers, Agilent Technologies Inc.,  
Middelburg, The Netherlands

## Peak identification

1. C<sub>6</sub>
2. C<sub>7</sub>
3. dimethyl ether (DME)
4. C<sub>8</sub>
5. C<sub>9</sub>
6. C<sub>10</sub>
7. methyl-tertiary-butyl-ether (MTBE)
8. ethyl-tertiary-butyl-ether (ETBE)
9. C<sub>11</sub>
10. C<sub>12</sub>



[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01637



**Agilent Technologies**