



# Biodiesel according to ASTM D6584

## Application Note

Energy & Fuels

### Authors

Agilent Technologies, Inc.

### Introduction

Analysis of biodiesel (fatty acid methyl esters) using Agilent VF-5ht with very low bleed at 380°C.




**Agilent Technologies**

## Conditions

Technique : GC-capillary  
Column : Agilent FactorFour VF-5ht, 0.32 mm x 15 m fused silica (df = 0.10  $\mu$ m) (Part no. CP9047)  
Temperature : 50 °C (1 min)  $\rightarrow$  180 °C, 15 °C/min, 7 °C  $\rightarrow$  230 °C, 30 °C  $\rightarrow$  380 °C  
Carrier Gas : Helium, 60 kPa  
Injector : Cold on-column, 1  $\mu$ L  
Detector : FID  
T = 380 °C  
Sample Size : 100 mg per 10 mL hexane  
Courtesy : C. Duvekot, Agilent application laboratory, Middelburg, The Netherlands

## Peak identification

1. C16:0
2. C18:0 18:1 18:2
3. C20:0 20:1
4. C22:0 22:1
5. C24:0 24:1



[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A02347



**Agilent Technologies**