
Application Bulletin

Of interest to: General analytical laboratories; Food

A 1, 7

Automatic determination of the formol number in fruit and vegetable juices

Summary

The formol number represents a further parameter for the characterization of fruit and vegetable juices. As this is merely an index (the formol number does not deal with the molecular size, nor with the quantity of amino acids), the conditions of the titration can be adapted to practical needs. This concerns mainly the pH value of the endpoint of the SET titration (pH = 8.5, pH = 9.0, pH = 9.2, etc.).

Instruments and accessories

- 702 SET/MET Titrino, 716 DMS Titrino, 719 SET Titrino, 736 GP Titrino, 751 GPD Titrino or 785 DMP Titrino or 726 or 796 Titroprocessor with 700 Dosino or 685 Dosimat
 - 765 or 776 auxiliary Dosimat for adding the formaldehyde solution
 - 2.728.0040 Magnetic Stirrer
 - 6.3014.223 Exchange Units
 - 6.0255.100 combined LL double-junction pH glass electrode with 6.2104.020 electrode cable
-

Reagents

- Titrant: $c(\text{NaOH}) = 0.1 \text{ mol/L}$
 - Formaldehyde solution: $w(\text{HCHO}) = 35\%$, adjusted to pH = 8.5 with NaOH
-

Analysis

Pipet 25.0 mL sample solution into a glass beaker and titrate in a first SET titration to pH = 8.5 with $c(\text{NaOH}) = 0.1 \text{ mol/L}$. By means of the activate pulse the second SET titration is then automatically started under the following conditions:

Addition of 15 mL formaldehyde solution (auxiliary Dosimat), 60 s waiting time, titration to pH = 8.5 with $c(\text{NaOH}) = 0.1 \text{ mol/L}$.

Calculation

The formol number is equivalent to the consumption of $c(\text{NaOH}) = 0.1 \text{ mol/L}$ per 100 mL sample solution.

Formol number = $\text{EP1} * \text{C01}$

EP1 = titrant consumption for the second SET titration in mL

C01 = 4 (for a sample volume of 25 mL)

Figures

```
'pa
719 S Titrino      OP1/110  719.0020
date 1999-05-11  time 09:34   26
SET pH           Formal1
parameters
>SET1
  EP at pH        8.50
  dynamics        0.5
  max.rate        10.0 ml/min
  min.rate        10 µl/min
  stop crit:      drift
  stop drift      20 µl/min
>SET2
  EP at pH        OFF
>titration parameters
  titr.direction: auto
  pause 1         0 s
  start V:        OFF
  pause 2         0 s
  extr.time       0 s
  meas.input:     1
  temperature     25.0 °C
>stop conditions
  stop V:         abs.
  stop V          99.99 ml
  filling rate    max. ml/min
>statistics
  status:         OFF
>preselections
  conditioning:   OFF
  req.ident:      OFF
  req.smpl size:  OFF
  activate pulse: OFF
=====
```

Fig. 1: Parameter settings on the 719 SET Titrino for the first SET titration.

```
'pa
719 S Titrino      OP1/110  719.0020
date 1999-05-11  time 09:34   26
SET pH           Formal2
parameters
>SET1
  EP at pH        8.50
  dynamics        0.5
  max.rate        10.0 ml/min
  min.rate        10 µl/min
  stop crit:      drift
  stop drift      20 µl/min
>SET2
  EP at pH        OFF
>titration parameters
  titr.direction: +
  pause 1         60 s
  start V:        OFF
  pause 2         0 s
  extr.time       0 s
  meas.input:     1
  temperature     25.0 °C
>stop conditions
  stop V:         abs.
  stop V          99.99 ml
  filling rate    max. ml/min
>statistics
  status:         OFF
>preselections
  conditioning:   OFF
  req.ident:      OFF
  req.smpl size:  OFF
  activate pulse: OFF
=====
```

Fig. 2: Parameter settings for the second SET titration.

```

'cr
719 S Titrino      OP1/110  719.0020
date 1999-05-11   time  08:26   17
meas.input:      1      CAL  *****
cal.date         1999-05-11
                    pH      U/mV
buffer 1         7.00     15
buffer 2         4.00     185
cal.temp         22.5 °C
slope(rel)       0.968     pH(as)  7.26
                    =====

'fr
719 S Titrino      OP1/110  719.0020
date 1999-05-11   time  09:41   27
pHc(init)        3.85     SET pH  Formal1
EP1              31.420 ml  8.51
                    =====

'fr
719 S Titrino      OP1/110  719.0020
date 1999-05-11   time  09:46   28
pHc(init)        7.04     SET pH  Formal2
EP1              4.842 ml  8.50
FZ              19.368
                    =====
    
```

Fig. 3: Result report for the determination of the formol number in orange juice.

Literature

- Schweizerisches Lebensmittelbuch, Kapitel 28
Frucht- und Gemüsesäfte, Fruchtnektare, Fruchtsirupe, Konzentrate und Pulver
Abschnitt 9.4 Bestimmung der Formolzahl (1990).