Thermo. Titr. Application Note No. H-107

Title: Determination of Ferrous Ion in Hydrometallurgical Leach Liquors

Scope:	Determination	of	the	ferrous	ion	content	of
	hydrometallurgical leach liquors						

Principle:	A measured amount of acidic hydrometallurgical leach liquor is further acidified with sulfuric acid, prior to being titrated with standard potassium dichromate solution to an exothermic endpoint.
	$\begin{array}{l} Cr_2O_7^{2-} + \ 14H^+ + \ 6e \leftrightarrow 2Cr^{3+} + \ 7H_2O \\ \underline{[Fe^{2+} \leftrightarrow Fe^{3+} + e] \ x \ 6} \\ Cr_2O_7^{2-} + \ 6Fe^{2+} + \ 14H^+ \leftrightarrow 2Cr^{3+} + \ 6Fe^{3+} + \ 7H_2O \end{array}$
	Thus, 1 mol $K_2Cr_2O_7 \equiv 6$ mol Fe^{2+}

Reagents:	Titrant: Standard potassium dichromate solution		
	$c(K_2Cr_2O_7) = 0.1 \text{ mol/L}$		
	Acid: ~25% v/v sulfuric acid		

Method:	Basic Experimental Parameters:				
	Titrant delivery rate (mL/min.)	4			
	No. of endothermic endpoints	1			
	Data smoothing factor (DSF)	50			
	Stirring speed (802 stirrer)	8			
	Delay before start of titration (secs.) 10			
	Basic titration procedure. A 15mL aliquot of acidic process liquor is pipetted by volumetric glass pipette into a PP titration tube, and 5mL 25% v/v H ₂ SO ₄ solution plus 10mL DI water added. The sample solution is swirled to mix prior to being placed in the sample rack.				

Example:

Acidic hydrometallurgical leach liquor, containing Fe(II), Fe(III)), Mg, Al, Mn, Cr, Cu, Co and Ca.

A Metrohm

Fe(II) = 10.66±0.02 g/L (n=5)

Calculation:

Fe, g/L = $((EP \text{ vol., mL- Blank, mL}) \times c(K_2Cr_2O_7) \text{ mol/L} \times AW \text{ Fe} \times 6)$ Sample vol., mL

