

IC Application Note No. M-4

Title: Traces of bromide and bromate in drinking water by IC-MS, determination of the Method Detection Limit (MDL)

Summary: Determination of bromide and bromate in drinking water applying anion chromatography with MS detection.

Sample: Drinking water

Sample Preparation: Direct injection

Column: 6.1006.510 Metrosep A Supp 5 – 100


Eluent: 25 mmol/L lithium hydroxide
25% acetonitrile

Suppressor: Metrohm Suppressor Module (MSM, 100 mmol/L sulfuric acid, 2% acetonitrile)


Flow: 0.7 mL/min

Injection Volume: 100 μ L

MS detection: Extracted ion chromatogram for bromate ($m/z = 127$, red) and bromide ($m/z = 79$, blue)



MS detection: Total ion chromatograms used for the determination of the Method Detection Limit



Results:	Bromate µg/L	Bromide µg/L
MDL – 1	0.1630	0.1790
MDL – 2	0.1688	0.1800
MDL – 3	0.1720	0.1800
MDL – 4	0.1687	0.1885
MDL – 5	0.1720	0.1848
MDL – 6	0.1679	0.1854
MDL – 7	0.1700	0.1799
Average	0.16891	0.18252
Standard Deviation	0.00306	0.00368
Method Detection Limit (MDL)	0.0096	0.011