


Speciation of Fe(II) and Fe(III) in soil applying IC-ICP/MS according to EPA SW-846 Method 6800



The speciation analysis of iron is crucial because the oxidation state determines the environmental behavior. Iron uptake by organisms as well as its transport and storage depend to a great extent on its oxidation state. Here, the two species are separated on a Metrosep A Supp 10 S-Guard/4.0 column. For detection, the IC-ICP/MS with speciated isotope dilution methodology is applied.

Results

	Fe(III) [$\mu\text{g}/\text{kg}$]	Fe(II) [$\mu\text{g}/\text{kg}$]
Standard solution	250	250

Sample

Standard

Sample preparation

Filtration (0.45 µm)

Columns

Metrosep A Supp 10 S-Guard/4.0 6.1020.510

IC Solutions

Eluent 4.0 mmol/L dipolonic acid
20.0 mmol/L ammonium nitrate, pH = 4.3

Parameters

Flow rate 0.8 mL/min

Injection volume (MiPT) 100 µL

P_{max} 25 MPa

Recording time 5 min

Parameters ICP/MS

RF power 1550 W

Plasma gas flow rate 15 L/min

Carrier gas flow rate 0.95 L/min

Makeup gas flow rate 0.15 L/min

Collision gas (He) flow rate 6.0 L/min

Sampling depth 8.0 mm

Spray chamber temperature 2 °C

Tuning solution 1 µg/L Li, Co, Y, Ce
Ti in 2% HNO₃
solution

Acquisition mode Spectrum and time
resolved analysis

Monitoring mass Fe 56 amu

Analysis

ICP/MS detection

Instrumentation

940 Professional IC Vario ONE 2.940.1100

ICP-MS Agilent 7700

858 Professional Sample Processor 2.858.0020

Remote box 6.2148.010

Remote cable Professional IC - MS-
Detector (Agilent) 6.2141.380

