

IC Application Note No. O-40

Title: Borate in seawater

Summary: Determination of borate in seawater using ion exclusion chromatography with suppressed conductivity detection after Inline Matrix Elimination.

Sample: Seawater

Sample Preparation: Dilution 1:20 with ultrapure water, the borate is injected onto an Organic Acids Guard column, the matrix (mainly chloride) is washed from the column with ultrapure water.


Column: 6.1005.200 Metrosep Organic Acids – 250

Eluent: 0.2 mmol/L perchloric acid / 50 mmol/L mannitol

Suppressor: Metrohm Suppressor Module (MSM, 100 mmol/L Li₂CO₃)

Flow: 0.6 mL/min

Injection Volume: 10 µL onto a 6.1005.250 Metrosep Organic Acids Guard



Results:	Borate as boron mg/L
	10.6