

Application of single particle ICP-MS for silver nanoparticles characterization

Chiu-Hun Su^{1*}, Chia-Chin (Donna) Hsu², Chun-Hua Chen¹, Meng-Tzu Lai¹, Yen-Ying (Brian) Laio², Ching-Heng (Jones) Hsu³

¹ Industrial Technology Research Institute, Hsinchu, 31040, Taiwan

² Agilent Technologies Taiwan Ltd., Taoyuan City 32450, Taiwan

³ BASF, Taoyuan City, 32853, Taiwan.

Corresponding author email: ch.su@itri.org.tw; donna.hsu@agilent.com

ABSTRACT


The nanomaterials are widely used in various products of human life such as water purification equipment, anti-bacteria clothes and wearing. In this work, to understanding the behavior and fate of nanoparticles in environmental, we investigated the particle size distribution and particle concentration of silver nanoparticles (Ag NPs) in four kinds of acidic matrix and environmental waters via spICP-MS. Comparison of the particle size stability of Ag NPs that spiked in HCl, CH₃COOH, HNO₃ and H₂SO₄ solutions (the pH value is about 5) for 72 hrs, the relative standard deviation (RSD) of particle size are 6.3%、9.2%、9.6% and 15.4%, respectively. It indicates that Ag NPs stored in H₂SO₄ matrix are the most unstable, followed by HNO₃ matrix. And the disintegration rate of Ag NPs increase when pH value decreases. Furthermore, we spiked Ag NPs in industrial waste water, rain water (pH value is 5.29), artificial acid rain water (mixture of H₂SO₄, HNO₃ and NH₄OH and pH value is 5.24), tap water and ultrapure water. In short- and long-term stability tests, we found that Ag NPs spiked in industrial waste water are the most unstable, and the short-term (16~18 hrs) stability of Ag NPs is good in ultrapure water and tap water. In conclusion, with the present study we have demonstrated that Ag NPs in environmental waters cannot exist for a long-time, they might collapse into Ag ions or transformed in different types of particles, namely Ag₂S and AgCl.

Single Nanoparticle Analysis


Agilent 7900 ICP-MS



Principle of spICP-MS



Workflow of spICP-MS



Experimental

Chemicals and Nanomaterials

Silver nanoparticles of 40 and 100 nm (Sigma-Aldrich) were acquired in stock suspensions at a nominal concentration of 20 mg Ag/L. Nano-Ag suspensions were made by diluting the stock solutions with 18.2 M-ohm Nanopure water (Millipore to final concentrations, ranging from 1.25 to 100 ppt. Aqueous Ag standards (High-Purity Standards), used for calibration, were diluted in 1% nitric acid to concentrations 1 ppb. Nitric acid (HNO₃) and Hydrochloric acid (HCl) (36%) from BASF, as well as Sulfuric acid (H₂SO₄) (98%) and Acetic acid (CH₃COOH) (99.8%) from J.T. Barker were used to prepare the sample.

Instrumentation


An Agilent 7900 with Micromist nebulizer was used for single particle (SP) analysis, run at 1550 W RF power.



Dilution of AgNPs


NP samples should give 500 ~ 2000 peaks/min

Conc.	Particle Conc. (particles/mL)	# of Particles	Mean Size (nm)	Mean Size RSD(%)
1.25 ppt	2503	68	39	0.6%
2.5 ppt	4483	122	40	0.6%
5 ppt	9479	259	39	0.6%
15 ppt	27970	763	38	1.0%
50 ppt	94609	2580	39	1.7%
100 ppt	180692	4928	40	0.3%




Repeatability test of AgNPs

The RSD of the Mean Size and the Particle Concentration at five replicated runs are 1.1% and 3.1%, respectively.




Stability testing of AgNPs (pH effects)

The disintegration rate of Ag NPs increase when pH value decreases.




Short-term stability test of AgNPs

The RSD of the Mean Size and the Particle Concentration in 51 minutes are 0.3% and 2.8%, respectively.




Stability testing of AgNPs (Matrix)

Ag NPs stored in H₂SO₄ matrix are the most unstable, followed by HNO₃ matrix.




Stability testing of AgNPs in environmental waters




Comparison Results

Ag NPs spiked in industrial waste water are the most unstable.

Mean size of AgNPs that spiked in Ultrapure water, Rain water, Tap water, Industrial waste water and Artificial acid rain water



of Particles of AgNPs that spiked in Ultrapure water, Rain water, Tap water, Industrial waste water and Artificial acid rain water



Reference

1. S. Wilbur, M. Yamanaka and S. Sannac, *Agilent publication*, **2015**, 5991-5516EN.
2. Ariel R. Donovan, Craig D. Adams, Yinfa Ma, Chady Stephan, Todd Eichholz, Honglan Shi, *Chemosphere*, **2016**, 144, 148–153.
3. Vera V. Pinto, Maria José Ferreira, Ricardo Silva, Hélder A. Santos, F. Silva, Carlos M. Pereira, *Colloids and Surfaces A: Physicochem. Eng. Aspects*, **2010**, 364, 19–25.