

SBSE/SFE-SFC-MS/MS System to Analyze Migrating Plastic Additives in Medical Solutions and Biological Matrices



Benjamin Caux^{1,2}, Clément De Saint Jores¹, Ramy Abou Naccoul², Shinnosuke Horie³, Toshiyuki Yamashita⁴, Takeshi Bamba⁴, Caroline West¹


¹University of Orleans, CNRS UMR7311, ICOA, rue de Chartres, BP 6759, 45067 Orléans cedex 2, France

²Shimadzu France, Le lizzard 2, Bat A, Bd Salvador Allende Noisiel, 77448 Marne-la-Vallée, France

³Shimadzu Europa GmbH, Albert-Hahn-Straße 6-10, 47269 Duisburg, Germany

⁴Division of Metabolomics, Medical Institute of Bioregulation, Kyushu University, Fukuoka, Japan

The use of **plastic additives** in the formulation of plastic materials is essential to obtain specific physico-chemical properties. However, certain compounds such as plasticizers can **migrate** from a plastic container to the content and can be **harmful for health** especially in the context of medical devices. Stir Bar Sorptive Extraction (SBSE) well described in literature as a technique to adsorb lipophilic compounds such as phthalates, which are common plasticizers. Thermo-desorption or solvent extraction is usually employed to desorb the stir-bar, followed by GC or LC analysis. The use of **supercritical CO₂ (scCO₂)** is an innovative way to desorb analytes from the stir-bar, followed by **online analysis by SFC-MS/MS**, offering several advantages.




Nexera UC




Design of Experiments

Equilibrium between the solution and the PDMS phase is an important factor to optimize in SBSE. The addition of methanol, salt, and stirring time can be adjusted to modify the equilibrium and enhance the extraction yield. A **Box-Behnken design** was performed using these **three factors**. For many plasticizers, such as DiBP and DEHP, neither methanol nor NaCl was required, but a minimum stirring time of 60 minutes was necessary to obtain optimal adsorption. An **aqueous solution of glucose** was spiked with plastic additives and used as sample to mimic a **hospital nutritive bag**.




SBSE Rinsing method

Normally, stir-bars are "cleaned" with solvents or at high temperatures to completely desorb residual contaminants. Here, it is proposed to use the **online SFE-SFC-MS/MS** system to directly **clean the adsorbent**, followed by an analysis for blank checking. To clean the stir bars, scCO₂ with methanol followed by **pure scCO₂** was used (for less than one hour). The eluent composition and cleaning time would need to be optimized depending on the contaminants.



Saliva


The SBSE method was applied to biological matrices such as saliva. The objective was to analyze the composition of orthodontic elastic bands using the online **SFE-SFC-MS/MS** method with pure scCO₂ to extract plastic additives and to compare this composition with the **migrated compounds** extracted and analyzed using the **SBSE-SFE-SFC-MS/MS** method. Two orthodontic bands were immersed in saliva to simulate long-term use in real-life conditions.



% relative area is calculated in function of targeted compounds

Human plasma

The SBSE method was applied to **human plasma** as a biological matrix. Some compounds were adsorbed onto the stir-bar and then desorbed using **pure scCO₂**. However, the signal area was reduced due to **matrix effects**. Plasma macromolecules can directly interact with analytes and modify the adsorption equilibrium, or macromolecules can get adsorbed onto the stir-bar and compete with the analytes.



The use of **SBSE** allows the adsorption of **migrated plastic additives** in liquid matrices such as **medical solutions** (nutritive bags) or **biological matrices** (saliva, plasma). The use of **supercritical CO₂ (scCO₂)** is an innovative approach to desorb analytes from the stir-bar, followed by **online analysis using SFC-MS/MS**. The next objective is to reduce matrix effects from biological samples in order to **improve extraction efficiency** and **sensitivity**.