

Analysis on Cephalosporins on a Polar Endcapped C18 3 μ m Particle Column

Valeria Barattini, Thermo Fisher Scientific, Runcorn, Cheshire, UK

Key Words

Synchronis, aQ, cephalosporins, polar endcapping, antibiotics

Abstract

The analysis of species containing polar groups using reversed phase HPLC can be challenging. The polar groups present on the analytes can result in secondary interactions with silanols on the silica surface. Furthermore, the overall polarity of the molecule sometimes calls for the use of mobile phases with very high aqueous content to ensure the analyte solubility, retention and elution, at the expense of stationary phase stability. In this application note we demonstrate the analysis of a mix of four analytes from the cephalosporin antibiotics family on a Thermo Scientific™ Synchronis™ aQ 3 μ m particle column. The analysis yields excellent peak shape and resolution in short run times, as well as outstanding reproducibility.

Introduction

One of the key goals for the chromatographer is to achieve a consistent, reproducible separation. The selection of a highly reproducible HPLC column is essential if this goal is to be attained. The Synchronis column range has been engineered to provide exceptional reproducibility due to its highly pure, high surface area silica, dense bonding and double endcapping, all controlled and characterized through the use of rigorous testing. The aQ phase endcapping group has been further optimized to provide a controlled interaction mechanism that retains and resolves polar analytes, with added media stability in 100% aqueous mobile phase.

Cephalosporins are a class of antibiotics constituting a class of β -Lactam antibiotics first discovered and isolated from the fungus *Cephalosporium Acremonium*. The core structure of these antibiotics (Figure 1) features a number of hydrogen bond donors and acceptor groups, as well as an acid function. Together these groups contribute towards the polarity of this class of compounds.





Figure 1: Core structure of cephalosporins



A group of four cephalosporins, listed in Figure 2, was chosen to be separated using a Synchronis aQ column with a 3 μ m particle. Two of the components (cephalexin and cefradine) differ only in the presence of a phenyl group in cephalexin compared to a 1,4-cyclohexadiene in cefradine.




Figure 2: Cephalosporins mix used for this analysis

Thermo
SCIENTIFIC

Experimental Details

Consumables		Part Number
Cephalexin, cephradine, cefazolin sodium salt and cephaloridine were all of a suitably high purity (97% or above)		
Fisher Scientific analytical reagent grade ammonium acetate		A/3440/53
Fisher Scientific OPTIMA grade acetic acid		A465-250
Fisher Scientific HPLC grade acetonitrile		A/0627/17
Thermo Scientific Smart2Pure® purified water (18.2 MΩ)		
Sample Handling Equipment		Part Number
Liquid handling hardware:	FinnPipette Kit 1	4700870
Vials and closures:	8 mm standard opening screw thread vial Convenience kit	60180-600 pack of 100
Separation Conditions		Part Number
Instrumentation:	Thermo Scientific Dionex Ultimate 3000 HPLC system	
Column:	Synchronis aQ, 100 x 3.0 mm, 3 µm	97303-103030
Mobile phase A:	25 mM ammonium acetate, adjusted to pH 5 with acetic acid	
Mobile phase B:	Acetonitrile	
Gradient:	Time (min)	%B
	0.0	10
	3.0	40
	4.0	40
	4.1	10
	6.0	10
Flow rate:	0.5 mL/min	
Run time:	6.0 min	
Column temperature:	35 °C	
Injection details:	1 µL of cephalosporins mix	
Injection wash solvent:	Water:acetonitrile 9:1 (v/v)	
UV detector wavelength:	254 nm	
Solutions		
Individual solutions of cephalexin, cephradine, cefazolin and cephaloridine were prepared by diluting 20 mg of compound (± 1 mg) in 20 mL 2:8 (v/v) of acetonitrile : water (acidified with 0.1% acetic acid). Each standard solution was sonicated for 10 minutes. A mixture of all four cephalosporins was prepared by diluting 100 µL of each standard in 600 µL of water (acidified with 0.1% acetic acid).		

Results

Under the conditions described above, separation of the cephalosporins mixture above was achieved with a total run time of six minutes, with an elution window of 1.8 minutes for all four analytes. The separation chromatogram is shown in Figure 3. Good retention and separation of these polar compounds is achieved, demonstrating the controlled interaction mechanism and selectivity provided by the Synchronis aQ polar endcapping group.




Figure 3: Separation chromatogram of a mixture of cephalosporins on a Synchronis aQ 3 μ m particle column

In Figure 4 the chromatograms obtained from six replicate injections are overlaid and % RSD values for the retention time of each analyte listed in the table inset. Excellent reproducibility was found, with % RSD values not exceeding 0.19%.




Figure 4: Overlay of six replicate injections of a mixture of four cephalosporins on Synchronis aQ 3 μ m particle column

Outstanding low asymmetry values were found for each analyte (Table 1), demonstrating the superior interaction mechanism and selectivity provided by Synchronis aQ polar endcapping group and lack of unwanted secondary interactions with the silica substrate.

	Average Asymmetry
Cephalexin	1.083
Cephadrine	1.057
Cefazolin	1.075
Cephaloridine	1.035

Table 1: Average peak asymmetry value measured for each cephalosporin analyte

Conclusion

- Synchronis aQ 3 μm provides excellent resolution power within short run times
- The highly pure, high surface area silica combined with polar endcapping provides good retention and superior peak shape for polar analytes
- Exceptional reproducibility of results is achieved.

thermoscientific.com/chromatography

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

USA and Canada +1 800 332 3331
France +33 (0)1 60 92 48 34
Germany +49 (0) 2423 9431 20 or 21
United Kingdom +44 (0)1928 534110
Japan +81 3 5826 1615

China +86 21 68654588 +86 10 84193588
+86 20 83145199 800 810 5118
India +91 22 6742 9494 +91 27 1766 2352
Australia 1 300 735 292 (free call domestic)
New Zealand 0800 933 966 (free call domestic)
All Other Enquiries +44 (0) 1928 534 050

Technical Support
North America +1 800 332 3331
Outside North America +44 (0) 1928 534 440

Thermo
 SCIENTIFIC
 Part of Thermo Fisher Scientific