

HILIC Method for the Separation of Diquat and Paraquat

Monica Dolci, Thermo Fisher Scientific, Runcorn, Cheshire, UK

Key Words

- Paraquat
- Diquat
- Quaternary amines
- Synchronis HILIC
- EPA

Abstract

This application note demonstrates the use of the Thermo Scientific Synchronis HILIC column for the fast analysis of the quaternary ammonium herbicides paraquat and diquat.

Introduction

Paraquat and diquat are non-selective contact herbicides, widely used to control crop and aquatic weeds [1]. Paraquat and diquat have high water solubility, making them easily absorbed into the soil and potentially into drinking water supplies. Because of their significant toxicity to humans, these compounds have been banned, or their use restricted, in several European countries and in Japan.

The World Health Organization considers paraquat and diquat moderately hazardous pesticides and they must be routinely monitored [2].

The structures of these compounds are shown in Figure 1.




Figure 1: Structure of paraquat (A) and diquat (B)

These highly charged quaternary amines cannot be retained in conventional reversed-phase.

An alternative separation and retention method is offered by hydrophilic interaction liquid chromatography (HILIC), a technique suitable for the retention of very polar and hydrophilic compounds. In HILIC the stationary phase is hydrophilic and often charged; the mobile phase typically consists of 60-95% acetonitrile in water, or a volatile buffer. In HILIC, water (or the aqueous buffer) acts as the strong eluent, necessary to establish a water-rich liquid layer near the stationary phase surface. Separation is achieved by partitioning of analytes from the mobile phase to the water-rich layer. When using charged HILIC stationary phases, secondary electrostatic interactions can also contribute to retention and selectivity. Zwitterionic stationary phases are successfully employed in HILIC because the overall effect they provide is of weak electrostatic interactions and lower buffer concentrations are required to disrupt these interactions.



Experimental Details

The separation of paraquat and diquat was achieved on a Synchronis HILIC column specifically designed for the retention of polar and hydrophilic compounds.

Chemicals and Reagents	Part Number
Fisher Scientific HPLC grade water	W/0106/17
Fisher Scientific HPLC grade acetonitrile	A/0626/17
Thermo Scientific 2 mL clear vial and Si/PTFE seal	60180-600

Sample Preparation

A 1400 µg/mL standard solution of paraquat and 1800 µg/mL standard solution of diquat were prepared in mobile phase; these solutions were then mixed in a 50:50 ratio and used for the analysis.

Separation Conditions	Part Number
Instrumentation:	Thermo Scientific Accela HPLC/UHPLC
Column:	Synchronis HILIC 5 µm, 100 x 2.1mm 97505-102130
Column Temperature:	45 °C
Injection volume:	1 µL
Flow rate:	0.3 mL/min
UV detection:	254 nm (for paraquat) and 308 nm (for diquat)
Mobile phase:	20:80 (v/v) water/acetonitrile + 0.1% TFA

Synchronis™ HILIC columns are based on highly pure 100 Å silica, with a surface area of 320 m²/g. The zwitterionic modified stationary phase results in total charge equalisation and therefore a neutral but highly polar surface, which offers enhanced retention of polar and hydrophilic analytes.

This application note shows an efficient and reproducible HILIC method for the analysis of paraquat and diquat, with excellent peak shape and baseline resolution.

Results

The herbicides paraquat and diquat were retained and separated on the Synchronis HILIC column, using an isocratic method with 0.1% TFA added to the mobile phase. Figure 2 shows the chromatogram obtained employing Synchronis HILIC 5 µm, 100 x 2.1 mm column, with paraquat eluting at 1.30 min and diquat eluting at 1.47 min.

The following chromatographic parameters were monitored: retention time t_R , peak area, and resolution R_s . Mean and % RSD values for the above parameters were based on data derived from six replicate injections, and they are reported in Table 1.




Figure 2: Chromatogram of paraquat and diquat, separated on a Synchronis HILIC 5 µm, 100 x 2.1 mm column

Conclusions

The use of a Synchronis HILIC column allowed the successful retention and separation of the quaternary ammonium compounds paraquat and diquat. Synchronis HILIC columns are an excellent choice for the HILIC separation of paraquat and diquat, affording highly reproducible analysis

References

- [1] Determination of diquat and paraquat in drinking water by Liquid-Solid Extraction and High Performance Liquid Chromatography with Ultraviolet Detection; U.S. EPA Method 549.2, Revision 1.0, Environmental Protection Agency: Cincinnati, OH, 1997.
- [2] <http://www.who.int/whosis/en/> World Health Organization.

Paraquat	t_R (min)	Peak Area
Mean	1.30	7361563
SD	0.01	98638.00
% RSD	0.43	1.34

Diquat	t_R (min)	Peak Area	R_s
Mean	1.47	7414798	1.61
SD	0.01	104719.97	0.01
% RSD	0.50	1.41	0.66

Table 1: Average and Method Precision (%RSD) of chromatographic parameters, derived from the analysis of paraquat and diquat on a Synchronis HILIC 5 µm, 100 x 2.1 mm column (data calculated from six replicate injections)

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

**North America
USA and Canada**
+1 800 332 3331

**Europe
France**
+33 (0)1 60 92 48 34

Germany
+49 (0) 2423 9431 -20
or -21

United Kingdom
+44 1928 534110

**Asia
Japan**
+81 3 5826 1615

China
+86-21-68654588
or +86-10-84193588
800-810-5118

India
+91-22-6742 9494

**Thermo Fisher
Scientific Australia
Pty Ltd**
1300 735 292 (free call
domestic)

**Thermo Fisher
Scientific New
Zealand Ltd**
0800 933 966 (free call
domestic)

All Other Enquiries
+44 (0) 1928 534 050

Technical Support

North America
800 332 3331

**Outside North
America**
+44 (0) 1928 534 440

www.thermoscientific.com/chromatography

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

ANCCSDIQAQSYN 0112