

## **Impurities test for Amoxicillin (EP-8.0 method):**

### **SAMPLE PREPARATION:**

**Test Solution (a):** Dissolve 30 mg of Metronidazole in Mobile phase A and dilute to 50ml with mobile phase A.

**Test Solution (b):** Dissolve 30 mg of Metronidazole in Mobile phase A, and dilute to 20ml with mobile phase A.

**Reference solution (a):** Dissolve 30 mg of Metronidazole in Mobile phase A and dilute to 50ml with mobile phase A.

**Reference solution (b):** Dissolve 4mg of cefadroxil CRS in mobile phase A, , dilute to 50 ml with mobile phase A, To 5ml of this solution add 5ml of reference solution (a) and dilute to 100 ml with mobile phase A.

**Reference solution (c):** Dilute 2ml of reference solution (a) to 20 ml with mobile phase A; dilute 5ml of this solution to 20 ml with mobile phase A.

### **CHROMATOGRAPHIC CONDITIONS:**

**Instrument:** UltiMate 3000 LC

**Column:** Hypersil Gold (4.6\*250mm, 5um, p/n 25005-254630, lot no.: 13285)

#### **Mobile phase:**

**Mobile phase A:** Mix 1 volumes of Acetonitrile and 99 volumes of 25 % v/v solution of 0.2M potassium Dihydrogen phosphate adjusted to pH 5.0 with dilute sodium hydroxide.

**Mobile phase B:** Mix 20 volumes of Acetonitrile and 80 volumes of 25 % v/v solution of 0.2M potassium Dihydrogen phosphate adjusted to pH 5.0 with dilute sodium hydroxide.

**Separation Mode:** Gradient

Time (min)	Mobile phase A (% v/v)	Mobile phase B (% v/v)
0	92	8
10	92	8
35	0	100
50	0	100
52	92	8
65	92	8

**Column temperature:** 25°C

**Flow rate:** 1.0 mL/min

**Injection Volume:** 50 µl

**Detector wavelength:** UV 254 nm


**Run Time:** 65min

#### **System Suitability Results:**

Sr. No.	Parameters	USP Criteria	Obtained Results
1	Resolution b/w Amoxicillin and Cefadroxil peak in reference solution b	NLT 2.0	10.01
2	RRT of Impurity C and Impurity J with respect to Amoxicillin peak	<ul style="list-style-type: none"><li>• Imp C- 3.4</li><li>• Imp-J- 4.1</li></ul>	<ul style="list-style-type: none"><li>• 3.2</li><li>• 3.8</li></ul>

CHROMATOGRAMS:

System Suitability:



Impurity Mix:

