

Impurities test for Caffeine (EP-8.0 method):

SAMPLE PREPARATION:

Test Solution: Dissolve 0.1gm of Caffeine in Mobile phase, dilute to 50 ml, and dilute 1ml of this solution to 10ml with mobile phase.

Reference solution (a): Dilute 2ml of test solution to 100ml with mobile phase, dilute 1ml of this solution to 10ml with mobile phase.

Reference solution (b): Dissolve 5mg of caffeine for system suitability CRS (containing impurity A, C, D and F) in mobile phase and dilute to 5ml, Dilute 2ml of this solution to 10ml with mobile phase.

CHROMATOGRAPHIC CONDITIONS:

Instrument: UltiMate 3000 LC

Column: Acclaim 120 C18 (4.6*150mm, 5um, p/n 059148, lot no.:018-01-171)

Mobile phase: Acetonitrile: Tetrahydrofuran: Buffer (25:20:955)

Buffer: 0.82 g/L of anhydrous sodium acetate, adjust pH to 4.5 with Glacial acetic acid.

Separation Mode: Isocratic

Column temperature: 25°C

Flow rate: 1.0 mL/min

Injection Volume: 10 µl

Detector wavelength: UV 275nm


Run Time: 15min

System Suitability Results:

Sr. No.	Parameters	USP Criteria	Obtained Results
1	<ul style="list-style-type: none">Resolution between the impurities C & D.Resolution between the impurities F & A	<ul style="list-style-type: none">NLT 2.5NLT 2.5	<ul style="list-style-type: none">4.44.2
2	Tailing Factor for caffeine	NMT 2.0	1.1

CHROMATOGRAMS:

System Suitability:



Impurity Mix:

