

Impurities test for Metronidazole (EP-8.0 method):

SAMPLE PREPARATION:

Test Solution: Dissolve 0.05gm of Metronidazole in Mobile phase and dilute to 100ml with mobile phase

Reference solution (a): Dilute 1ml of test solution to 100ml with mobile phase, dilute 1ml of this solution to 100ml with mobile phase.

Reference solution (b): Dissolve 5mg of Metronidazole impurity A CRS in mobile phase, add 10 ml of test solution, dilute to 100ml with mobile phase, and dilute 1ml to 100 ml with mobile phase.

CHROMATOGRAPHIC CONDITIONS:

Instrument: UltiMate 3000 LC

Column: Synchronis C18 (4.6*250mm, 5um, p/n 97105-254630, lot no.: 11811)

Mobile phase: 30:70 (Methanol: 1.36 g/L solution of potassium Dihydrogen phosphate.

Separation Mode: Isocratic

Column temperature: 25°C

Flow rate: 1.0 mL/min

Injection Volume: 10 µl

Detector wavelength: UV315 nm


Run Time: 21min

System Suitability Results:

Sr. No.	Parameters	USP Criteria	Obtained Results
1	Resolution b/w Metronidazole and Impurity-A	NLT 2.0	6.9
2	Tailing Factor for Metronidazole peak	NMT 2.0	1.0

CHROMATOGRAMS:

System Suitability:



Impurity Mix:

