

Impurities test for Hydrochlorothiazide (EP- method):

SAMPLE PREPARATION:

Reference solution-a: Dissolve 3mg of Impurity-A and 3mg of Hydrochlorothiazide in 5ml of mixture of equal volumes of acetonitrile and methanol if necessary dilute to 20ml with phosphate buffer solution pH 3.2. Dilute 5ml of this solution to 100ml with solvent mixture.

Test Solution-a: Dissolve 30mg of Hydrochlorothiazide in 5ml of mixture of equal volumes of acetonitrile and methanol if necessary dilute to 20ml with phosphate buffer solution pH 3.2.

Reference solution-b: Dilute 1ml of test solution (a) to 100ml with the solvent mixture. Dilute 1 ml of this solution to 10ml with solvent mixture.

Solvent Mixture: Dilute 50ml of mixture of equal volumes of acetonitrile and methanol to 200ml of phosphate buffer solution pH 3.2.

CHROMATOGRAPHIC CONDITIONS:

Instrument: UltiMate 3000 LC

Column: Hypersil Gold (4.6*100mm, 3 μ m, p/n 25003-104630, lot no.:12337)

Mobile phase A: Buffer: Methanol: Tetrahydrofuran (940:60: 10)

Mobile phase B: Buffer: Methanol: Tetrahydrofuran (500:500: 50)

Diluent: Solvent mixture and phosphate buffer pH 3.2

Phosphate Buffer: 35.8 g/l of disodium hydrogen phosphate in water, adjust with phosphoric acid to a pH of 3.2. Dilute 100 mL of this solution to 2000 mL with water.

Separation Mode: Gradient

Time (min)	Mobile phase A (% v/v)	Mobile phase B (% v/v)
0	100	0
17	55	45
30	55	45
30.1	100	0
45	100	0

Column temperature: 25°C

Flow rate: 0.8 mL/min

Injection Volume: 10 μ l



Detector wavelength: UV 224nm

Run Time: 45 min

System Suitability Results:


Sr. No.	Parameters	USP Criteria	Obtained Results
1	Resolution between Hydrochlorothiazide & impurity-A in reference solution (a)	NLT 2.5	5.74
2	RRT of Impurity A, B and C with reference to Hydrochlorothiazide	Impurity-A= about 0.9 min Impurity-B= about 0.7 min Impurity-C= about 2.8 min	Imp-A=0.88 Imp-B= 0.67 Imp-C =2.94

CHROMATOGRAMS:



System Suitability:

Impurity Mix:



Reference Solution-b: