

## **Impurities test for Ampicillin Sodium (EP-8.0 method):**

### **SAMPLE PREPARATION:**

**Reference solution (a):** Dissolve 27.0 mg of Ampicillin in Mobile phase A and dilute to 50ml with mobile phase A.

**Reference solution (b):** Dissolve 2mg of cefradine CRS in mobile phase A, , dilute to 50 ml with mobile phase A, To 5ml of this solution add 5ml of reference solution (a).

**Reference solution (c):** Dilute 1.0 ml of reference solution (a) to 20 ml with mobile phase A.

**Reference solution (b):** To 0.2gm of Ampicillin add 1 ml of water. Heat the solution to 60°C for 1 hour. Dilute 0.5 ml of this solution to 50 ml with mobile phase A.

### **CHROMATOGRAPHIC CONDITIONS:**

**Instrument:** UltiMate 3000 LC

**Column:** Acclaim120-C18 (4.6\*250mm, 5um, p/n 059149, lot no.: 018-01-152)

#### **Mobile phase:**

**Mobile phase A:** Mix 0.5 ml of dilute acetic acid, 50ml of 0.2M potassium Dihydrogen phosphate and 50 ml of acetonitrile, then dilute to 1000ml with water.

**Mobile phase B:** Mix 0.5 ml of dilute acetic acid, 50ml of 0.2M potassium Dihydrogen phosphate and 400 ml of acetonitrile, then dilute to 1000ml with water.

**Separation Mode:** Solution (b) and solution (c) with Isocratic elution at Initial mobile phase composition for 20 minutes.

Solution (d), Impurity Mix solution and individual impurities as per the gradient described below:

Time (min)	Mobile phase A (% v/v)	Mobile phase B (% v/v)
0	90*	10*
11	90	10
41	0	100
56	0	100
57	90	10
70	90	10

**\*Note:** As described in Monograph the initial mobile phase composition of A and B can be adjusted to achieve the resolution, so the composition of A and B was finalized as 90 and 10 respectively.

**Column temperature:** 25°C

**Flow rate:** 1.0 mL/min

**Injection Volume:** 50 µl

**Detector wavelength:** UV 254 nm


**Run Time:** 70min

#### **System Suitability Results:**


Sr. No.	Parameters	EP Criteria	Obtained Results
1	Resolution b/w Ampicillin and Cefradine	Minimum 3.0	3.23
2	RRT with reference to ampicillin and ampicillin dimer in reference solution (d)	About 2.8	2.84

**CHROMATOGRAMS:**

**System Suitability:**



**Impurity Mix:**



**Reference solution (d):**

