

Impurities test for Cefradine (EP- method):

SAMPLE PREPARATION:

Reference solution-a: Dissolve 3.0 mg of Impurity-B in Mobile phase A and dilute to 50 ml with mobile phase A.

Reference solution-b: Dissolve 3.0 mg of Cefradine and Cefalexin monohydrate in Mobile phase A and dilute to 25 ml with mobile phase A.

Reference solution-c: Dilute 1ml of test solution to 100 ml with mobile phase A

Test Solution: Dissolve 0.3 g of Cefradine in mobile phase A and dilute to 50 ml with mobile phase A.

Reference solution-e: Dissolve content of vial of cefradine impurity mixture (Impurity A and G) in 1 ml mobile phase.

CHROMATOGRAPHIC CONDITIONS:

Instrument: UltiMate 3000 LC

Column: Hypersil Gold (4.6*150mm, 5 um, p/n 25005-154630, lot no.:12750)

Mobile phase A: 2.72 g/L solution of potassium Dihydrogen phosphate adjusted to pH 3 with OPA.

Mobile phase B: Methanol

Separation Mode: Gradient

Time (min)	Mobile phase A (% v/v)	Mobile phase B (% v/v)
0	99.5	0.5
2.5	97	3
11	75	25
13	60	40
16	60	40
19	20	80
19.1	99.5	0.5
25	99.5	0.5

Column temperature: 30°C

Flow rate: 1.0 mL/min

Injection Volume: 25 µl



Detector wavelength: UV 220nm

Run Time: 25min

System Suitability Results:


Sr. No.	Parameters	EP Criteria	Obtained Results
1	Resolution between Cefalexin and Cefradine	Minimum 4.0	7.32
2	Retention Time of Cefradine	About 15 minutes	15.057 min
3	RRT with reference to Cefradine	Imp-A= about 0.27, Imp-B= about 0.32 min, Imp-C= about 0.53 min, Imp-D=about 0.63 min, Imp-E= about 0.80 min, Imp-F= about 0.92 min, Cefalexin= about 0.95, 4',5' dihydrocefradine= about 1.06 min, Imp-G= about 1.32 min	Imp-A= 0.23 min, Imp-B= 0.28 min Imp-C= 0.53 min, Imp-D= 0.59 min Imp-E= 0.78 min, Imp-F= 0.98 min Cefalexin= 0.95 min 4',5' dihydrocefradine= 1.06 min Imp-G= 1.18 min

CHROMATOGRAMS:



Reference solution –b (System Suitability)

Impurity Mix



Reference solution –e