

## Impurities test for Betamethasone (EP- method)

### SAMPLE PREPARATION:

**Test solution:** Dissolve 25mg of Betamethasone in a mixture of equal volumes of Acetonitrile and Methanol and dilute to 10 ml with same solution.

**Reference solution-a:** Dissolve 2 mg of Betamethasone and 2 mg of Methylprednisolone in mobile phase A then dilute to 100 ml with mobile phase A.

**Reference solution-b:** Dilute 1ml of test solution to 100ml with Mobile phase A.

### CHROMATOGRAPHIC CONDITIONS:

**Instrument:** UltiMate 3000 LC

**Column:** Synchronis C18 (4.6\*250mm, 5 um, p/n 97105-254630, lot no.:11811)

**Mobile phase A:** Water: Acetonitrile (250:750)

**Mobile phase B:** Acetonitrile

**Separation Mode:** Gradient

Time (min)	Mobile phase A (% v/v)	Mobile phase B (% v/v)
0	100	0
15	100	0
40	0	100
41	100	0
46	100	0

**Column temperature:** 45°C

**Flow rate:** 2.5 mL/min

**Injection Volume:** 20 µl

**Detector wavelength:** UV 254nm


**Run Time:** 46 min

**System Suitability Results:**


Sr. No.	Parameters	EP Criteria	Obtained Results
1	Resolution between Methylprednisolone and Betamethasone peak	Minimum 1.5	1.55*

\* The obtained results are keeping all the chromatographic parameters same as per requirement of USP monograph, however as per monograph, the required resolution can further be increased by changing the acetonitrile concentration in mobile phase A.

### CHROMATOGRAMS:



**System Suitability:**



**Impurity Mix:**