

# Multidraw Infinity II Bio LC Upgrade Kit (G7137-68711)

## Installation Note

In this note we describe how to install the Multidraw Infinity II Bio LC Upgrade Kit (G7137-68711) into an Agilent 1290 Infinity II Bio Multisampler (G7137A).

## Contents

**General Information** 2

**Delivery Checklist** 3

**Installing the Multidraw Infinity II Bio LC Upgrade Kit** 4

**Configuration of the Controller** 5

OpenLab CDS ChemStation C.01.10 Update 3 and OpenLab 2.5 and above 6  
Lab Advisor B.02.16 and above 6



## General Information

The Multidraw Infinity II Bio LC Upgrade Kit can be installed in the Agilent 1290 Infinity II Bio Multisampler.

With the kit, you can add a maximum of 500  $\mu\text{L}$  or 1500  $\mu\text{L}$  to the injection volume (100  $\mu\text{L}$  analytical head) of your injector.

The maximum injection volume will vary, depending on the analytical head and sample loop installed.

### NOTE

The delay volume of your Multisampler is extended when using the extended seat capillaries from the multidraw kit. When calculating the delay volume of the Multisampler, double the volume of the extended capillaries (used in partial filling mode). The delay volume can be reduced by bypassing the Multisampler once the sample has reached the head of the column, see your *User Manual* for more information.

---

When using the Multisampler in multidraw mode, the syringe ejects an equivalent of the injected volume into waste. Therefore it is recommended to verify the correct installation of the waste tube to the waste outlet of the Multisampler.

# Delivery Checklist


**NOTE**

Do not connect a waste outlet directly to the plastic waste tube (connected to port 4 of the injection valve). This can lead to siphoning effects that will influence the precision of the Multisampler.


Make sure that all parts and materials have been delivered with the upgrade kit. Report missing or damaged parts to your local Agilent Technologies sales and service office.

p/n	Description
G7137-68307	Extended seat capillary, MP35N, 500 $\mu$ L, 0.5 mm id
G7137-68308	Extended seat capillary, MP35N, 1500 $\mu$ L, 0.9 mm id
5023-2625	Union MP35N

Extended seat capillary  
500  $\mu$ L



Extended seat capillary  
1500  $\mu$ L  
(larger external diameter  
than 500  $\mu$ L loop)







Union

**NOTE**


Both extended seat capillaries have SL fittings, compatible with VICI and IDEX/Rheodyne Valves.

# Installing the Multidraw Infinity II Bio LC Upgrade Kit


Tools required	p/n	Description
	Wrench, 1/4 inch	¼ inch wrench (one supplied in the autosampler accessory kit)

<p><b>1</b> Open the front door.</p> 	<p><b>2</b> Disconnect the seat-capillary fitting from the injection valve (port 5).</p> 
<p><b>3</b> Install the 500 or 1500 µL seat extension loop (longer side) to the injection valve (port 5).</p> 	<p><b>4</b> Install the union between the seat capillary and the other side of the seat extension loop. Store the extension loop in the leak tray.</p> 

**5** Connect the waste tube to the waste guidance.



**6** Close the front door.



**7** Connect the waste tube (shipped with the start up kit of your Multisampler).


## Configuration of the Controller

The configuration of your controller is necessary to enable the multidraw mode.

When setting an injection larger as the configured injection volume the multidraw mode is active. Multiples of the injection syringe volume are stored in the extended seat capillary prior to switching the injection valve.

## OpenLab CDS ChemStation C.01.10 Update 3 and OpenLab 2.5 and above

- 1 Select **Capillaries** in the instrument function.




- 2 In the **Configuration** menu, change seat capillary to the value of the installed extended seat capillary and click **OK**.

## Lab Advisor B.02.16 and above

- 1 Connect to the stack in **System Overview**.
- 2 Select **Instrument Control** and open **Controls** of the Multisampler.

- 3 On the **Configuration** tab, select the volume for your installed extended seat capillary (500  $\mu$ L or 1500  $\mu$ L) from the **Set Capillary** drop-down menu.



The screenshot shows a configuration window titled "Set Capillaries". It contains two dropdown menus. The "Seat Capillary" dropdown is currently set to "High-Pressure Seat Bio-Compatible 0.12 (G7137-87012)". The "Loop Capillary" dropdown is currently set to "High-Pressure Seat Bio-Compatible 0.12 (G7137-87012)". A "Send" button is located to the right of the dropdown menus.

Field	Selected Option	Available Options
Seat Capillary:	High-Pressure Seat Bio-Compatible 0.12 (G7137-87012)	High-Pressure Seat Bio-Compatible 0.12 (G7137-87012) Multi-draw Seat Capillary MP35N 500ul (G7137-68307) Multi-draw Seat Capillary MP35N 1500ul (G7137-68308)
Loop Capillary:	High-Pressure Seat Bio-Compatible 0.12 (G7137-87012)	High-Pressure Seat Bio-Compatible 0.12 (G7137-87012) Multi-draw Seat Capillary MP35N 500ul (G7137-68307) Multi-draw Seat Capillary MP35N 1500ul (G7137-68308)

**Configuration of the Controller**  
Lab Advisor B.02.16 and above



Part Number: G7137-90110 Rev. A  
Document Number: D0004814 Rev. A

Edition: 07/2020  
Printed in Germany

© Agilent Technologies, Inc 2020

Agilent Technologies, Inc  
Hewlett-Packard-Strasse 8  
76337 Waldbronn  
Germany