

HPLC, LC/MS Columns

InertSustain C18



The First Choice ODS Column





The First Choice ODS Column

It is high inertness, low adsorption, and wide pH range.

InertSustain C18

InertSustain C18 inherits the advantages of all current Inertsil HPLC columns (e.g., extremely low operating back pressure, superior inertness to almost all analytes, high efficiency, and compatibility with a wide range of solvents) while additionally enabling wide pH analysis with consistent column-to-column and lot-to-lot performance.

In general, silica based columns are mechanically stable and highly efficient but cannot be used under alkaline conditions as their residual silanol groups tend to adsorb organic bases. InertSustain C18 employs a radically new type of silica that is uniquely surface-modified for precise control of the silica properties.


- Base Material : High Purity ES Silica Gel
- Particle Size : 2 μm , 3 μm , 3.5 μm , 4 μm , 5 μm , 7 μm , 10 μm
- Surface Area : 350 m^2/g
- Pore Size : 100 Å (10 nm)
- Pore Volume : 0.85 mL/g
- Functional Group : Octadecyl
- End-capping : Yes
- Carbon Loading : 14 %
- USP Code : L1
- pH Range : 1 - 10




Reliable Reproducibility, Performance and Quality

Rigorous quality control of physical properties and strict chromatographic tests for inertness and selectivity, contribute to the production of InertSustain C18 with an outstanding reproducibility and long column lifetime.


Strong Basic Compound Test




Same Retention Behavior between Various Particle Sizes



Strong Acidic Compound Test




Strong Chelating Compound Test




Wide pH Compatibility with Long Column Lifetime

As shown in the experiment below, due to the introduction of Evolved Surface Silica, InertSustain C18 maintained high efficiency and peak shape for 300 hours while other "wide pH" column brands failed.



Column Performance Comparison


Basic Compounds



Conditions

Column Size : 3 μ m, 150 $\times$ 2.1 mm I.D.
 Eluent : A) CH₃CN
 B) 25 mM K₂HPO₄ (pH 7.0, KH₂PO₄)
 A/B = 30/70, v/v
 Flow Rate : 0.2 mL/min
 Col. Temp. : 40 $^{\circ}$ C
 Detection : UV 230 nm
 Sample : 1. Uracil
 2. Pyridine
 3. Phenol
 4. Berberine chloride
 5. Dextromethorphan


Acidic Compounds



Conditions

Column Size : 3 μ m, 150 $\times$ 2.1 mm I.D.
 Eluent : A) CH₃CN
 B) 0.1 % H₃PO₄
 A/B = 25/75, v/v
 Flow Rate : 0.2 mL/min
 Col. Temp. : 40 $^{\circ}$ C
 Detection : UV 254 nm
 Sample : 1. Brilliant Blue FCF
 2. Phenol
 3. Salicylic acid

Chelating Compounds



Conditions

Column Size : 5 μ m, 150 $\times$ 4.6 mm I.D.
 Eluent : A) CH₃CN
 B) 0.1 % H₃PO₄
 A / B = 40 / 60, v / v
 Flow Rate : 1.0 mL / min
 Col. Temp. : 40 $^{\circ}$ C
 Detection : UV 254 nm
 Sample : 1. Hinokitiol

Please refer to our website for complete list of our columns.

https://www.glsciences.com/product/lc_columns/c18_column/01708.html



Contact us or your local GL Sciences representative.

<https://www.glsciences.com/contactus/>

Authorized distributor:

<https://www.glsciences.com/company/distributor.html>

For analytical use only. Information subject to change without prior notice.



<https://www.glsciences.com>
 E-Mail world@glsc.co.jp