

# Determination of Mineral Oil in Water According to ISO 9733-2:2000 Using ISOLUTE® Na<sub>2</sub>SO<sub>4</sub>/FL Solid Phase Extraction Columns

## Introduction

This application note describes the extraction and clean-up of mineral oil (petroleum oil) from water samples according to ISO 9377-2:2000. This method involves a liquid-liquid extraction step followed by clean up using an SPE column containing layers of sodium sulfate and florisil, to remove halogenated and other polar interferences from the extract.

## Sample Preparation Procedure

### Format

ISOLUTE® Na<sub>2</sub>SO<sub>4</sub>/FL 4 g/6 mL, Part Number 976-0400-C

### Pre-Treatment

1. Cool sample (250–1000 mL) to ~ 10 °C (if required, to prevent volatilization).
2. Acidify to pH 2 with HCl. This step may be performed in the field at collection.
3. Add magnesium sulfate to the sample (80 g/900 mL, as required (see Additional Information section)).
4. Add hexane (50 mL), close bottle, and shake or stir vigorously for 30 mins. Further treatment such as centrifugation may be required if an emulsion is formed.
5. Remove hexane layer and transfer to a pre-solvated SPE column for clean-up.

### Column Solvation

Rinse column with hexane (5 mL)

### Sample Application

Apply the hexane extract to the column, and allow to flow under gravity. Low vacuum or pressure may be used (optimize flow rate as required).

### Analyte Elution

Elute hydrocarbons with hexane (2 x 10 mL), under gravity or low vacuum/pressure (optimise flow rate as required).

## Post Extraction

1. Concentrate the extract using a rotary evaporator with controlled vacuum to a volume of approximately 6 mL.
2. Concentrate further to just below 1 mL using a gentle flow of nitrogen.
3. Make up to exactly 1 mL with hexane, or calculate the exact volume of the concentrated extract by weighing.
4. Transfer an aliquot of the final extract to a septum vial for GC analysis. If a high HC oil index is expected or if a large volume injection system is used, concentration of the extract to 1 mL is not required, but a known volume should be used.

## Chemicals and Reagents

- » Hexane
- » Magnesium Sulfate

## Additional Information

1. Hexane is suggested as extraction solvent, however other non-branched alkanes may be suitable.
2. Samples may need to be preserved in the field by adjustment with mineral acid (e.g. HCl) to pH 2. If high levels of humic acids are present however, acidification should be avoided to prevent precipitation of the humic substances.
3. Addition of magnesium sulfate is only required if the samples are known to emulsify.
4. For full details of the method, see ISO 9377-2:2000 Determination of hydrocarbon oil index. Part 2: Method using solvent extraction and gas chromatography.





## Ordering Information

Part Number	Description	Quantity
976-0400-C	ISOLUTE® Na <sub>2</sub> SO <sub>4</sub> /FL 4 g/6 mL	30

### EUROPE

Main Office: +46 18 565900  
 Toll Free: +800 18 565710  
 Fax: +46 18 591922  
 Order Tel: +46 18 565710  
 Order Fax: +46 18 565705  
 order@biotage.com  
 Support Tel: +46 18 56 59 11  
 Support Fax: + 46 18 56 57 11  
 eu-1-pointsupport@biotage.com

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
 Toll Free: +1 800 446 4752  
 Fax: +1 704 654 4917  
 Order Tel: +1 704 654 4900  
 Order Fax: +1 434 296 8217  
 ordermailbox@biotage.com  
 Support Tel: +1 800 446 4752  
 Outside US: +1 704 654 4900  
 us-1-pointsupport@biotage.com

### JAPAN

Tel: +81 3 5627 3123  
 Fax: +81 3 5627 3121  
 jp\_order@biotage.com  
 jp-1-pointsupport@biotage.com

### CHINA

Tel: +86 21 68162810  
 Fax: +86 21 68162829  
 cn\_order@biotage.com  
 cn-1-pointsupport@biotage.com

### KOREA

Tel: +82 31 706 8500  
 Fax: +82 31 706 8510  
 korea\_info@biotage.com  
 kr-1-pointsupport@biotage.com

### INDIA

Tel: +91 11 45653772  
 india@biotage.com

Distributors in other regions  
 are listed on [www.biotage.com](http://www.biotage.com)

### Literature Number: IST1077.V.1

© 2020 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.