

FTIR SPECTROSCOPY REFERENCE GUIDE

The Measure of Confidence

INTRODUCTION

Infrared spectroscopy is the study of the interactions between infrared electromagnetic energy and matter. The technique of infrared spectroscopy measures the vibrations of molecules, allowing for qualitative and quantitative measurements of samples. A Fourier transform infrared (FTIR) spectrometer is an ideal tool for the identification of unknown organic and inorganic samples whether they exist in the form of a gas, liquid or a solid.

ELECTROMAGNETIC SPECTRUM

Radiation in the infrared region is commonly referred to in terms of a unit called a wavenumber ($\bar{\nu}$), rather than wavelength (λ). Wavenumbers are expressed as reciprocal centimeters (cm^{-1}) and are the preferred unit as they are directly proportional to energy. A higher wavenumber corresponds to a higher energy.


Convert a wavenumber ($\bar{\nu}$) to a wavelength (λ) by using the following formula:

$$\bar{\nu} (\text{cm}^{-1}) = \frac{1}{\lambda (\text{cm})}$$

Convert a wavenumber ($\bar{\nu}$) to Energy (E) by multiplying it by the speed of light (c, in centimeters per second) and Planck's constant (h):

$$E (\text{J}) = h (\text{J}\cdot\text{s}) \times c (\text{cm}\cdot\text{s}^{-1}) \times \bar{\nu} (\text{cm}^{-1})$$

High ← Energy → Low
High ← Wavenumber ($\bar{\nu}$) → Low



Short ← Wavelength (λ) → Long

A portion of the electromagnetic spectrum showing the relationship of the infrared region to other types of radiation (not shown to scale). The relationship between energy, wavenumber, and wavelength is also highlighted.

COLLECTING AN FTIR SPECTRUM

There are three quick, simple steps involved in obtaining a spectrum of a sample:

STEP 1 Record a spectrum with no sample present. (Known as a 'background'.)

STEP 2 Insert the sample into the spectrometer.

STEP 3 Record a second spectrum and interpret the data.


The key collection parameters are:

- The number of scans (background and sample)
- The scan range
- The spectral resolution.

In general increasing the number of scans that are co-added improves the signal-to-noise ratio of the spectrum and can assist in the analysis of weakly absorbing samples. The spectral resolution is user-defined and helps to distinguish closely spaced absorption peaks, and is expressed in wavenumbers.

The simplified correlation table on the right allows users to extract structural information from IR spectra. Computer-based search programs are also available for assisting in compound identification.

CARY 600 SERIES SPECTROMETER OPTICS




INTERPRETING SPECTRA

Bond	Type of Vibration	Wavenumber Range (cm^{-1})	Bond	Type of Vibration	Wavenumber Range (cm^{-1})	
C-H	Alkane (stretch)	3000 – 2850	C-O	Alcohols, esters, ethers, carboxylic acid, anhydrides	1300 – 1000	
	-CH ₃ (bend)	1450 & 1375		O-H	Alcohols, phenols	3650 – 3600
	-CH ₂ (bend)	1465	Free		3400 – 3200	
	Alkene (stretch)	3100 – 3000	H-Bonded		3400 – 2400	
	C=C	Alkene (out-of-plane bend)	1000 – 650	N-H	Primary & secondary amines & amides (stretch)	3500 – 3100
Aromatic (stretch)		3150 – 3050	(bend)		1640 – 1550	
C≡C	Alkyne (out-of-plane bend)	900 – 600	C-N	Amines	1350 – 1000	
	Alkyne (stretch)	~ 3300		C=N	Imines & oximes	1690 – 1640
C=O	Aldehyde	1740 – 1720	C≡N	Nitriles	2260 – 2240	
	Ketone	1725 – 1705	N=O	Nitro (R-NO ₂)	1550 & 1350	
	Carboxylic acid	1725 – 1700		S-H	Mercaptans	2550
	Ester	1750 – 1730	C-X		Halides	
	Amide	1680 – 1630			Fluoride	1400 – 1000
Anhydride	1810 & 1760	Chloride		785 – 540		
			Bromide, iodide	< 650		

SPECTRAL RANGE COVERAGE

Spectral range coverage for common sources, beamsplitters, and detectors. The combination of these components generate a spectrometer's working scan range.




THE COMPLETE SOLUTION



Several different accessories facilitate sample preparation and spectral acquisition, including:

- Attenuated Total Reflectance (ATR)
- Diffuse Reflectance
- Specular Reflectance
- Grazing Angle Reflectance
- Microscopy and Chemical Imaging
- ATR Chemical Imaging (Micro and Macro)
- Fiber Optic Probes
- PM-IRRAS
- Photoacoustic Spectroscopy
- TGA-FTIR
- GC-FTIR
- GPC-FTIR

A HISTORY OF COMMERCIAL 'FIRSTS'



This information is subject to change without notice. © Agilent Technologies, Inc. 2013. Published in Agilent, February 1, 2013. 59900-0000

