

# Manual and Automated Hybrid Searches in NIST26 Video/Handout

James Little


Mass Spec Interpretation Services

April 25, 2026

[mzinterpretation.com](http://mzinterpretation.com)

See **Full Course** on NIST26 with new **Integrated** Deconvolution/Library Searching for **EI GC-MS** and **LC-MS/MS**!

## Mass Spec ( $m/z$ ) Interpretation Services Organic Mass Spectrometry



## Hybrid MS/MS Search: Extending Library Coverage

Conventional MS/MS searching depends on exact spectral matches, limiting identification when compounds are absent from libraries. Hybrid MS/MS search overcomes this by recognizing ***structurally related compounds***, expanding effective library coverage.

In the NIST Mass Spectral Search Program, the hybrid search is available both ***manually*** in the Library Search window and ***automatically*** within the Chromatogram workflow.


The use of  **$\Delta$ form (dform)** highlights mass differences corresponding to plausible structural modifications, guiding interpretation. These capabilities enable more confident exploration of unknowns beyond exact matches.

# Additional Information

## Novel NIST Hybrid Search for Unknown Identifications

[See Link](#)

### Example Identifying Compound with NIST Hybrid Search



[Excel DeltaMass Table](#)

[PDF DeltaMass Table](#)

NIST initially developed the hybrid search and I have listed their publication

#### Videos

[YouTube: Advanced NIST Hybrid Search of EI Spectra](#)

[YouTube: Advanced NIST Hybrid Search for MS/MS \(Tandem\) Spectra](#)

[YouTube: Application NIST Hybrid in EI in Seized Drug Analyses](#)

#### Associated Handouts

[Handout: Advanced NIST Hybrid Search of EI Spectra](#)

[Handout: Advanced NIST Hybrid Search of MS/MS \(Tandem\) Spectra](#)

[Handout: Application NIST Hybrid in EI in Seized Drug Analyses](#)


Supporting Software Courses for Utilizing Hybrid and Adaptive Search

[NIST EI Search Software](#)


[NIST MSMS Search Software](#)

# Manually Sending One Spectrum to Lib Search from Chromatogram

- Perform standard search in Chromatogram
- Find, select, and display NoID's
- Sort by Abundance Relative
- Step through with keyboard arrows and butterfly plots
- Send ones of interest manually to Lib Search for Hybrid search to find candidates



# Many related to PEG type Materials, Note Precursor for $(M+NH_4)^+$ and loss of $NH_3$ , $MH^+ -17$



| #  | Scan | RT      | Score | Score (Unfiltered) | Abund.Rel. | Prec. m/z | dPPM | Prec. |
|----|------|---------|-------|--------------------|------------|-----------|------|-------|
| 1  | 1907 | 9.1726  | 0     | 0                  | 2.88       | 511.4714  | 0.0  |       |
| 2  | 1909 | 9.2624  | 0     | 0                  | 2.35       | 528.4980  | 0.0  |       |
| 3  | 1908 | 9.1736  | 0     | 0                  | 2.01       | 497.4558  | 0.0  |       |
| 4  | 1420 | 7.9990  | 0     | 0                  | 1.87       | 576.2362  | 0.0  |       |
| 5  | 2801 | 11.6457 | 0     | 0                  |            |           |      |       |
| 6  | 2660 | 11.3833 | 0     | 0                  |            |           |      |       |
| 7  | 2694 | 11.3542 | 0     | 0                  |            |           |      |       |
| 8  | 1626 | 8.4805  | 0     | 0                  |            |           |      |       |
| 9  | 1170 | 6.9208  | 0     | 0                  |            |           |      |       |
| 10 | 2026 | 9.6909  | 0     | 0                  |            |           |      |       |
| 11 | 1399 | 7.5871  | 0     | 0                  |            |           |      |       |
| 12 | 2883 | 11.9782 | 0     | 0                  |            |           |      |       |
| 13 | 2769 | 11.5863 | 0     | 0                  |            |           |      |       |
| 14 | 2658 | 11.3514 | 0     | 0                  |            |           |      |       |
| 15 | 2738 | 11.4931 | 0     | 0                  |            |           |      |       |
| 16 | 1844 | 8.9222  | 0     | 0                  |            |           |      |       |

- Library Search
- Library Search Options
- Send To
- Copy Selected Hits to Clipboard
- Export Selected Hits to Text File
- Show Selected
- Show All
- Properties

Library Search Options

**Search**

MS/MS

Libraries

Automation

Spectrum Search Type

Identity  Similarity

MS/MS Hybrid

Precursor Ion m/z

In spectrum

Spectrum Scoring Options


Method: Full Spectrum Search (Score)

Compound Ubiquity Correction


- Most reasonable candidate is a dodecyl ether of PEG (polyethylene glycol)
- Sent Hit 1 to MS interpreter to understand fragmentation
- dForm showing loss of three PEG groups for Delta Mass of 3 x EG group
- Observed dForm 132.079, theoretical 132.079
- Sent proposed structure to drawing program and removed 3 EG repeat groups
- Sent to modified structure to MS Interpreter for confirmation

Definition Below\*

| # | Library        | S... | DotProd | DeltaMass | dForm      | pctRelForm | Prec. Type           | Instr. Type | Energy | DBs | Name                               |
|---|----------------|------|---------|-----------|------------|------------|----------------------|-------------|--------|-----|------------------------------------|
| 1 | hr_msms_nist#2 | 871  | 908     | -132.079  | C-6H-12O-3 | 0.25       | [M+NH4] <sup>+</sup> | HCD         | 20%    | 25  | Nonaethylene glycol monododecyl... |




Best Candidate



Reasonable Structure

\*pctRefForm Percent: Probability that a dForm(ula) is correct. Based on statistics of occurrences of each dForm relative to the most common (CH<sub>2</sub> which is assigned 100%) – for high scoring hits (>900). The higher, the most likely correct – values like 0.01 are so rare, they are very likely accidental and wrong for example

# Proposed Structure, Changed Imported Structure to Protonated, Default is Ionized, Loss of ammonia *1.2 ppm difference*



# General Comments for Ammonium in MSMS

- **More (M+H)<sup>+</sup> than (M+NH<sub>4</sub>)<sup>+</sup> adducts for species in NIST MS/MS Library**
- NIST **does not purposely add** ammonium when acquiring reference spectra
- Often better when user sees loss of 17 for NH<sub>3</sub>
- Change specified ion to the (M+H)<sup>+</sup> value noted for fragmentation
- Thus, for the previous example, unclick “In spectrum” and insert value for (M+H)<sup>+</sup> noted in spectrum
- Top two hits still monodecyl ethers
- **The major fragments in MS/MS of (M+NH<sub>4</sub>)<sup>+</sup> and (M+H)<sup>+</sup> are the same**

Library Search Options

Search

MS/MS

Libraries

Automation

Limits

Constraints

Presearch

Spectrum Search Type

Identity  Similarity

MS/MS Hybrid

Precursor Ion m/z


In spectrum 451.3624

Spectrum Scoring Options


Method Full Spectrum Search (Score)

Compound Ubiquity Correction

| # | Library        | Score | DotProd | DeltaMass | dForm      | pctRelForm | Prec. Type         | Instr. Type | Energy | DBs | Name                                   |
|---|----------------|-------|---------|-----------|------------|------------|--------------------|-------------|--------|-----|--|
| 1 | hr_msms_nist#2 | 838   | 899     | 44.0257   | C2H4O      | 5.1        | [M+H] <sup>+</sup> | HCD         | 20%    | 12  | Pentaethylene glycol monododecyl ether |
| 2 | hr_msms_nist#2 | 812   | 871     | -132.079  | C-6H-12O-3 | 0.25       | [M+H] <sup>+</sup> | HCD         | 15%    | 25  | Nonaethylene glycol monododecyl ether  |






1



2

## Running Hybrid Search in Automated Chromatogram Approach

- Select Input File with Hybrid Search checked
- Select **“Run in Background”**
- Note will take **significantly more time** than standard method
- Note user alert shows **“Running Background Task”**
- Alerted when finished
- **Properties** in menu **automatically displayed**
- See results sorted by NoID, selected, sorted by abundance
- **Same hit shown** for component done manually for comparison
- Very useful to **send to Lib Search** Window for manual search **for structures**



| # | Scan | RT      | Score | Score (Unfiltered) | Abund.Rel. | Prec. m/z | dPPM | Prec. Type | HScore | HdMass    | HdForm     | HpctRelForm | HPrec. Type          | HName                                 | Lib ID      |
|---|------|---------|-------|--------------------|------------|-----------|------|------------|--------|-----------|------------|-------------|----------------------|---------------------------------------|-------------|
| 1 | 1907 | 9.1726  | 0     | 0                  | 2.88       | 511.4714  | 0.0  |            | 683    | 196.2184  | C14H28     | 0.014       | [M+H] <sup>+</sup>   | Dihexyl adipate                       | NoID_511.47 |
| 2 | 1909 | 9.2624  | 0     | 0                  | 2.35       | 528.4980  | 0.0  |            | 573    | 213.2450  | -          |             | [M+H] <sup>+</sup>   | Dihexyl adipate                       | NoID_528.50 |
| 3 | 1908 | 9.1736  | 0     | 0                  | 2.01       | 497.4558  | 0.0  |            | 719    | 182.2028  | C13H26     | 0.086       | [M+H] <sup>+</sup>   | Dihexyl adipate                       | NoID_497.46 |
| 4 | 1420 | 7.9990  | 0     | 0                  | 1.87       | 576.2362  | 0.0  |            | 30     | 196.0950  | C14H12O    | 0.19        | [M+H] <sup>+</sup>   | Licofelone                            | NoID_576.24 |
| 5 | 2801 | 11.6457 | 0     | 0                  | 1.64       | 468.3887  | 0.0  |            | 871    | -132.0794 | C-6H-12O-3 | 0.25        | [M+NH4] <sup>+</sup> | Nonaethylene glycol monododecyl ether | NoID_468.39 |