



Volatile amines are an occupational health hazard. There are several methods published by the U.S. Occupational Health and Safety Administration (OSHA). Typical methods use some kind of trapping reagent in an air-sampling tube. In this case, 4-chloro, 7-nitro, 1,2,3-benzoxadiazole (NBD-chloride) reacts with amines in the trap, to produce a highly-colored, fluorescent derivative. OSHA methods describe the use of normal-phase LC, but reversed-phase LC can provide equivalent or better results. As shown here, the Acclaim C30 resolved four common volatile amines with simple, isocratic reversed-phase conditions. The unreacted reagent elutes at 5 minutes and does not interfere with the analytes of interest. This new method is compatible with four existing methods (34, 36, 40 and 41), while providing improved peak shape and faster analysis time. Linearity is excellent from 0.1 to 75 µg/mL, with detection limits below 0.05 µg/mL.