

# Comprehensive and Confident Identification of Narcotics, Steroids and Pharmaceuticals in Urine

David E. Alonso<sup>1</sup>, Petra Gerhards<sup>2</sup>, Charles Lyle<sup>1</sup> and Joe Binkley<sup>1</sup> | <sup>1</sup>LECO Corporation, St. Joseph, MI; <sup>2</sup>LECO European LSCA Centre, Moenchengladbach, Germany

## Introduction

Monitoring of patients in hospitals and clinics has traditionally relied on targeted methods of analysis. These screening methods are not comprehensive and result in an incomplete picture of a patient's activities. Gas chromatography high resolution time-of-flight mass spectrometry (GC-HRT) provides a fast and convenient method for analysis of urine samples. The objective of this study was to develop a comprehensive, fast and accurate method for compound characterization and effective molecular profiling of urine samples.

## Why Urine?

- Used since ancient times for medical diagnosis (Urology → Urinalysis)
- In ancient times considered "a divine fluid and window to the body"



Fig. 1: A Physician Examining Urine through a Matula (Left). A Urine Color Chart (Right).

- Urine-based measurements exhibit advantages over other biological fluids:

- Samples are easily collected in large quantities
- Contain relatively high concentrations of pharmaceuticals, illicit drugs, and steroids
- Compounds are typically detectable over extended periods of time

## Instrument: LECO Pegasus® GC-HRT

- High Quality Spectral Data
  - Comprehensive
  - Search Against Well-Established Databases (e.g., NIST, Wiley)
- Excellent Mass Accuracy Values (< 1 ppm) = Robust Formulas for Fragment, Molecular & Adduct Ions
- High Resolution Deconvolution™ (HRD™)
- High Resolving Power (up to 50,000)

## Compound Characterization Strategy




Fig. 2: The Utilization of EI- and CI-HRT Technology for Confident Compound Characterization.

## Experimental

### Samples

- Obtained from a collaborating European hospital
- 52 patient monitoring samples
- Sample preparation




Fig. 3: SPE Urine Sample Preparation.

### Acquisition Parameters

GC	Agilent 7890B
Column	Agilent J&W VF-DA (12m x 0.20 mm x 0.33µm)
Carrier Gas, Flow	He, 1.0 mL/min Constant Flow
Injection/Volume	Splitless/2 µL
Injector Temp.	280 °C
Temp. Program	70 °C (1 min) to 320OC at 25 °C/min (5 min)
MS	LECO Pegasus® GC-HRT
Transfer Line Temp.	300 °C
Ion Source Temp.	250 °C (CI 200 °C)
Ionization	El (70 eV); CI (140 eV)
Mass Range	45 – 520 (CI 60 – 800, Reagent Gas = 5% NH <sub>3</sub> in CH <sub>4</sub> )
Acquisition Rate	10 cps
Mass Calibration	PFTBA (Internal)

## Results: Illustrative Urine Samples:




Fig. 4: Illustrative Urine Samples A, B, C.

## Sample A (Traditional Drugs)

### Representative Compounds



Fig. 5: Analytical Ion Chromatogram (AIC) – Sample A.

Table 1: Compounds in Sample A.




Fig. 6: Peak True, Deconvoluted (Top) and Library Mass Spectra (Bottom) for Caffeine in Sample A.




Fig. 7: Peak True Mass Spectrum and Table of Accurate Mass Ions for Cholesterol in Sample A.

### Traditional Drugs




Fig. 9: Peak True Mass Spectra for Meconin (Top) and Lupanine (Bottom) in Sample A.




Fig. 10: Peak True Mass Spectra for Codeine (Top) and Morphine (Bottom) in Sample A.

## Sample B (Synthetic Drugs)

### Representative Compounds




Fig. 11: AIC and Table Listing Compounds in Sample B.

### Synthetic Drugs




Fig. 13: An XIC showing an unknown synthetic drug or metabolite not present in the NIST or Wiley libraries (A). Peak True mass spectrum for the unknown (B). The ChromaTOF-HRT® Spectra Table (C) showing the formula C<sub>23</sub>H<sub>21</sub>NO<sub>2</sub> for m/z = 357.13595 (mass accuracy = 0.02 ppm) and the initiation of a spectral database search resulting in the potential hit: 4-(3-(1-naphthyl)1H-indol-1-yl), a JWH 073 N-butanolic acid metabolite.

## Sample C (Amphetamine + Pharmaceuticals + Metabolites)




Fig. 14: AIC Displaying Representative Compounds in Sample C.

Table 2: Representative Compounds in Sample C.




Fig. 15: AIC and Table Listing Pharmaceuticals and Metabolites in Sample C.

## Summary

- GC-HRT analysis provides an effectual "molecular profile" of samples
- Confident compound identification through spectral similarity searches & robust formula determinations for fragment, molecular and adduct ions

LECO, Pegasus, and ChromaTOF-HRT are registered trademarks of LECO Corporation. High Resolution Deconvolution and HRD are trademarks of LECO Corporation.

Fig. 16: EI-HRT Peak True (Top) and Library Mass Spectra (Bottom) for Quetiapine (anti-psychotic). CI-HRT Peak True Spectrum for Quetiapine (Right).