

New Column Technologies for Better and Faster Separations in UHPLC

Xiaoli Wang, Wu Chen, Jason Link, William Barber
Agilent Technologies
March 14, 2011



Outline

- Review new column technologies
 - Properties
 - Theoretical considerations
- Applications of superficially porous particles
 - Fast separations
 - “Common” separations
 - High efficiency separations
- Practical considerations
 - Instrument impact
 - Method transfer
- Conclusions




Column technologies for small molecules

Particle	Particle size	Special need	Optimum applications
Traditional fully porous	2.5-5.0 μm	n/a	Widely used in many established methods
Sub-2- μm fully porous	1.5-2.0 μm	UHPLC; low extra-column broadening	Fast separation or high efficiency?
Sub-3- μm superficially porous	2.6-2.7 μm	low extra-column broadening	Fast separation or high efficiency?
Sub-2- μm superficially porous	1.7 μm	UHPLC; low extra-column broadening	Fast separation or high efficiency?
Monolith	n/a	n/a	High efficiency separations



Modern superficially porous particles



- Physical-chemical properties
 - Small particle size distribution: $d_{90}/d_{10} < 1.20$
 - Available particle size: 2.7, 2.6, 1.7 μm
 - Small particle porosity
- Chromatographic properties
 - Smaller h_{min} compared to fully porous
 - Smaller A term: better packed bed homogeneity
 - Smaller B term: less porosity thus less longitudinal diffusion
 - Similar C term for small molecules
- What are the optimum applications for superficially porous particles?

Gritti *et al.* *J Chrom A*, **2010**, 1217, 3819



Theory to optimize efficiency

- Three-parameter optimization (Knox-Saleem limit)
 - Simultaneously optimize particle size, column length and flow rate
 - Useful to assess the *performance potential* of different *particle types*

$$d_p^* = \left[\frac{\Phi \eta^{B/C}}{P_{\max}} \right]^{1/4} (\lambda t_0)^{1/4} D_m^{1/2} \quad L^* = \left[\frac{P_{\max}^{B/C}}{\Phi \eta} \right]^{1/4} (\lambda t_0)^{3/4} D_m^{1/2} \quad u_e^* = \left[\frac{P_{\max}^{B/C}}{\Phi \eta} \right]^{1/4} (\lambda t_0)^{-1/4} D_m^{1/2}$$

- Two-parameter optimization (Pope or kinetic plots)
 - Fix particle size, then optimize column length and flow rate
 - Useful to compare the *real performance* of different *particles*

$$L^* = \sqrt{\frac{P_{\max} \lambda t_0}{\Phi \eta}} d_p = \psi \lambda t_0^{1/2} d_p \quad u_e^* = \sqrt{\frac{P_{\max}}{\Phi \eta \lambda t_0}} d_p = \psi d_p t_0^{-1/2}$$

Carr et al. Anal. Chem., 2009, 81, 5342



Simple optimization procedure

Goal – To achieve the highest N in desired time

- Step 1
 - Conduct 3-parameter optimization to find optimal particle size, column length and flow rate
- Step 2
 - Choose the closest particle size available (**compromise**)
 - Conduct 2-parameter optimization to find optimal column length and flow rate
- Step 3
 - Chose the closest column length available (**compromise**)
 - Operate at the flow rate that gives the desired time

Wang et al. *LC-GC North America*, **2010**, 28, 932



Applications of superficially porous particles

- Ultrafast separations
- “Common” separations
- High efficiency separations



Ultrafast LC for dissolution test

	UV analysis	Traditional HPLC analysis
Easy to operate	Yes	Yes
Speed per sample	0.5-2 min	5-10 min
Speed per run	0.5 -2 hrs	3-5 hrs
On-line analysis	Readily available	Not readily available
Sample volume	8-10 mL	< 100 µL
Specificity	Moderate	Excellent
Sensitivity	Moderate	Excellent
Linearity range	Moderate	Excellent
Bubble problem	More likely to occur	Less likely to occur
Information	-	Rich info

Wang et al. LC-GC North America, 2010, 28, 932



Calculation for fully porous particles

- Optimization goal
 - $t_0 = 4$ sec or $t_R \sim 30$ sec
 - Use sub-2- μm fully porous particles at 1100 bar and 80 °C

Table V: Stepwise optimization for $t_0 = 4$ s separation on ZORBAX RRHD SB-C18 at 1100 bar and 80 °C*

Step	t_0 (s)	d_p (μm)	L (mm)	u_e (cm/s)	ΔP (bar)	N	% Loss
1	4	0.9	38	1.35	1100	18,500	—
2	4	1.8	77	2.70	1100	11,800	36
3	4	1.8	75	2.65	1057	11,700	37

* Conditions: 10% acetonitrile, $D_m = 1 \times 10^{-5} \text{ cm}^2/\text{s}$ for ciprofloxacin at 80 °C.


Wang et al. LC-GC North America, 2010, 28, 932



Ultrafast separation

Mobile phase A: 0.1% TFA in water, mobile phase B: ACN, **80 °C**, 0.3 µL injection, 278 nm

Zorbax RRHD SB-C18, 1.8 µm, 2.1x75 mm



Wang et al. LC-GC North America, 2010, 28, 932



Calculation for superficially porous particles

- Optimization goal
 - $t_0 = 4$ sec or $t_R \sim 30$ sec
 - Used sub-3- μm superficially porous particles at 550 bar and 80 °C

Table VI: Stepwise optimization for $t_0 = 4$ s separation on Poroshell 120 SB-C18 at 550 bar and 80 °C*

Step	t_0 (s)	d_p (μm)	L (mm)	u_e (cm/s)	ΔP (bar)	N	% Loss
1	4	1.2	37	1.33	550	15,500	—
2	4	2.7	82	2.90	550	9200	41
3	4	2.7	75	2.67	466	8900	43

* Conditions: 10% acetonitrile, $D_m = 1 \times 10^{-5} \text{ cm}^2/\text{s}$ for ciprofloxacin at 80 °C.


Wang et al. LC-GC North America, 2010, 28, 932



Ultrafast separation

Mobile phase A: 0.1% TFA in water, mobile phase B: ACN, **80 °C**, 0.3 µL injection, 278 nm


Poroshell 120 SB-C18, 2.7 µm, 2.1x75 mm



Wang et al. *LC-GC North America*, **2010**, 28, 932




Ultrafast separation for dissolution test




- The ultrafast LC methods provide equivalent dissolution results but is
 - 3-5 times faster than traditional HPLC methods
 - similar to automated UV methods in speed
- Such ultrafast separations can be useful for other quick assay applications, such as content uniformity and cleaning validation

Food and beverage additives



Eclipse Plus C18
1.8 μm , 3.0x100 mm




Poroshell 120 EC-C18
2.7 μm , 3.0x100 mm

Agilent Application Note 5990-6082en



Pharmaceutical screening

14 polar APIs
2.1x50 mm columns
0.8 mL/min, 3 min gradient
from 40% to 90% ACN



**Sub-3- μ m
superficially porous
Vendor 1**

**Sub-3- μ m
superficially porous
Vendor 2**


**Sub-3- μ m
superficially porous
Vendor 3**

**Sub-2- μ m fully porous
Vendor 4**


Fekete *et al. Talanta*, **2011**, 84, 416

High efficiency separation


Poroshell 120 SB-C18, 2.7 μ m, three 4.6x150 mm columns in series



N = 83,000
1 mL/min
316 bar



N = 103,000
1.5 mL/min
478 bar




N = 115,000
1.8 mL/min
573 bar

Agilent Application Note 5990-5602en



High efficiency separation

Poroshell 300 SB-C18, 5.0 μm , four 2.1x150 mm columns in series



Mobile phase A: 0.1% TFA in water with 50 mM NaClO_4 , mobile phase B: 0.1% TFA in 80/20 ACN/ H_2O , gradient from 0% to 42% B in 120 min, 0.50 mL/min, **70 °C**, 214 nm

Wang et al. *J Chrom A*, **2006**, 1107, 139




Practical considerations

- Instrument effect
- Method transfer from fully porous to superficially porous particles




Instrument broadening

Extra-column broadening on typical UHPLC instruments



Effect of extra-column broadening on 2-parameter optimum of 2.7 μm superficially porous at 600 bar



Gritti *et al.* *J Chrom A*, **2010**, 1217, 7677


Method transfer

- Consider the case where one wants to transfer an HPLC method to a faster method with 2.7 μm superficially porous particles, but keep the same resolution
- Example – Transfer USP Naproxen method (isocratic)

Method	Original HPLC Method	Quick Transfer #1	Quick Transfer #2
Particles	5 μm fully porous	2.7 μm superficially porous	2.7 μm superficially porous
L (mm)	150	75	100
F (mL/min)	1.2	2.4	2.6
Pressure (bar)	70	300	500



USP Naproxen method transfer



Conclusions

- Superficially porous particles have *great potential* to generate best performance, mainly due to lower $h_{\min}$
- Current superficially porous particles generate best performance *in a very wide analysis time range ($t_0 > 30$ sec or $t_R > 3$ min)*
- Various successful applications confirm the theoretical predictions
- Instrument performance can have strong impact on the performance of superficially porous particles, especially for fast separations



Acknowledgements

- Bill Long
- Anne Mack
- Maureen Joseph
- Ron Majors

