

Food Contact Material (FCM) Migration Study using HR-LCMS and Novel Software Database Suite



Charles Yang¹, Allan Brown³, Kate Comstock², Michal Godula⁴

^{1,2}Thermo Fisher Scientific, San Jose, CA, ³Scholle Packaging, IL, USA, ⁴Thermo Fisher Scientific, Prague, Czech Republic

Introduction

To ensure consumer safety, food/beverage packaging companies are required to conduct Food Contact Material (FCM) migration studies following regulatory guidelines^[1]. FCM migration study is the key component of marketing approval document "Food Contact Notification (FCN)" by authorities around the world.

FCM migration study requires components identification, characterization, and quantitation for safety assessment. Because of sample complexity and the unknown, unexpected nature of some components, advanced analytical instruments, combined with good software and databases, can significantly alleviate the challenge of FCM migration study.

This study presents a migration study workflow for food packaging bags, both non-gamma-irradiated and gamma-irradiated, using HR-LCMS and database search.

LCMS Analysis

Sample Preparation

Beverage bags non-irradiated and gamma irradiated at 10 kGy were extracted using 3% acetic acid and 50% ethanol food simulants for 12 days at 40°C [2].

Liquid Chromatography

LC separations were carried out on the Thermo Ultimate™ 3000 RS UHPLC system consisting of: DGP-3000RS pump, WPS-3000RS sampler, TCC-3000RS column compartment, and DAD-3000RS UV detector.

Column: Thermo Accucore C18, 2.1x100 mm 2.6 um. Column Temp: 35°C.

Injection Volume: 10 ul

LC Mobile phase: A: H₂O/0.1FA B: ACN/0.1% FA, flow rate: 400ul/min.

Gradient: Time (min.) 0 2 28 39 39.1
Mobile B (%) 5 15 90 90 5


Mass Spectrometry

MS analyses were carried out on Thermo Scientific™ Q Exactive™ Plus mass spectrometer using electrospray ionization (ESI). High resolution accurate mass (HRAM) full scan MS and top 3 data-dependent HCD MS/MS data were collected at resolving power of 70,000 and 17,500 (FWHM m/z200) with polarity switching. The scan range is 150-1500 amu. Stepped HCD normalized collision energy: 20, 40, and 60. Quantitation was conducted using SIM scan at isolation window 5 amu and resolution 70,000.

Results

The HRAM full scan data allow confident component identification and elemental composition assignment. The information-rich HCD MS/MS fragments provide valuable data for structure elucidation. Rapid positive/negative polarity switching gives additional information and confidence in component detection and characterization.

FIGURE 2. Base Peak Chromatograms of 50% EtOH Extraction Samples ESI (+)



Q Executive High Resolution Accurate Mass Measurement

Q Exactive High Resolution Accurate Mass (HRAM) data provides ultimate confidence for qualitative and quantitative analysis. The very sensitive, rapid polarity switching ensure detection of structurally diverse compounds at all levels.

SIEVE Software for Component Extraction and Differential Analysis

SIEVE, a differential analysis software, was used for component extraction. Figure 4 shows the base peaks alignment of extracted components. The extracted components are then searched against "ChemSpider" or Thermo E&L Compound Database to identify each component. A hit list is generated. A good database can improve hit rate and quickly identify the known compounds through HRAM data and ms/ms fragment ions matching.

Figure 1. HR-LCMS Analysis Workflow for FCM Migration Study




FIGURE 3. High Mass Accuracy for Positive/Negative Mode with Polarity Switching





FIGURE 4. SIEVE software for component extraction and Differential Analysis



Mass Frontier and MzCloud for Structure Elucidation

For unknown components, multiple possible structures were obtained for each component through ChemSpider database or other database searching, in order to determine the correct structures. "Mass Frontier™", a small molecule structure analysis software, was used. The "HighChem Fragmentation Library™" in Mass Frontier 7.0 has extensive published literature references. For each proposed structure, the "Fragments and Mechanisms" feature in Mass Frontier was used to generate predicted "fragments and mechanisms" through HighChem Library search. A high degree of correlation between predicted and experimental fragments confirmed the proposed structure. Mass Frontier then automatically annotated the matching fragments based on library search results, see Figure 5. Mass Frontier can build customized libraries.

FIGURE 5. Mass Frontier Auto Annotation for Fragment Ions



mzCloud Spectral Database Searching

A search was also conducted with "mzCloud", a high resolution spectral database. mzCloud provides several search criteria for small molecule structure identification using tandem mass spectra, including spectra, fragments, precursor ions, etc., all of which can be very useful for unknown structure elucidation. Figure 6 shows identification of Erucamide using the ms/ms spectrum search feature.

FIGURE 6. mzCloud Spectral Database Searching for Erucamide




FIGURE 7. Base Peak Chromatogram of 50% EtOH extraction for γ-Irradiated Bag





TABLE 1. Major Components Identified from 50% EtOH Extraction of γ-Irradiated Bag

Peak ID	RT (min.)	Measured	Calculated	Elemental Composition	Error (ppm)
P1	18.76	371.1006(H ₂ O)	371.1012	C ₁₂ H ₁₆ O ₂ S ₂	-1.3
		389.1112(H ₂ O)	389.1118	C ₁₂ H ₁₆ O ₂ S ₂	-1.9
P2	20.50	411.0001(H ₂ O)	411.0007	C ₁₂ H ₁₆ O ₂ S ₂	-1.7
		392.3006(H ₂ O)	392.3010	C ₁₂ H ₁₆ O ₂ N	-1.2
P4	22.13	445.1186(H ₂ O)	445.1190	C ₁₂ H ₁₆ O ₂ S ₂	-3.3
		463.1303(H ₂ O)	463.1306	C ₁₂ H ₁₆ O ₂ S ₂	-3.5
P5	22.90	485.1125(H ₂ O)	485.1125	C ₁₂ H ₁₆ O ₂ S ₂ Na	-0.1
		282.2787(H ₂ O)	282.2791	C ₁₂ H ₁₆ O ₂ N	-1.5
P6	25.19	513.1386(H ₂ O)	513.1388	C ₁₂ H ₁₆ O ₂ S ₂	-2.4
		537.1490(H ₂ O)	537.1494	C ₁₂ H ₁₆ O ₂ S ₂	-1.8
P7	25.55	564.1758(H ₂ O)	564.1759	C ₁₂ H ₁₆ O ₂ S ₂	-0.7
		473.2023(H ₂ O)	473.2026	C ₁₂ H ₁₆ O ₂ P	-3.6
P8	25.80	310.3099(H ₂ O)	310.3104	C ₁₂ H ₁₆ O ₂ N	-1.7
		700.4307(H ₂ O)	700.4300	C ₁₂ H ₁₆ N ₂ O ₂	-1.8
P10	26.30	336.3254(H ₂ O)	336.3261	C ₁₂ H ₁₆ O ₂	-2.0
		693.1667(H ₂ O)	693.1676	C ₁₂ H ₁₆ O ₂ S ₂	-1.6
P11	27.78	611.1672(H ₂ O)	611.1682	C ₁₂ H ₁₆ O ₂ S ₂	-1.6
		628.1841(H ₂ O)	628.1847	C ₁₂ H ₁₆ O ₂ S ₂	-1.0
P12	28.43	338.3419(H ₂ O)	338.3417	C ₁₂ H ₁₆ O ₂ N	-1.1
		338.3409(H ₂ O)	338.3417	C ₁₂ H ₁₆ O ₂ S ₂	-3.6
P14	30.01	667.1759(H ₂ O)	667.1759	C ₁₂ H ₁₆ O ₂ S ₂	-2.0
		685.1836(H ₂ O)	685.1836	Weak peak	-
P15	31.27	366.3725(H ₂ O)	366.3730	C ₁₂ H ₁₆ O ₂ N	-1.6
		741.1945(H ₂ O)	741.1952	C ₁₂ H ₁₆ O ₂ S ₂	-1.8
P17	38.75	776.2123(H ₂ O)	776.2123	C ₁₂ H ₁₆ O ₂ S ₂	-1.5
		693.4530(H ₂ O)	693.4537	C ₁₂ H ₁₆ O ₂ P	-1.3

FIGURE 8. Proposed Structures of Identified Compounds (Partial List)



High Resolution Quantitation of Antioxidant Irganox 1035

Quantitation is required for certain FCM. In this study, Antioxidant Irganox 1035 in food simulant 50% EtOH was used to demonstrate the Quan capability of Q Exactive Plus.

3,5-Bis(1,1-dimethylethyl)-4-hydroxybenzenepropanoic acid thiodi-2,1-ethanediyl ester

CAS# 41484-35-9

Molecular formula: C₃₈H₅₈O₆S

Formula weight: 642.3954

Stock solution: 1mg/mL in IPA. Working solutions were prepared by serial dilution of stock using food simulant 50%EtOH, see calibration curve for concentrations. The calibration curve was generated over the range of 0.05-100ppb, linear regression and 1/x weighting. The results show an excellent LOD 90.05ppb and linearity, see Figure 9 and Table 2. Negative mode data not shown

FIGURE 9. LCMS Chromatogram of IPA Reflux of Sample A




TABLE 2. Calibration Table for Antioxidant Irganox 1035

FileNames	Level	RT	Response	Specified Cal	Calculated C	% DB	% RSD	% CV	Integration Type	Exclude
1	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
2	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
3	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
4	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
5	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
6	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
7	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
8	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
9	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
10	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
11	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
12	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
13	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
14	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
15	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
16	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
17	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
18	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
19	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
20	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
21	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
22	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
23	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
24	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
25	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
26	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
27	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
28	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
29	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
30	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
31	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
32	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
33	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
34	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
35	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
36	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
37	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
38	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
39	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	
40	1	4.10	246.352	0.000	0.000	0.00	10.71	0.17	Manual Settings	

Conclusion

This study demonstrated an extractable analysis workflow for food contact material migration study for identification and quantitation. The UHPLC/HRAM full MS/HCD MS² with polarity switching on the fly data acquisition, coupled with novel database search, significantly increase the confidence and throughput of routine extractable analysis, in particular for unknown components identification and structure characterization.

References

- FDA CFR 21.94, CFR 600.11(b) and 600.11(h), CFR 211.160
- FDA: Guidance for Industry: Preparation of Premarket Submissions for Food Contact Substances; Chemistry Recommendations April 2002; December 2007

