


Analysis of Phosphopeptides by Automated IMAC in Gallium Form



Column: Thermo Scientific™ Acclaim™ 300
C18 3 μ m, 4.6 \times 150 mm

Pump: Thermo Scientific™ Dionex™ Summit
P680A DGP6

Mobile Phases: (A) 70:30 acetonitrile:water (v/v)
(B) water
(C) 0.40% TFA (v/v) +
1.0 mM $\text{Na}_4\text{P}_2\text{O}_7$

Gradient:	Times	0	15	25
%A	14	75	75	
%B	61	0	0	
%C	25	25	25	

Flow: 1.0 mL/min

Temperature: TCC-100 thermostat, 30 °C

Injection: ASI-100 autosampler; 30 μ L

Detector: UVD-340U; UV at 214nm

Note: See Thermo Scientific Dionex
Technical Note 705 for a complete
description of the automation and
IMAC conditions.

Samples: Fractionated on Thermo Scientific™
ProPac™ IMAC-10 2 \times 50 mm column
in gallium form:
(A) Unretained fraction
(B) Fraction retained by IMAC

Sample Prep: β -Casein, 2 mg/mL, reduced and
carboxymethylated, digestion with
immobilized trypsin, diluted 1:1 in
200 mM formic acid.

Peaks: 1. FQSpEEQQTEDELQDK
2. RELEELNVPGEIVSpLSpSpSpEESITR

24334

Immobilized metal affinity chromatography (IMAC) is a powerful tool for enriching phosphopeptide isolates from protein digests. Gallium and iron are typically used as the immobilizing metals in IMAC. An automation package consisting of the Thermo Scientific Dionex dual-gradient LC system, ProPac IMAC-10 (for enrichment) and Acclaim 300 analytical column is the perfect solution for this technique.