


Development and application of an exact mass LC-MS/MS library for the screening of mycotoxins and fungal metabolites in food and feed




Elisabeth Varga¹, Thomas Glauner², Emma Rennie³, Bernhard Wüst^{2,3}, Michael Sulyok¹, Rainer Schuhmacher¹, Rudolf Krška¹, Franz Berthiller¹

¹ University of Natural Resources and Life Sciences, Vienna (BOKU), Department IFA-Tulln, Center for Analytical Chemistry and Christian Doppler Laboratory for Mycotoxin Metabolism Konrad-Lorenz-Str. 20, 3430 Tulln, Austria

² Agilent Technologies, Chemical Analysis Group Hewlett-Packard-Str. 8, 76337 Waldbronn, Germany

³ Agilent Technologies, Life Sciences and Chemical Analysis 5301 Stevens Creek Blvd, Santa Clara, CA 95051, USA



University of Natural Resources and Life Sciences, Vienna

Department for Agrobiotechnology (IFA-Tulln)

Overview and Introduction

Purpose:

- creation of an exact mass HR-MS/MS library for mycotoxins and other fungal metabolites
- application for the screening and confirmation of mycotoxins in food and feed commodities

Methods:

- HR-MS/MS spectral library → flow injection of single analyte solutions
- classical screening approach and All-Ions MS/MS approach

Results:

- HR-MS/MS database (> 400) and library (~ 150) for mycotoxins and fungal metabolites
- All-Ions MS/MS yields fragments without precursor selection and is a suitable tool for screening

Mycotoxins are secondary **fungal metabolites** capable of causing **various toxic effects** including hepatotoxicity, mutagenicity, carcinogenicity or estrogenicity [1].


- produced by e.g. *Fusarium* spp., *Aspergillus* spp., *Penicillium* spp.



- detected in a **broad range of food and feed** e.g. cereals, beer, milk, spices, coffee, nuts, dried fruits [2]



- **400+ known** – many more unknown
- several are **regulated** in 100+ countries and standards are available for a limited number of substances [3]




Accurate mass screening for food contaminants is of growing interest due to the complexity of the samples and the increasing number of relevant analytes (masked and emerging mycotoxins).

Sample Preparation and LC-MS-Method

Due to the huge **differences in chemical and physicochemical properties** of mycotoxins, a **general sample preparation procedure** is required. Extraction was performed according to [4] and no clean-up was performed to avoid discrimination of certain mycotoxins.

Method Details:

LC-MS system: 1290 Infinity UHPLC 6550 iFunnel QTOF (Agilent Technologies), correction of mass axis with reference masses during the whole run
column: Zorbax SB-C18 RRHD (150 x 2.1 mm; 1.8 µm particle size)
eluent: water and methanol containing 0.1% formic acid
gradient: from 10% to 100% B in 19.5 min




Scheme 1: Applied analytical scheme

All-Ions MS/MS Approach

Alternatively the **"All-Ions MS/MS"** acquisition was applied for elimination of potential false positives, which used **fragmentation without precursor selection**. The presence of mycotoxins was confirmed by the **co-elution of characteristic fragment ions**


1) HR-MS/MS All-Ions scan

Measurement with "low energy channel" (no collision energy) and at least one "high energy channel" (collision energy e.g. 20 eV) (a)



2) Database search and library confirmation

- I) Database search: applying "Find by Formula" algorithm on low channel
 - II) Extraction of ion chromatograms (EICs) in high channel mode
 - III) Alignment EICs of fragments with parent EICs and correlate qualified by library confirmation
- The coelution score is the major function to qualify compounds based on similarity of peak shapes.
- (b) summary list; (c) fragment list; (d) overlaid EIC chromatograms; (e) co-elution plot; (f) MS-scan results, (g) All-Ions MS/MS results



Scheme 3: Example of an identification pathway using the All-Ions approach



Creation of Database & MS/MS-Library

- acquired **HR-MS/MS library spectra** for about 150 mycotoxins and other fungal metabolites → injection of individual standards

collision energies: typically 10, 20 and 40 eV
acquisition rate: MS: 5 spectra/s; MS/MS: 5 spectra/s
isolation width: narrow (1.3 amu)
ion species: pos. mode: [M+H]⁺, [M+Na]⁺; neg. mode: [M-H]⁻, [M+HCOO]⁻

- import **spectra and structure formulae** into user defined database and library
- **curation of fragment masses** based on fragment formulas and structures using **in silico fragmentation** (MassHunter Molecular Structure Correlator)
- only **explainable fragments** are included in the database with their **exact m/z values**


Figure 1: Library spectrum and structure of mycophenolic acid



Classic Screening Approach


1) HR-MS full scan

total ion chromatogram (TIC) of a naturally contaminated hazelnut sample (a)




2) Database search

extracted ion chromatogram (EIC) (b) and extracted MS spectrum (c) of mycophenolic acid




3) MS/MS scan

measurement principle of a targeted MS/MS scan (d)



4) Library confirmation

extracted MS/MS spectra at collision energies of 10 eV (e), 20 eV (f) and 40 eV (g) and library hit confirming the presence of the immunosuppressive mycotoxin mycophenolic acid (h)



Mycophenolic acid
C₁₇H₂₀O₆
produced by several *Penicillium* spp.
low acute toxicity, immunosuppressive

Scheme 2: Example of an identification pathway using a classic screening approach

Validation and Application

Validation:

- **45 analytes** were selected to cover a **wide range of different properties**: ionisation mode (neg./pos.), molecular mass (100 to 1300 amu), polarity, regulated mycotoxins and representatives of the most important groups.
- Neat standard solution and three matrices (**maize, hazelnut and wine**) were spiked on several concentration levels after extraction to cover three orders of magnitude.
- All samples were measured in negative and positive ionisation mode using both approaches. Preliminary results show, that **high library match scores** are achieved and false positive results are completely avoided. While the method is yet not sensitive enough to allow the detection of e.g. aflatoxins at the regulated levels, at a spiking level of 100 µg/kg the vast majority of the analytes were identified in the MS scans and confirmed by their MS/MS spectra. Thus, the method offers a possibility for the screening and **unambiguous confirmation** of a wide range of different mycotoxins in food and feed **even with a lack of reference standards**.

Application:

Additionally, the two approaches were applied to naturally contaminated samples. While in the initial TOF screen several contaminants were suspected, **applying the exact mass library** to the MS/MS data efficiently **eliminated false positives**.

Conclusion

Accurate mass screening using a UHPLC-QTOF in combination with an **exact mass library** is a **suitable technique for routine screening and confirmation** of a wide range of mycotoxins in food and feed. The application of the All-Ions approach reduces the measurement time (only one injection instead of two) and allows also **post-acquisition evaluation including MS/MS** information. The database currently comprises **more than 400 mycotoxins** and fungal metabolites. One third of the toxins are already included with their MS/MS spectra.

References

- [1] Bennett JW, Klich M (2003) Clin Microbiol Rev 16:497–516.
- [2] Zöllner P, Mayer-Helm B (2006) J Chromatogr A 1136:123–169.
- [3] Van Egmond HP, Schothorst RC, Jonker MA (2007) Anal Bioanal Chem 389: 147-157.
- [4] Sulyok M, Berthiller F, Krška R, Schuhmacher R (2006) Rapid Commun Mass Spectrom 20:2649–2659.



Pictures: *Fusarium* sp.: Marc Lemmens (BOKU); *Aspergillus* sp.: <http://www.eapri.eu/>; *Penicillium* sp.: <http://website.nbm-mnb.ca/mycologywebpages/Moulds/Penicillium.html>; Maize: Marc Lemmens (BOKU); Cheese: <http://www.yumsgar.com/Burning-Question-When-OK-Eat-Moldy-Food-4391316>; Pistachio: <http://www.buhlergroup.com/northamerica/en/process-technologies/optical-sorting/nut-sorting/pistachio-sorting.htm>; Chromatograms and spectra: Mass Hunter qualitative and quantitative analysis (Agilent Technologies)



Correspondence should be addressed to:
DI Elisabeth Varga, elisabeth.varga@boku.ac.at

University of Natural Resources and Life Sciences, Vienna
Department for Agrobiotechnology (IFA-Tulln), Center for Analytical Chemistry
Christian Doppler Laboratory for Mycotoxin Metabolism