

Tips and Tricks of Faster LC Analysis Without Capital Investment

It may be easier than you think!

Agilent Technologies, Inc.

Fall 2008



There Are Many Reasons to Reduce Analysis Time

- Run More Samples/Day
- Reduce Mobile Phase Use
- Reduce Mobile Phase Waste Disposal
- Better Utilize Existing Resources in Instruments and Personnel
- Implement “HB5”* Work Schedule

(* “HB5” **H**ome **B**y **5**, Ray Lombardi, Agilent Technologies, Inc. 2007)

At What Cost and Difficulty?

That Depends on What You Want to Accomplish!

How Fast Do You Want To Go?

Your Answer Will Determine the Cost of Increased Speed

- 2 to 3 X increase in Speed – Easy
- 5-10 X Increase – May Require Instrument Tweak
- 10 X + Increase – Will Require Instrument Upgrade and Tweak

What Are the Paths to LC Method Speed Gains?

- Most Speed From Column and Particle Size Optimization
- Instrument Optimization Will Enable or Boost Those Gains

The Following Discussion is Designed to


Simplify the Process - Often at No Cost!


How Do Small Particles Help Speed Up Methods?

- **Increase Kinetics of the separation**
- **Increase Efficiency (N) in Same Length Column**
- **Allow Shorter Column to Match Efficiency of Longer Column**
- **Provide Expanded Flow Rate Range Without Loss of Efficiency**

Trick: Use Smaller Particles For More Efficient Separations Across a Wide Flow Rate Range

Small Particle Columns including Monolith





ISOCRATIC ELUTION



Tip: Short Column/Smaller Particle Maintains Rs

Trick: Tailor Column Dimensions and Particle Size to Analysis Needs

Column Length (mm)	Resolving Power N(5 µm)	Resolving Power N(3.5 µm)	Resolving Power N(1.8 µm)	Typical Pressure Bar (1.8 µm)	Analysis Time*
150	12,500	21,000	32,500	>400	
100	8,500	14,000	24,000	>400	Analysis Time -33%
75	6000	10,500	17,000	320	Peak Volume -50%
50	4,200	7,000	12,000	210	Solvent Usage -67%
30	N.A.	4,200	6,500	126	
15	N.A.	2,100	2,500	55	Solvent Usage -90%

* Reduction in analysis time compared to 150 mm column


- pressure determined with 60:40 MeOH/water, 1ml/min, 4.6mm ID

Goal: Run Current Method 2X Faster

Trick: $\frac{1}{2}$ Column Length and Smaller Particle = $\frac{1}{2}$ Time

StableBond SB-C18


4.6 x 150 mm, 5 μ m



Rapid Resolution


StableBond SB-C18

4.6 x 75 mm, 3.5 μ m



Goal: Run Current Method 3X Faster


Trick: 1/3 Column Length and Smaller Particle = 1/3 Time



How Long Can the Column Be with MeOH?

System Pressure With MeOH/Water

4.6 x 50mm Column with 1.8um particles at 1 ml/min, 30°C




Percentage of MeOH (%)	Pressure (bar)
0	142
10	169
20	203
30	229
40	249
50	252
60	240
70	216
80	184
90	145
100	98

- Use P of 4.6 x 50mm as baseline
- Multiply P by L_{col2} / L_{col1}
- If P is Lower than Max P of Instrument you can run that Length

How Long Can the Column Be with ACN?

System Pressure with ACN/Water

4.6 x 50mm Column with 1.8um particles at 1 ml/min, 30°C




Percentage of ACN (%)	Pressure (bar)
0	142
10	154
20	159
30	156
40	150
50	135
60	120
70	104
80	88
90	70
100	59

- Use P of 4.6 x 50mm as baseline
- Multiply P by L_{col2} / L_{col1}
- If P is Lower than Max P of Instrument you can run that Length

More Resolution in Original Method than Needed

Tip: 1.8u, Very Short Column, Faster Flow 5-10X Faster


Isocratic on Instrument with 400bar Pressure Limit



Columns: ZORBAX SB-C18 Mobile Phase: 50% 20 mM NaH₂PO₄, pH 2.8: 50% ACN Flow Rate: 1 mL/min Temperature: RT
Detection: UV 230 nm Sample: 1. Estradiol 2. Ethynylestradiol 3. Dienestrol 4. Norethindrone

Tip: Constant Linear Velocity Maintains Efficiency

Trick: Reduce Flow Rate When Reducing Column Diameter



Quick Reference for Changing to Common Column Diameters

Maintains Equivalent Linear Velocity for Different Column ID's:

Column Type	Column ID	Flow Rate
Analytical	4.6 mm	1.0 mL/min
Solvent Saver	3.0 mm	0.42 mL/min
NarrowBore	2.1 mm	0.21 mL/min
MicroBore	1.0 mm	47 μ L/min
Capillary	0.5 mm	12 μ L/min
Capillary	0.3 mm	4.2 μ L/min
Nano	0.1 mm	472 nL/min
Nano	0.075 mm	266 nL/min

Flow rate column 2 = (diameter column 2)² / (diameter column 1)² x Flow rate column 1

- **Maintain equivalent mobile phase linear velocity when changing column diameter.**

Isocratic Tips and Tricks Summary

- Increased Isocratic Speed Does Not Necessarily Require Sub-Two Micron Particles
- Particle Diameter and Column Length are Major Factors in Resolving Power
- Chose New, Shorter Column to Approximate Original Column Efficiency - Chose Length Based on Time Savings Goal

Example 3X: 150mm, 5 μ m (~12,000N) ~ 50mm, 1.8 μ m (~12,000N)

- Usually Only Need to Reduce Col. Length and Particle Size – Not Necessary to Change Flow, Injection Vol. or Organic Content of M.P.
- Altered Column Diameter Requires Change in Flow and Injection Volume






GRADIENT ELUTION



Reducing Particle Size Improves Kinetics

Reduces Peak Width, Increases Height- It Does NOT Change the Gradient!



➤ Chromatograms Offset for Better View

Resolution Equation for Gradient Elution

Relationship of k^* to Gradient Time, Flow and Column Dimensions

$$R \approx \frac{\sqrt{N}}{4} \propto k^*$$


$$k^* \propto \frac{t_g F}{S (\Delta\%B) V_m}$$

$\Delta\%B$ = difference between initial and final % B
 S = constant (≈ 4 for 100 - 500 Da)

F = flow rate (mL/min.)
 t_g = gradient time (min.)
 V_m = column void volume (mL)



$$\text{Zorbax } V_m = \pi \times (\text{Col Internal Radius})^2 \times \text{Length} \times 0.6$$

Expected: Use of a Longer Gradient Time Increases Gradient Retention




• Increased gradient retention improves resolution of several peak pairs – 1,2 and 4,5.

Tip: A Shorter Column (smaller V_m) Also Increases Gradient Retention and R_s Increases Gradient Retention, Increases Overall Resolution



Column: ZORBAX SB-C8
Gradient: 20 – 60% B
Mobile Phase: A:H₂O with 0.1% TFA, pH 2
B: Acetonitrile
Flow Rate: 1.0 mL/min
Temperature: 35°C
Sample: Herbicides

Tip: Improve Resolution By Using Short Column Length (Vm) for BioMolecules



Mobile Phase:

A: 95% Water : 5 % ACN, 0.1% TFA
B: 5% Water : 95% ACN, 0.085% TFA
Gradient: 10-60% B in 30 min.

Flow Rate:

1.0 mL / min.

Temperature:

Ambient

Sample: 1. Gly-Tyr

2. Val-Tyr-Val

3. [Gln¹¹] Amyloid-β-
Protein Fragm 1-16

4. (TYR8) Bradykinin

5. Met-Enk

6. Leu-Enk

7. Angiotensin II

8. Kinetensin

9. RNase

10. Insulin (Eq.)



Tip: Maintain k^* To Keep Relative Peak Position in a Chromatogram Unchanged While Reducing Time

Any Decrease in

V_m-Column length

V_m-Column (i.d.)

Can be Offset by a Proportional

Decrease in t_G or F

Decrease in t_G or F

$$k^* \propto \frac{t_G \cdot F}{S \cdot \Delta\%B \cdot V_m}$$

Trick: Change V_m and t_G by Same Proportion

Two Chromatograms Both Having the Same Gradient Steepness


Sample: 1. Tebuthiuron 2. Prometon 3. Prometryne 4. Atrazine 5. Bentazon 6. Propazine 7. Propanil
8. Metolachlor

Column: **StableBond SB-C8**
4.6 x 150 mm, 5 μ m

Gradient
Time: **30 min.**

Flow Rate: 1.0 mL/min

**Analysis
Time: 24 min**




Column: **Rapid Resolution
StableBond SB-C8**
4.6 x 75 mm, 3.5 μ m

Gradient
Time: **15 min.**

Flow Rate: 1.0 mL/min

**Analysis
Time: 12 min**



001784S1.PPT

Very Fast LC on Conventional 1100 HPLC

G1379 Degasser

G1311 Quaternary pump

G1313A ALS autosampler

G1316A column compartment

G1314A VWD (standard cell G1314-60086, 10mm, 14uL)

The Instrument

Acetophenone

Diethyl phthalate

Benzophenone

Butyrophenone

Valerophenone

Hexanophenone


Heptanophenone

Octanophenone

The Sample



Conventional Column - 4.6 x 150mm, 5µm, SB-C18



Flow Rate 1.0 ml/min
Injection Volume 15µL
Temperature 30°C
Wavelength 246nm
Sample rate 2.5 Hz

Time (min)	% Acetonitrile
0	50
10	90
13.5	90
13.6	50
15	50


Initial Pressure: 69 bar

Final Pressure: 38 bar

Tip: Shorten Column and Gradient Time by Same Factor

1/3 Column Length- 1/3 Gradient Time

RRHT Column – 4.6 x **50mm**, 1.8 μ m, SB-C18



Flow Rate 1.0 ml/min
Injection Volume 5 μ L
Temperature 30°C
Wavelength 246nm
Sample rate **13.74 Hz**


Time (min)	% Acetonitrile
0	50
3.33	90
4.5	90
4.53	50
5	50

Initial Pressure: 132 bar

Final Pressure: 74 bar

The Comparison of chromatogram with the same RRHT column in 1100 and 1200 SL

Overlay the two Chromatograms




The red one is 1200, the blue one is 1100.

Tip: Increase Column Flow-Reduce Gradient Time

Double Flow (2mL/min) – ½ Gradient Time

RRHT 4.6 x 50mm, 1.8µm, SB-C18



Flow Rate **2.0 ml/min**
Injection Volume 5µL
Temperature 30°C
Wavelength 246nm
Sample rate **13.74 Hz**

Time (min)	% Acetonitrile
0	50
1.67	90
2.25	90
2.27	50
3.34	50

Initial Pressure: 266 bar
Final Pressure: 146 bar

Gradient Tips and Tricks Summary

$$k^* \propto \frac{t_g F}{S (\Delta\%B) V_m}$$

Gradient Retention Relationship Formula Is the Key to -

- Improving Resolution

- Longer Gradient Time (t_g)
- Smaller Column Volume (V_m) Length or Inner Diameter
- *Keeping All other Parameters Constant*

- Increasing Speed Without Losing Resolution

- Shorter Column with Proportionally Shorter Gradient Time
- Shorter Gradient Time with Proportionally Faster Flow Rate



Tip: Method Translator Makes Changes Easy

Trick: Let Agilent Method Translator Do the Math

Basic Mode for Easy Transfer of Conventional Method to RRLC

Original Method

System Info

- Agilent 1100/1200 Series Binary Pump
- Agilent 1100/1200 Series Quat. Pump

Column Info

Column ID [mm]: 4.6

Column Length [mm]: 150

Particle Size [µm]: 5.0

Method Info

Flow Rate [mL/min]: 1.50

Injection Vol. [µL]: 5.00

Pressure [bar]: 60.92

Solvent: Water / Acetonitrile

Temperature [°C]: 40.00

max. Solvent Visc. [cP]: 0.75

	Time [min]	%B	Flow [ml/min]
Initial	0.000	10.0	1.500
Gradient	30.000	90.0	1.500
Hold to	40.000	90.0	1.500
Return to	41.000	10.0	1.500
Stop	50.000	10.0	1.500

New Method

System Info

Agilent 1200 Series RRLC

Column Info

Column ID [mm]: 2.1

Column Length [mm]: 50

Particle Size [µm]: 1.8

Method Info

Flow Rate [mL/min]: 0.313

Injection Vol. [µL]: 0.30

Pressure [bar]: 156.70

Detector Settings: >0.03 min (0.5 s)

Time-Saving Factor: 3.0

fast ————— ultra-fast

- Simple Conversion
- Speed Optimized
- Resolution Optimized

	Time [min]	%B	Flow [ml/min]
Initial	0.000	10.0	0.313
Gradient	10.000	90.0	0.313
Hold to	13.333	90.0	0.313
Return to	13.667	10.0	0.313
Stop	16.667	10.0	0.313

- Injection Volume Conversion
- Detector Settings recommendation
- Gradient and Isocratic Method Conversion (auto-detected)

Agilent Method Translator – Advanced Mode

More Detailed Information, But Still Easy to Use

Detailed Input Detailed Output

Original Method

System Info

Max. System Pressure [bar] 400 600
 Allowed Pressure [%]
 Max. System Flow [mL/min]
 System Dispersion [µL]
 System Delay Volume [µL]

Column Info

Column ID [mm]
 Column Length [mm]
 Particle Size [µm]
 Flow Resistance Factor
 Porosity
 Red. Plate Height

Method Info

Max. Solvent Viscosity [cP]
 Flow Rate [mL/min]
 Injection Volume [µL]

Calculated Values

Pressure [bar]
 Linear Velocity [mm/s]
 Estimated k* (grad. k')
 Column Volume [mL]
 Plate Numbers
 Effective Plate Numbers
 Plate Number Yield [%]
 Peak Volume [µL]
 Peak Width [min]
 Peak Width [s]
 Peak Capacity

Method Info

Max. Solvent Viscosity [cP]
 Flow Rate [mL/min]
 Injection Volume [µL]

	Time [min]	%B	Flow [ml/min]
Initial	0.000	10.0	1.500
Gradient	30.000	90.0	1.500
Hold to	40.000	90.0	1.500
Return to	41.000	10.0	1.500
Stop	50.000	10.0	1.500

Comment

New Method

System Info

Max. System Pressure [bar] 400 600
 Allowed Pressure [%]
 Max. System Flow [mL/min]
 System Dispersion [µL]
 System Delay Volume [µL]

Column Info

Column ID [mm]
 Column Length [mm]
 Particle Size [µm]
 Flow Resistance Factor
 Porosity
 Red. Plate Height

Method Info

Max. Solvent Viscosity [cP]
 Flow Rate [mL/min]
 Injection Volume [µL]

Calculated Values

Pressure [bar]
 Linear Velocity [mm/s]
 Estimated k* (grad. k')
 Column Volume [mL]
 Plate Numbers
 Effective Plate Numbers
 Plate Number Yield [%]
 Peak Volume [µL]
 Peak Width [min]
 Peak Width [s]
 Peak Capacity

Time-Saving Factor:


Simple Conversion
 Speed Optimized
 Resolution Optimized
 Correct for System Delay Volume

	Time [min]	%B	Flow [ml/min]
Initial	0.000	10.0	0.313
Gradient	10.000	90.0	0.313
Hold to	13.333	90.0	0.313
Return to	13.667	10.0	0.313
Stop	16.667	10.0	0.313

Original Method → **New Method**

Does it work? - Example

Analysis of impurities of an active pharmaceutical ingredient by conventional HPLC (4.6mm ID x 250 mm, 5.0 μm):



Converting to a 4.6 mm ID x 100 mm, 1.8 μ m column:

Agilent Method Translator 1.0

Basic Mode | **Advanced Mode** | Viscosity Table | Help

Original Method

System Info

Agilent 1100/1200 Series Binary Pump
 Agilent 1100/1200 Series Quat. Pump

Column Info

Column ID [mm]
 Column Length [mm]
 Particle Size [μ m]

Method Info

Flow Rate [mL/min] Solvent
 Injection Vol. [μ L] Temperature [$^{\circ}$ C]
 Pressure [bar] max. Solvent Visc. [cP]

	Time [min]	%B	Flow [ml/min]
Initial	0.000	5.0	1.400
Gradient	20.000	90.0	1.400
Hold to	23.000	90.0	1.400
Return to	23.010	5.0	1.400
Stop	30.000	5.0	1.400

Comment

New Method

System Info

Agilent 1200 Series RRLC

Column Info

Column ID [mm]
 Column Length [mm]
 Particle Size [μ m]

Method Info

Flow Rate [mL/min]
 Injection Vol. [μ L]
 Pressure [bar]
 Detector Settings



Time-Saving Factor:

4.6



fast ultra-fast

- Simple Conversion
- Speed Optimized
- Resolution Optimized


	Time [min]	%B	Flow [ml/min]
Initial	0.000	5.0	2.585
Gradient	4.333	90.0	2.585
Hold to	4.983	90.0	2.585
Return to	4.986	5.0	2.585
Stop	6.500	5.0	2.585

Does it work? Yes!



Does it work?




4.6 mmID x 100 mm, **1.8 μ m** Zorbax SB C18

Speed Optimized

0.00 min	5%	B
4.33 min	90%	B
4.98 min	90%	B
4.99 min	5%	B
6.5 min	5%	B

ZORBAX Selectivity Choices

Options for Improved Speed and Resolution





Tip: Bonded Phase Can Change Selectivity Improve Peak Shape and Shorten HPLC Assay Time

- Each ZORBAX phase available in matching 1.8 μ m, 3.5 μ m and 5 μ m choices, most in 7 μ m!
- Over 140 RRHT, 1.8 μ m 600 bar column choices available
- 14 Different Column Chemistries – 13 bonded phases and silica for HILIC use
- 7 column lengths (250*, 150, 100, 75, 50, 30 and 20 mm long) * 1.8 μ m as custom column
- 3 internal diameters (4.6, 3.0, 2.1 mm and Prep ID)



Tip: C8 Often Yields Same Peak Order in < Time Than C18

Trick: Evaluate Different Phases with Same M.P.




Less retention can save significant time – the C8 is a good choice here.

The RRHT column is delivering the efficiency and resolution expected, but the C8 bonded phase may be the best choice.

Tip: Polar Phases Can Save Analysis Time

Trick: Evaluate Non-Polar and Polar Phases with Same M.P.




Mobile Phase 40 % ACN 60 % 25 mM Sodium Phosphate Buffer pH 2.4 Flow Rate: 1.5 ml/min 4.6 x 50mm UV 210 nm
2 μ l Elution order for Eclipse Plus Phenyl Hexyl: (1) Piroxicam, (2) Sulindac, (3) Tolmetin, (4) Naproxen, (5) Ibuprofen, (6) Diclofenac, (7) Celebrex (equal portions of approximately 1 mg/ml solutions)

Tip: Higher Temperature Can Improve Speed and Resolution






Higher Temperature should always be considered as a parameter during speed optimization

- Provides more rapid mass transfer:
 - Improves Efficiency – **enhances resolution**
 - Decreases analysis time – **faster separations** with no loss in resolution
- Decreases Mobile Phase Viscosity
 - Lowers backpressure – allows for higher flow rates, **faster separations**, greater efficiency
- Can change selectivity – **optimize resolution**
- Not Necessary to Go to the Limit-Midrange Offers Benefits



Tip: Temperature Changes Can Alter Selectivity

Trick: Optimize Col Temp For Best Speed/Rs





Tips on When You Should Move to the Higher Pressure and Performance 1200 SL (RRLC) System

- Very Fast Methods with High Velocity Mobile Phase That Increase Backpressure Above 400bar
- Resolution of Difficult Samples Requires Long Columns With 1.8 μ m Particles (>400bar Backpressure)
- Labs Required to Run Both Traditional and New, Higher Performance LC Methods
- Need Faster Data Acquisition for Narrow, High Efficiency Peaks Generated in Higher Speed Methods

Tip: 1.8µm Particle Increases Resolution (How?)

Trick: Substitute 1.8µm Particles in Same Length Column



Conditions for both experiments

Pumps

- Solvent A: H₂O + 0.1% TFA
- Solvent B: ACN + 0.1% TFA
- Gradient: 10% to 95% ACN in 40min, hold for 1min
- Flow Rate: 0.4ml/min

Autosamplers

- Injection volume: 3µl

Thermostatted Column Comp.





- Temperature: 50°C

Detectors

- DAD 2µl cell and **20Hz**, 220nm,


Ref: Appl. Note 5989-4506 by Edgar Naegele

Tip: Small Particles Improve Detection of Low Level Impurities



Particles Reveal More Information and Improve Detection and Integration

Customer Sample, Translation of Isocratic Impurity Methods, Zoom Critical Time Range (t = 7min)




Up to 60% higher resolution
without loss in sensitivity

Example of What Is Possible!

Shorter Column, Smaller Particle, > Flow Rate, High Temp


F = 1.20ml/min
T = 40°C
Analysis Time = 11min
Solvent Cons. = 13.2ml



Easy!

High Resolution:
4.6mm x 150mm 5.0µm


F = 4.80ml/min
T = 40°C
Analysis Time = 1.05min
Solvent Cons. = 5.1ml



Easy!

High Speed:
4.6mm x 50mm 5.0µm


F = 1.00ml/min
T = 40°C
Analysis Time = 1.1min
Solvent Cons. = 1.1ml



Easy!

High Speed & Resolution:
2.1mm x 50mm 1.8µm

F = 2.40ml/min
T = 95°C
Analysis Time: 0.4min
Solvent Cons. = 1.0ml




Easy, But You Will Need > 400bar!

Max Speed at T = 95°C
2.1mm x 50mm 1.8µm

Tip: Higher Data Improves Response and Rs

Trick: Optimize Data Rate for Accurate Profile



80Hz versus 10Hz Data Rate

- Peak Width: **- 55%**
- Resolution: **+ 90%**
- Peak Capacity: **+ 120%**
- App. Column Eff.: **+ 260%**

Data Rate	Peak Width	Resolution	Peak Capacity
80 Hz	0.300	2.25	60
40 Hz	0.329	2.05	55
20 Hz	0.416	1.71	45
10 Hz	0.666	1.17	29
5 Hz	1.236	0.67	16

Sample: Phenones Test Mix
 Column: Zorbax SB-C18, 4.6x30, 1.8um
 Gradient: 50-100%ACN in 0.3min
 Flow Rate: 5ml/min

Moving From Conventional LC to RRLC:

*“Doesn't this
cost a fortune ?”*

Moving to RRLC is simple and easy !


- Same software (e.g. ChemStation, EZChrom)
 - ⇒ **No** additional training effort
- Limit capital investment by modular step-wise upgrade
 - ⇒ **Leverage** existing 1100 modules, existing method and stationary phase
- Configurable delay volume in 1200 SL pump
 - ⇒ **Easy** to also run legacy LC methods
- Simple LC to RRLC method translator
 - ⇒ **No** tedious method development



*“No - not
at all !”*

Stepwise Scale-up to Rapid Resolution LC

From 1100 to 1200 RRLC in two steps -
Add What You Need When You Need It!



Analysis Time	> 5 min	> 1.5min	> 0.2min
Cycle times	> 6 min	> 2min	> 0.4min
Peak Width	> 3 sec	> 1.5 sec	> 0.2 sec
N	5 - 12,000	5 - 30,000	5 - 60,000
Column ID	4.6 mm	2.1 - 4.6 mm	2.1 - 4.6 mm
Column Length	50 mm	50 - 150 mm	20 - 150mm
Flow rates	0.2 - 10ml/min	0.05 - 5 ml/min	0.05 - 5 ml/min
Temperature	80 C	80 C	100 C
Pressure	400bar	600 bar	600bar

SUMMARY

- Increasing Speed of HPLC Separations Need Not be Difficult or Costly
- Changes in Columns and Equipment Depend on Your Goals
- Majority of Time Saved Will Be Based on Column and Particle Size
- Isocratic and Gradient 2X , 3X and 5X Times Savings are Easy
- Greater Time Savings May Require Instrument Modifications
- 10-20X Increase in Speed May Require Instrument Upgrade
- Long Columns with 1.8 μ m Particles May Require >400 bar Pump and Autosampler Capability

