


Aggregation of Antibody-Drug Conjugates: The Light Scattering Toolbox for Screening and Characterization


Chris Broomell and Daniel Some, Wyatt Technology Corporation, 6300 Hollister Ave, Santa Barbara, CA 93117

Abstract



Antibody-Drug Conjugates (ADCs) show much promise as effective therapeutics for cancer and other diseases. However, they often exhibit an increased aggregation propensity compared to their unmodified counterparts due to non-specific interactions arising from attached drug and linker moieties. Light scattering offers multiple techniques for addressing the challenges of formulation screening, and characterizing both aggregates and propensity for aggregation. We present these tools in the context of ADCs and the dependence of aggregation profiles on linker chemistry.

High-Throughput Dynamic Light Scattering





Dynamic light scattering (DLS) analyzes Brownian motion via the intensity fluctuations of light scattered by macromolecules and nanoparticles in solution. DLS determines directly diffusion coefficients D_t . Particle radius R_h is calculated from the Stokes-Einstein equation:

$$R_h = \frac{kT}{6\pi\eta D_t}$$

DLS determines size distributions without fractionation, providing polydispersity estimates. High-throughput screening via the DynaPro® Plate Reader II is a powerful tool for characterizing multiple aspects of biotherapeutic stability:


- Degree and type of aggregation
- Thermal conformational stability (T_M , T_{onset})
- Chemical conformational stability
- Colloidal stability



Conclusions


- High-throughput screening via the DLS Plate Reader is essential for assessing stability and aggregation of ADCs, and the impacts of linker chemistry on the onset of aggregation.
- FFF-MALS yields high-resolution characterization of aggregates, leading to a better understanding of the aggregation process.
- In a separate study, MALS-UV-RI analysis determined drug-antibody ratio (DAR).
- The light scattering toolkit contributes to complete characterization of ADCs, leading to an efficient and rapid development timeline.

Linker-Induced Instability Studied by DLS



Protein stability is typically studied by temperature of melting (T_M) or onset of aggregation. ADC₂ unfolds and aggregates at a lower temperature, and is therefore less stable than ADC₁. These data suggest distinct aggregation pathways.

Identical mAb and drug, different linkers



Size Distributions by FFF-MALS

Multi-angle, static light scattering (MALS) analyzes molar mass M_w from first principles, via the average intensity of scattered light I and concentration c :


$$I_{scattered} \propto M_w c$$

This analysis is independent of questionable reference standards and column calibration.



ADC drug and linkers molecules are often hydrophobic, making them incompatible with fractionation via size-exclusion chromatography (SEC). Field flow fractionation (FFF) separates macromolecules and nanoparticles by size without a stationary phase, eliminating most of the non-ideal surface interactions prevalent in SEC.


Coupling an Eclipse™ FFF device to a Wyatt DAWN® HELEOS® II MALS detector creates a versatile system for accurate and robust characterization of molar mass and size distributions from 1 nm to 1000 nm.




In an Asymmetric-Flow FFF separation channel, macromolecules and nanoparticles are gently pushed against a semipermeable membrane by crossflow. Smaller particles diffuse back up towards the center of the channel. Laminar channel flow induces a parabolic flow velocity profile, causing smaller particles to elute earlier.

Aggregation Analyses of Stressed and Unstressed ADCs by FFF-MALS

Typical Fractogram – Unstressed Sample




Typical Fractogram – Stressed Sample



* Small shoulder suggestive of fragmentation. This is present to some extent in all samples.

Average M_w and Mass Fraction ($n=2$)



	Peak1		Peak 2		Peak 3	
	%	M_w	%	M_w	%	M_w
Unstressed	0.75	148.4	97.75	231.5	1.5	
Stressed	1.40	149.6	95.1	297.1	3.5	

M_w is determined by the ratio of light scattering signals to UV. Mass fraction is determined by integrating UV peaks. FFF-MALS is highly sensitive to small quantities of aggregates due to the dependence of scattered intensity on M_w .