

Evaluation of Amine Solutions for CO₂ Absorption by TOC/TN

M. Tanaka

User Benefits

- ◆ Separation and recovery processes for the greenhouse effect gas (GHG) CO₂ by amine solutions can be evaluated by using the TOC-L analyzer.
- ◆ The concentration of inorganic carbon, such as the concentrations of dissolved CO₂ and hydrogen carbonate ions, can be obtained by inorganic carbon (IC) measurement of amine solutions.
- ◆ Concentration control of amine solutions is possible by TOC/TN measurement.

Introduction

Global warming, which has become an issue in recent years, is caused by the increasing concentration of greenhouse effect gases (GHG) in the atmosphere. In particular, it is feared that the effects of CO₂ may affect ecosystems at the global scale. As a countermeasure for this problem, research on technologies such as CCS (carbon capture and storage) and CCUS (carbon capture, utilization, and storage) for separation and recovery of CO₂ and reuse of CO₂ by conversion to valuable substances is underway with the aim of realizing carbon neutrality.


CO₂ recovery techniques include chemical absorption methods, physical absorption methods, and membrane separation methods. The optimum method is selected based on the physical properties of the discharge source, such as the CO₂ concentration and pressure. Among these various approaches, chemical absorption methods utilizing amine-based solutions, which bond strongly with CO₂ by chemical reaction, have the advantage of applicability to low partial pressures such as combustion off-gas and other sources of CO₂ because of their high reactivity, and various research efforts are in progress.

This article introduces an example of measurement of CO₂ capture by amine-based solutions and TN measurement and TOC measurement of the amine solutions by adding a TNM-L total nitrogen measurement unit to Shimadzu TOC-L total organic carbon analyzer.

Measurement Method

Two types of amine solutions were prepared using 2-amino-2-methyl-1-propanol [abbreviation: AMP, chemical formula: (CH₃)₂C(NH₂)CH₂OH] or monoethanolamine [abbreviation: MEA, chemical formula: HOCH₂CH₂NH₂]. The AMP and MEA were diluted with pure water to prepare solutions having concentrations of 20 wt% of AMP and 30 wt% of MEA.

It is known that there are two paths in the reaction of amine and CO₂, as shown in ① and ② below (for the case of Class 1 amine R-NH₂).



MEA is widely used because it is inexpensive and has high reactivity with CO₂. However, since reaction ① is the main reaction, the stoichiometric ratio of the amine molecule and CO₂ molecule is 2:1. As a result, the heat of reaction with CO₂ is large, and the energy consumption required in CO₂ recovery is also large.

On the other hand, when using AMP, reaction ① is difficult because a bulky substituent is bonded around the amine group, and reaction ② becomes the main reaction, but because the stoichiometric ratio of the amine molecule and CO₂ molecule is 1:1, energy consumption for CO₂ recovery can be reduced. In addition, AMP also has the advantage of a large CO₂ absorption rate.

In this experiment, samples were prepared by flowing CO₂ gas through these two types of amine solutions, as shown in Fig. 1, to cause CO₂ absorption.




Fig. 1 Experimental Setup for CO₂ Absorption in Amine Solutions

First, 20.6 mL of the amine solution was introduced into an impinger, and CO₂ gas was flowed for approximately 2 h at a flowrate of 100 mL/min from a 4.74 vol% CO₂ (in N₂) cylinder to absorb CO₂ in the solution. Due to reaction ②, part of the CO₂ captured by the amine solution exists in the form of hydrogen carbonate ions.

These samples were evaluated by using a TOC-L analyzer, which has a function that enables separate quantitation of the total organic carbon (TOC) and inorganic carbon (IC) in the sample. IC measurements are carried out by the acidification and sparging method, in which the sample is first acidified to convert HCO₃⁻ and CO₃²⁻ to the form of dissolved CO₂, and this CO₂ is then extracted (sparged) with a clean gas and quantified using an infrared CO₂ detector. This method makes it possible to determine the carbon content of HCO₃⁻ unaffected by carbon originating from the amine, even though the CO₂ capture solution contains a large amount of amine. If the optional TNM-L total nitrogen measurement unit is added to the system, it is also possible to measure the total nitrogen (TN) content, which can be used to control the concentration of the amine solution.

By using the TOC-L analyzer, it was possible to obtain the quantity of CO₂ captured by the amine solutions based on the IC measurements and the TN concentration and TOC concentration of the amine solutions from the TN and TOC measurements.

■ Measurement Conditions

The IC, TN, and TOC concentrations were measured before and after absorption of CO₂ by the amine solution. The samples for the IC measurements were diluted 100x and the samples for the TN and TOC measurements were diluted 500x with pure water. Table 1 shows the measurement conditions.

Table 1 Measurement Conditions

Instruments	: TOC-L _{CPH} total organic carbon analyzer, TNM-L total nitrogen measurement unit
Catalyst	: TOC/TN catalyst
Measurement items	: IC, TN, TOC (= TOC after acidification and sparging)
Calibration curves	: IC: One-point calibration curve by 200 mgC/L of sodium carbonate/sodium bicarbonate (aq) TN: One-point calibration curve by 200 mgN/L of potassium nitrate (aq) TOC: One-point calibration curve by 400 mgC/L of potassium biphthalate (aq)
Injection volume	: 40 μL
Samples	: AMP (2-amino-2-methyl-1-propanol), concentration 20 wt% solution MEA (monoethanolamine), concentration 30 wt% solution
CO ₂ absorption conditions	: Absorption of 4.74 Vol% CO ₂ (in N ₂) by approximately 20 mL of sample at flowrate of 100 mL/min for approximately 2 h
Sample dilution rates	: IC measurements: 100x TN measurements: 500x TOC measurements: 500x

■ Measurement Results

Table 2 shows the results of the measurements of the IC, TN, and TOC of the amine solution before and after CO₂ gas absorption. All measurement values have been corrected for the dilution rate.

Table 2 Measurement Results

Sample	IC measurements values (%C)	
	AMP solution	MEA solution
Before CO ₂ absorption	0.00118	0.0006
After CO ₂ absorption	1.26	1.07

Sample	TN measurements values (%N)	
	AMP solution	MEA solution
Before CO ₂ absorption	3.37	7.18
After CO ₂ absorption	3.34	7.08

Sample	TOC measurements values (%C)	
	AMP solution	MEA solution
Before CO ₂ absorption	11.7	12.2
After CO ₂ absorption	11.4	12.1

From Table 2, it can be understood that the IC concentrations of both the AMP solution and the MEA solution increased greatly, by more than 1000 times, as a result of absorption of CO₂ gas. From this, it is clear that the IC concentration increases because the amine solutions capture CO₂ gas and this CO₂ is dissolved in the solutions as hydrogen carbonate ions.

On the other hand, no large changes were detected in the TOC and TN concentrations, indicating that there was no significant effect on the concentration of the amine solutions in this CO₂ gas capture experiment.

■ Conclusion

This experiment demonstrated that the amount of CO₂ captured by amine solutions can be quantified and the concentration of the amine solutions can be controlled by using the TOC-L analyzer. At present, various research projects on low cost, energy saving CO₂ separation and recovery technologies are underway with the aim of realizing carbon neutrality. Among those efforts, chemical adsorption methods using amine solutions have attracted considerable attention, and the TOC analyzer is expected to be useful in that research.



Fig. 2 Shimadzu TOC-L Total Organic Carbon Analyzer and TNM-L Total Nitrogen Measurement Unit