

# Application

Note  
141

## Air Sampling of VOCs by SPME for Analysis by Capillary GC

The disposable SPME portable field sampler and the Carboxen/PDMS fiber are an ideal combination for monitoring VOCs in air samples. The ability to collect air samples in the lab or field with the sampler for analysis at a later time is convenient, simple, and reliable for compounds at trace levels.

### Key Words

- SPME • portable field sampler
- volatile organic compounds • air analysis

Solid phase microextraction (SPME) can be used efficiently in the field to extract volatile organic compounds (VOCs) from air samples and to store the samples for testing in the lab at a later time. Our SPME portable field sampler, which includes an SPME fiber coated with 75µm of Carboxen™/polydimethylsiloxane (PDMS) material, has a unique sealing mechanism that eliminates the need to ship bag samples or liter containers from the field to the lab.

With the SPME portable field sampler, volatile and semivolatile compounds extracted from air or water remain on the SPME fiber until they are thermally desorbed. The highly retentive Carboxen/PDMS coating ensures that the compounds remain on the fiber for up to 3 days after sampling, at -4°C. After sampling, the fiber is retracted into the protective outer needle, the outer needle is drawn within a replaceable sealing septum in the nosepiece, and the unit is transported for analysis. The lightweight disposable holder can be used for the life of the fiber.

We extracted US EPA Method 624 VOCs from a water sample, stored the analytes for 3 days under various conditions, and then analyzed the samples. We compared the response values with those obtained in an extraction in which we performed *immediate* desorption (Table 1). With the Carboxen/PDMS fiber, stored at -4°C, no loss of sample was observed. The slight differences were within expected deviations between two extractions. Results of fibers stored at ambient temperature were nearly identical to the fiber stored at -4°C. Even when the fiber was left unsealed for 3 days at ambient temperature, the average loss was less than 15%. The Carboxen/PDMS fiber proved to be much more suitable for storing VOCs than the 100µm PDMS fiber, which typically has been used for VOC extractions.

Using a mixture of 9 VOCs, we attempted to determine the capacity of the Carboxen/PDMS fiber, and to determine whether analytes with low distribution constants would be displaced by analytes with higher distribution constants. The analytes were obtained from a 1 ppm gas stream and spiked into a 125mL bulb with a final

**Table 1. Percent Change in Analyte Response After 3 Days Storage**

Method 624 Analytes and Naphthalene at 20ppb	75µm Carboxen/PDMS		100µm
	-4°C, sealed with septa	Ambient, not sealed	-4°C, sealed with septa
<b>PDMS</b>			
<b>Analyte</b>			
Chloromethane	2.9	-22.7	-100.0
Vinyl chloride	-13.9	22.6	-100.0
Bromomethane	0.0	-44.8	-100.0
Chloroethane	-4.7	-10.8	-100.0
Trichlorofluoromethane	1.8	-29.0	-100.0
1,1-Dichloroethene	-1.8	-12.4	-100.0
Methylene chloride	5.1	-17.8	-100.0
trans-1,2-Dichloroethene	5.3	-7.7	-100.0
1,1-Dichloroethane	6.6	-4.8	-100.0
Chloroform	0.0	-2.9	-100.0
1,1,1-Trichloroethane	-0.8	-23.7	-100.0
Carbon tetrachloride	1.9	-18.5	-100.0
Benzene	1.5	-14.1	-99.6
1,2-Dichloroethane	3.3	-9.7	-99.7
Trichloroethene	-0.6	-3.7	-98.7
1,2-Dichloropropane	2.2	-17.1	-94.6
Bromodichloromethane	6.5	-7.7	-97.7
cis-1,3-Dichloropropene	7.0	-28.8	-99.2
Toluene	-2.0	-5.0	-95.8
trans-1,3-Dichloropropene	-1.0	-40.9	-94.8
1,1,2-Trichloroethane	1.8	-33.5	-96.3
Tetrachloroethene	-2.4	-18.5	-92.1
Dibromochloromethane	-4.1	-17.3	-95.2
Chlorobenzene	-0.1	-41.3	-88.1
Ethylbenzene	-1.7	-16.7	-82.7
Bromoform	-1.6	20.5	-83.7
1,1,2,2-Tetrachloroethane	0.0	-16.9	-69.6
1,3-Dichlorobenzene	-0.4	-1.3	-51.1
1,4-Dichlorobenzene	-0.2	-14.5	-49.9
1,2-Dichlorobenzene	0.1	-1.4	-37.6
Naphthalene	-4.8	-0.2	-14.2
<b>Mean Change</b>	<b>0.2</b>	<b>-14.2</b>	<b>-88.4</b>


concentration ranging from 400ppt to 400ppb. Using an SPME portable field sampler containing a Carboxen/PDMS fiber, the analytes were extracted for 10 minutes at ambient conditions. Figure A shows the analytes at 10ppb.

Figure B shows the concentration plots for three analytes. Excellent linearity is indicated by the correlation coefficients of 0.98 or higher. Of the 9 VOCs, only 1,2-dichloroethane yielded a value less than 0.98 (Table 2). Vinyl chloride and 1,3-butadiene had the lowest distribution constants, while tetrachloroethene and trichloroethene had the highest.


Classical adsorption mechanism theory suggests that, due to limited adsorption sites, the lighter analytes will be displaced by the heavier analytes as concentration increases. At concentrations

**Figure A. VOCs in Air at 10ppb**

Sample: VOCs in air at 10ppb, 125mL bulb  
 SPME Fiber: Carboxen/PDMS, 75µm film  
 Cat. No.: 57318  
 Extraction: headspace, 10 min  
 Desorption: 5 min, 300°C  
 Column: SPB™-1 SULFUR, 30m x 0.32mm ID, 4.0µm film  
 Cat. No.: 24158  
 Oven: 40°C (2 min) to 150°C at 8°C/min  
 Carrier: helium, 35cm/sec  
 Inj.: splitless (closed 2 min), 0.75mm ID liner  
 Det.: GC/MS ion trap, m/z = 45 - 260



**Figure B. Linear Analyte Response Using Carboxen/PDMS Fiber**



Contact our Technical Service Department (phone 800-359-3041 or 814-359-3041, FAX 800-359-3044 or 814-359-5468) for expert answers to your questions.

**Table 2. Correlation Coefficients for VOCs**

Analyte	R <sup>2</sup>
Vinyl chloride	0.980
1,3-Butadiene	0.990
Acrylonitrile	0.980
Methylene chloride	0.986
Chloroform	0.984
1,2-Dichloroethane	0.953
Benzene	0.995
Trichloroethene	0.982
Tetrachloroethene	0.990

above 400ppt, the amount of analyte extracted would level off. However, the unique pore structure of Carboxen-1006 enables extraction of all analytes without displacement of the lighter analytes.

The portable field sampler and the Carboxen/PDMS fiber are an ideal combination for monitoring air samples. The ability to collect both air and water samples in the lab or field with the sampler and analyze at a later time is convenient, simple, and reliable for trace level analyses.

**Ordering Information:**

Description	Cat. No.
SPME Portable Field Sampler, pk. of 2 Permanent 75µm Carboxen/PDMS fiber**	504831
SPME Fiber Assembly, pk. of 3 75µm Carboxen/PDMS	57318
SPME Holder for Manual Sampling Order fiber assemblies separately.	57330-U
Thermogreen™ LB-2 Septa, 5mm pk. of 100	20638
SPME Septum Tool	504858
SPB-1 SULFUR Capillary GC Column 30m x 0.32mm ID, 4.0µm film	24158

\*Solid phase microextraction technology licensed exclusively to Supelco. US patent 5,691,206; European patent #0523092.

\*\*Note: Our original SPME unit and replaceable Carboxen/PDMS fibers can be used for this application if the analytes are analyzed immediately after the extraction.

Carboxen, SPB, and Thermogreen are trademarks of Sigma-Aldrich Co.

Note 141

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien  
 BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8  
 CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki  
 FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens  
 HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100  
 ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · Selangor  
 MEXICO · Sigma-Aldrich Quimica S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo  
 POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Quimica, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.  
 SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm  
 SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH  
 UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: [supelco@sial.com](mailto:supelco@sial.com)

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

BKF