

# TheReporter

Reprinted from Volume 14, No. 2, 1995

T295042

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

**ARGENTINA** • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119  
**AUSTRALIA** • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154  
**AUSTRIA** • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien  
**BELGIUM** • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem  
**BRAZIL** • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP  
**CANADA** • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8  
**CZECH REPUBLIC** • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8  
**DENMARK** • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallengbaek Strand  
**FINLAND** • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki  
**FRANCE** • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeaus Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex  
**GERMANY** • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen  
**GREECE** • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens  
**HUNGARY** • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest  
**INDIA** • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048  
**IRELAND** • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24  
**ISRAEL** • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100  
**ITALY** • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano  
**JAPAN** • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103  
**KOREA** • Sigma-Aldrich Korea • Samhan Camus Annex, 10<sup>th</sup> Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul  
**MALAYSIA** • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor  
**MEXICO** • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca  
**NETHERLANDS** • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht  
**NORWAY** • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo  
**POLAND** • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań  
**PORTUGAL** • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710  
**RUSSIA** • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062  
**SINGAPORE** • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529  
**SOUTH AFRICA** • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459  
**SPAIN** • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid  
**SWEDEN** • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm  
**SWITZERLAND** • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs  
**UNITED KINGDOM** • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH  
**UNITED STATES** • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email: [supelco@sial.com](mailto:supelco@sial.com)

H



Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

# Field Concentration of Toxic Organics in Air

M. Bruns and F. Li (MTI Analytical Instruments)

*This guest article describes analysis of airborne organics using MTI's new Trapper 2000 portable sample concentrator and custom adsorbent traps developed by Supelco. The Trapper 2000 concentrator, combined with the MTI high-speed micro gas chromatograph, allows for field concentration and analysis for on-site environmental assessment and remediation. This system accommodates two adsorbent traps with simultaneous operation, allowing for analytical duplication and increased sample throughput. The traps, with integral heaters, rechargeable 12-volt battery, and refillable purge gas supply, are field replaceable. The battery provides 6 hours of operation.*

The recent enactment of the 1990 Amendments to the US Clean Air Act and stricter enforcement of environmental regulations worldwide have challenged the development of portable instrumentation capable of field concentration (collection) and analysis of toxic organics in air (1). In response to this challenge, MTI Analytical Instruments (Fremont, CA, USA), along with the technical, manufacturing, and product support of Supelco, has developed a high-throughput, field-portable sample concentrator for use in conjunction with a portable gas chromatograph.

MTI's new Trapper 2000™ portable sample concentrator (Figure A) was developed for concentration of the 41 toxic volatile organic

compounds (VOCs) in the US Environmental Protection Agency's Compendium Method TO-14. This instrument, combined with the MTI high-speed micro gas chromatograph, allows field concentration and analysis for on-site environmental assessment and remediation (2). The system yields results comparable to laboratory purge-and-trap and GC/MS systems (3).

Air samples are drawn through and concentrated on the trap at ambient temperatures, using an internal pump. Immediately following adsorption, and prior to desorption of the sample, the trap undergoes a purge with dry helium to remove air and water, thus preventing oxidization of the sample or adsorbents. The trap is then rapidly heated to the user-defined desorption temperature, and the desorbed concentrated sample is backflushed with helium into a gas-tight syringe for chromatographic analysis. Finally, the entire system is flushed with helium, while at the desorption temperature, to ensure minimal sample carryover. The system is microprocessor controlled and provides onboard diagnostics, alerting the user of any potential problems. The user can modify up to 6 parameters per trap to optimize system operation.

The Trapper 2000 concentrator uses a custom-developed four-adsorbent bed thermal desorption tube containing Carbopack™ C, Carbopack B, Carboxen™ 1000, and Carboxen 1001 adsorbents to achieve a 100- to 1000-fold concentration of toxic VOCs in air.

**Figure A. MTI Trapper 2000 Concentrator and P200 Micro Gas Chromatograph**



Photo provided by MTI Analytical Instruments, 41762 Christy Street, Fremont, CA 94538 USA.

**Table 1. Analyte Recovery and Analytical Reproducibility**


Compound	Boiling Point (°C)	Original Concentration (ng/L)	Desorption Efficiency (%)	Percent RSD (n = 5)
Freon 12	-29.8	50	101	4.9
Vinyl chloride	-13.4	20	100	4.5
Ethyl chloride	12.3	120	92	5.7
Freon 11	23.7	20	92	6.6
Dichloromethane	39.8	60	90	8.1
Carbon tetrachloride	76.5	40	94	6.3
Benzene	80.1	40	90	4.2
Trichloroethylene	87.0	40	94	4.3
Toluene	110.6	40	96	4.6
Tetrachloroethylene	121.1	40	87	8.5
Chlorobenzene	132.0	40	90	2.7
Ethylbenzene	136.2	40	84	4.6
o-Xylene	144.4	80	83	4.8
o-Dichlorobenzene	180.5	30	87	5.0
Adsorption:		1000mL at 85mL/min at 30°C (sample saturated with moisture)		
Desorption:		2mL at 0.8mL/min at 350°C		
Adsorbents:		Carbopack C, Carbopack B, Carboxen 1000, Carboxen 1001		

## Figure B. Toxic VOCs in Air


Sample Conc.: adsorption — 1000mL at 85mL/min at 30°C  
 desorption — 2mL at 0.8mL/min at 350°C  
 adsorbents — Carboxen C, Carboxen B, Carboxen 1000,  
 Carboxen 1001  
 Column: OV<sup>®</sup>-1, 0.15mm ID, 1.2µm film  
 Oven: 53°C  
 Carrier: helium, column head pressure 21.7psi  
 Inj. Time: 100msec  
 Sample Time: 5 sec  
 Det.: micro thermal conductivity, medium gain

- |                                 |                                |
|---------------------------------|--------------------------------|
| 1. Ethylchloride (120ppm)       | 6. Trichloroethylene (40ppm)   |
| 2. Freon 11 (20ppm)             | 7. Toluene (40ppm)             |
| 3. Dichloromethane (60ppm)      | 8. Tetrachloroethylene (40ppm) |
| 4. Benzene (40ppm)              | 9. Chlorobenzene (40ppm)       |
| 5. Carbon tetrachloride (40ppm) | 10. Ethylbenzene (40ppm)       |


B1 — Concentrated Standard Mixture  
 of VOCs at ppm Level Prior to Dilution



B2 — First 2mL Desorption



B3 — Second 2mL Desorption



Figures provided by MTI Analytical Instruments, 41762 Christy Street, Fremont, CA 94538 USA.

Analyte adsorption is achieved without the use of cryogenics, which would hinder portability and use in the field. This particular adsorbent trap has been demonstrated to have a breakthrough volume of 5 liters of Freon<sup>®</sup> 12 (at 50ppm) in moisture-saturated helium at 40°C. It provides excellent analyte recovery and precision, yielding results comparable to full-scale laboratory techniques and instrumentation (Table 1).

Figure B1 illustrates a 500-fold dilution of a VOC sample, and Figure B2 a 500-fold concentration. The four-adsorbent bed trap effectively adsorbs and desorbs the trapped VOCs (Figure B2) to achieve the original concentration of the standard mixture (Figure B1). The absence of a detectable signal from the desorption volume in Figure B3 demonstrates the efficient desorption (with negligible carryover) of the adsorbent material with the first 2mL volume. (Note that the amplitude axis in Figure B3 is expanded.) No sample breakthrough (Figures B1 and B2) or carryover (Figures B2 and B3) from the initial desorption volume is detected.

If your work involves the analysis of airborne toxics, and field analysis would be helpful, consider the convenience of the MTI Trapper 2000 portable concentrator and the dependability of Supelco adsorbents. Supelco can customize adsorbent tubes according to your unique requirements.

### Ordering Information:

For more information about the Trapper 2000 Concentrator, the P200H Micro Gas Chromatograph, or other MTI products, contact MTI Analytical Instruments, 41762 Christy Street, Fremont, CA 94538 USA (phone 510-490-0900, FAX 510-651-2498).

For more information about custom packed adsorbent traps for VOCs analysis or any other application, phone Supelco at 800-247-6628 or 814-359-3441, FAX 800-447-3044 or 814-359-3044, or contact your nearest Supelco representative.

### References

1. Carney, K.R., R.L. Wong, E.B. Overton, M.A. Jackish, and C.F. Steele, *Sampling and Analysis of Airborne Pollutants*, E.D. Winegar, L.H. Keith, Eds., CRC Press, Inc., Boca Raton, FL, USA, 22-36 (1993).
2. Bruns, M.W., *Silicon Micromachining and High Speed Gas Chromatography*, in *Proceedings for the 1992 International Conference on Industrial Electronics, Control, Instrumentation, and Automation*, Institute for Electrical and Electronics Engineers: New York, NY, USA, Vol. 3, 1640-1644 (1992).
3. Carney, K.R., and E.B. Overton, article submitted to *American Laboratory* for publication in early 1995.

Reference 3 is available from MTI Analytical Instruments, 41762 Christy Street, Fremont, CA 94538 USA.

This article was submitted by Mark Bruns, Ph.D., MTI Analytical Instruments, 41762 Christy Street, Fremont, CA 94538 USA.

### Trademarks

Carboxen, Carboxen — Sigma-Aldrich Co.  
 Freon — E.I. du Pont de Nemours & Co., Inc.  
 OV — Ohio Valley Specialty Chemical Co.  
 Trapper 2000 — MTI Analytical Instruments