

TheReporter

Reprinted from Volume 15, No. 3, 1996

T296013

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119
AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154
AUSTRIA • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem
BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP
CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8
DENMARK • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallensbaek Strand
FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex
GERMANY • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen
GREECE • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens
HUNGARY • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest
INDIA • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048
IRELAND • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24
ISRAEL • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100
ITALY • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano
JAPAN • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103
KOREA • Sigma-Aldrich Korea • Samhan Camus Annex, 10th Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul
MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor
MEXICO • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca
NETHERLANDS • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht
NORWAY • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań
PORTUGAL • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710
RUSSIA • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062
SINGAPORE • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459
SPAIN • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid
SWEDEN • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm
SWITZERLAND • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs
UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.

 **SUPELCO**

ISO 9001
REGISTERED

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

Molecular Sieve PLOT Columns for the Capillary GC Separation of Permanent Gases

W. Betz, M. Keeler

PLOT capillary columns containing Carboxen carbon molecular sieves or molecular sieve 5A have sufficient capacity for directly injected, highly concentrated samples of permanent gases and C2 and C3 hydrocarbons. The pore size ranges and bonded layering of the particles provide effective separations of these analytes.

Porous solids, such as molecular sieves, are traditionally used in packed bed separation processes ranging from bulk-scale to bench-scale. Processes at the bench scale include gas-solid chromatography (GSC), using molecular sieves consisting of either carbon or aluminosilicate. Advancements in chromatographic techniques, including improvements in wall-coated capillary columns and molecular sieves, have enhanced efficiency and led to the development of a family of porous layer open tubular (PLOT) columns.

Carboxen-1006[®] and Molecular Sieve 5A PLOT capillary columns exhibit adsorption strengths that allow the separation of analytes with subambient boiling points. Table 1 illustrates the reproducibility of these two sets of columns.

The carbon molecular sieve particles in the PLOT layer (30 μ m) are spherical, and allow an extended path length through the intact particle. These particles also provide a large internal surface area


Table 1. Reproducibility Data for Molecular Sieve PLOT Columns

Description	k'	Standard Deviation
Carboxen-1006 (n = 24)	5.58 (CO ₂)	0.54
Mol Sieve 5A (n = 20)	4.97 (CO)	0.36

(725m²/g). A proprietary process is used to adhere these carbon particles to the walls of the column. This provides a stable porous solid surface, eliminates particle loss from the walls, and allows rapid temperature programming of the GC oven for fast-run analyses.

Macropores have a positive effect on the kinetic behavior of a porous solid such as Carboxen[™]-1006 adsorbent. The plots in Figure A indicate poor kinetic behavior of a micropores-only particle. The presence of macropores allows separation of analytes at linear velocities well above optimum.


Figure A. van Deemter Plot Shows the Effect of Porosity on Efficiency



796-0148

Figure B. Light Gases on a Carboxen-1006 PLOT Column

Column: Carboxen-1006 PLOT, 30m x 0.53mm ID
 Cat. No.: 25461
 Oven: 35°C (1 min) to 250°C at 24°C/min
 Carrier: helium, 10mL/min
 Det.: TCD, 250°C
 Inj.: 70 μ L, direct, 230°C



794-0618


Figure B illustrates the ability of the Carboxen-1006 molecular sieve to separate nitrogen, carbon monoxide, and other light gases. Note the excellent baseline resolution. Figures C and D show the use of the Carboxen-1006 column for the separation of formalin solutions and impurities in propylene.

The aluminosilicate PLOT column incorporates a molecular sieve 5A in the sub-micron particle size range. The process used to prepare this column is similar to that for the Carboxen-1006 column. Thus, similar optimization procedures are available. Figure E illustrates the use of the Molecular Sieve 5A column.

If you currently are using packed columns and traditional adsorbents for GSC separations of permanent gases or light hydrocarbons, consider the high efficiency of our PLOT capillary columns.

Figure C. Formalin Components on a Carboxen-1006 PLOT Column

Column: **Carboxen-1006 PLOT, 30m x 0.53mm ID**
 Cat. No.: **25461**
 Oven: 220°C
 Carrier: helium, 22cm/sec
 Det.: TCD, 230°C
 Inj.: 0.2µL formalin solution, direct, 230°C




1. Water
2. Formaldehyde
3. Methanol

794-0122

Figure D. Impurities in Propylene on a Carboxen-1006 PLOT Column

Column: **Carboxen-1006 PLOT, 30m x 0.32mm ID**
 Cat. No.: **24241-U**
 Oven: 200°C (4 min) to 250°C at 24°C/min
 Carrier: helium, 76cm/sec
 Det.: FID, 230°C
 Inj.: 20µL propylene, approx. 1% each contaminant, direct, 230°C




1. Methane
2. Acetylene
3. Ethylene
4. Ethane
5. Methylacetylene
6. Propylene
7. Propane

794-0359

Figure E. Small Molecules on a Molecular Sieve 5A PLOT Column

Column: **Mol Sieve 5A PLOT, 30m x 0.53mm ID**
 Cat. No.: **25463**
 Oven: 65°C
 Carrier: helium, 10mL/min
 Det.: TCD, 260°C
 Inj.: direct, 260°C



1. Hydrogen
2. Oxygen
3. Nitrogen
4. Methane
5. Carbon monoxide

795-0367

Ordering Information:

Description	Max. Temp. (°C)	Cat. No.
Carboxen-1006 PLOT Columns		
30m x 0.32mm ID	250	24241
30m x 0.53mm ID	250	25461
Mol Sieve 5A PLOT Columns		
30m x 0.32mm ID	300	24243
30m x 0.53mm ID	300	25463

■ US Pat. No. 4,839,331.

Carboxen is a trademark of Sigma-Aldrich Co.

