

# TheReporter

Reprinted from Volume 15, No. 6, 1996

T296036

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

**ARGENTINA** • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119  
**AUSTRALIA** • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154  
**AUSTRIA** • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien  
**BELGIUM** • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem  
**BRAZIL** • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP  
**CANADA** • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8  
**CZECH REPUBLIC** • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8  
**DENMARK** • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallensbaek Strand  
**FINLAND** • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki  
**FRANCE** • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex  
**GERMANY** • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen  
**GREECE** • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens  
**HUNGARY** • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest  
**INDIA** • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048  
**IRELAND** • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24  
**ISRAEL** • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100  
**ITALY** • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano  
**JAPAN** • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103  
**KOREA** • Sigma-Aldrich Korea • Samhan Camus Annex, 10<sup>th</sup> Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul  
**MALAYSIA** • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor  
**MEXICO** • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca  
**NETHERLANDS** • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht  
**NORWAY** • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo  
**POLAND** • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań  
**PORTUGAL** • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710  
**RUSSIA** • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062  
**SINGAPORE** • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529  
**SOUTH AFRICA** • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459  
**SPAIN** • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid  
**SWEDEN** • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm  
**SWITZERLAND** • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs  
**UNITED KINGDOM** • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH  
**UNITED STATES** • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.

 **SUPELCO**

ISO 9001  
REGISTERED

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

# Alumina PLOT Columns for Rapid and Consistent Separations of Light Hydrocarbons by Gas Solid Chromatography

W.R. Betz, M.J. Keeler, J.H. Whitford, S.B. Cole  
Supelco, Bellefonte, PA USA

*Incorporation of alumina, carbons, porous polymers, and other adsorbents into porous layer open tubular columns (PLOT columns) with relatively small diameters has considerably improved gas-solid chromatographic separations of many small analytes. 30-meter PLOT columns containing sub-micron, granular, activated alumina particles elute unsaturated hydrocarbons after the saturated hydrocarbons of the same carbon number. A patented adhesive technology immobilizes the adsorbent to the column wall – no particles are lost during routine analyses or rapid temperature programming. The high temperature limit of the adhesive, >360°C, makes the operating temperature limit of these columns a function of the particles used, not a function of the adhesive. The columns effectively separate mixtures of C1-C7 hydrocarbons.*

Carbons, porous polymers, alumina, and other adsorbents have been used for decades for gas-solid chromatographic separations of small analytes. The evolution of GSC to the degree that these adsorbents can be incorporated into porous layer open tubular columns (PLOT columns) with relatively small diameters has considerably improved resolution for many key analytes.

The choice of adsorbent depends on the properties of the analytes to be separated. Aluminum oxide provides a unique retention pattern for saturated and unsaturated light hydrocarbons, hence its application as a gas-solid adsorbent. We prepare 30-meter PLOT columns with sub-micron, granular, activated alumina particles with a surface area of 350m<sup>2</sup>/g. The use of sodium sulfate as a desiccant and activating agent ensures that the activity of the alumina surface is effective in eluting (distilling) unsaturated hydrocarbons after the saturated hydrocarbons of the same carbon number.


In preparing these PLOT columns we incorporate a proprietary, patented adhesive technology through which the adsorbent particles are immobilized, or cemented, to the wall of the capillary column. Because the adhesive also bonds the particles to one another, no particles are lost during routine analyses or rapid temperature programming. A significant feature of the adhesive is its high temperature limit, greater than 360°C, which makes the operating temperature limit of the columns a function of the particles used, not a function of the adhesive.

We also optimize column preparation to ensure that the adhesive does not interfere with the pores in the particles. This, in turn, ensures that analyte molecules are transferred, unrestricted, to the internal micropores in the particles, and that there is no adverse


effect on the external working surface. We have used this adhesive technology to prepare PLOT columns containing carbon, Tenax®, Porapak® N, U, and T, HayeSep® A, C, D, N, and T, polymer G-45, Mol Sieve 5X and 13X, and alumina/KCl adsorbents.

C1 to C4 hydrocarbons are separated effectively by a 30-meter alumina PLOT column, using a temperature program (Figure A). The elution pattern desired for analyses of trace impurities in unsaturated hydrocarbons – trace impurities eluting before major peaks – is a significant benefit of this column. The C4 isomers in

**Figure A. C1-C4 Hydrocarbons on an Alumina PLOT Column**



**Figure B. Natural Gas Standard (C4 and C5 Hydrocarbons)**



Reprinted from Volume 15, No. 6, 1996


American Society for Testing and Materials (ASTM) natural gas reference standard blend 6 also are effectively separated, using the same column and temperature program (Figure B).

Mixtures containing C1 to C7 hydrocarbons are analyzed by programming the column to 200°C. Figure C shows an analysis of a calorimetric natural gas reference standard comprised of two permanent gases (nitrogen and carbon dioxide, which are not shown), methane (the major component of the sample), and ethane, propane, isobutane, n-butane, isopentane, neopentane, n-hexane, and n-heptane.

Our PLOT columns prepared with carbon molecular sieves, zeolite molecular sieves, porous polymers, and activated alumina provide options for a wide range of analyses, including permanent gases, light hydrocarbons, Freon® gases, sulfur gases, and many polar compounds. Our proprietary adhesive process eliminates particle loss from the column wall and allows rapid temperature programs to be used. Further, our production processes ensure reproducible performance from column to column. We recommend our 30-meter alumina PLOT columns for GSC separations of small hydrocarbons.

### Figure C. Natural Gas Standard (C1-C7 Hydrocarbons)

Column: Alumina PLOT, 30m x 0.53mm ID  
 Cat. No.: 25464  
 Oven: 35°C (2.5 min) to 200°C at 5°C/min  
 Carrier: helium, 4mL/min  
 Det.: FID, 200°C  
 Inj.: natural gas standard (Cat. No. 3-03-101), direct injection, 200°C



96-0298

### Ordering Information:

Description	Cat. No.
<b>Alumina PLOT Columns</b>	
30m x 0.32mm	24244
30m x 0.53mm	25464
<b>Light Hydrocarbon Mixes*</b>	
C1-C4 Straight Chain	22566
C4 Saturates/Unsaturates	22567
<b>Natural Gas Standard**</b>	<b>303101</b>

For PLOT columns of other dimensions, please inquire.

For an extensive selection of other chemical standards for light hydrocarbons analyses, refer to the Supelco catalog.

\*10-20ppm each analyte in nitrogen in a 4-liter SCOTTY® I cylinder. Cat. No. 22566 includes the 8 compounds in Figure A; Cat. No. 22567 is similar to the mix in Figure B, but includes ethylacetylene (1-butyne), does not include n-pentane or isopentane. For 37-liter cylinders, please see our general catalog.

\*\*Includes the components in Figure C, plus helium, nitrogen, and carbon dioxide (which do not elicit a response from a flame ionization detector), in a 14-liter SCOTTY II cylinder. For mole percent composition, please see our general catalog.

### Trademarks

HayeSep – Hayes Separations Inc.  
 Porapak – Waters Associates, Inc.  
 SCOTTY – Scott Specialty Gases, Inc.  
 Tenax – Enka Research Institute Arnhem



Reprinted from Volume 15, No. 6, 1996