



Volatile amines

Application Note

Materials Testing & Research

Authors


Agilent Technologies, Inc.

Introduction

Gas chromatography using an Agilent CP-Wax 51 for Amines column separates nine volatile amines in 20 minutes.



Agilent Technologies



Conditions

Technique : GC-capillary
Column : Agilent CP-Wax 51, 0.32 mm x 10 m fused silica
WCOT CP-Wax 51 for amines (0.2 μ m)
(Custom-made)
Temperature : 70 °C (1 min) \rightarrow 150 °C, 5 °C/min
Carrier Gas : N₂, 40 kPa (0.4 bar, 5.8 psi), 9.3 cm/s
Injector : Splitter, 40 mL/min
Detector : FID, 8 x 10⁻¹² Afs

Peak identification

1. diisopropylamine
2. 2-aminobutane
3. diisobutylamine
4. 1-aminopentane
5. dibutylamine
6. pyridine
7. decylamine
8. dicyclohexylamine
9. 2,6-dimethylaniline

www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A00089



Agilent Technologies